Waters

Determination of Fluoroquinolone Residues in Bovine Kidney Using SPE and LC/MS/MS

Kevin M. Jenkins, Michael S. Young Waters Corporation, 34 Maple Street, Milford, MA 01757 USA

OVERVIEW

An LC/MS method is described for the trace quantification of nine fluoroquinolone (FLQ) residues (ciprofloxacin, enrofloxacin, sarafloxacin, enofloxacin, flumequine, norfloxacin, ofloxacin, lomefloxacin and danofloxacin) in bovine kidney. The drug residues were extracted from a 2 g kidney sample using an aqueous buffer and the sample extracts were cleaned and enriched using mixed-mode anion exchange solid phase extraction (Oasis® MAX) for sample clean-up and enrichment. All nine FLQ's can be resolved chromatographically using isocratic mobile phase conditions. Positive ion atmospheric pressure chemical ionization (APCI+) mass spectrometry was used to quantify and confirm the parent ion $[M+H]^+$ and daughter fragments for each analyte.

Fluoroquinolones are a synthetic class of antibiotics that are highly effective against gram-positive and gram-negative bacteria. This class of antibiotics was originally developed in 1962 for human use; however, in 1995 selected drugs from this family of antibiotics were approved for animal husbandry in North America. Now, there is concern because of reports that human diseases have developed antibacterial resistance to these drugs. Therefore, the FDA has recently banned the extralabel use of the fluoroquinolones in food-producing animals.

This presentation outlines a protocol for the effective clean-up and extraction of FLQ residues in bovine kidney and muscle. For regulatory purposes, the method was optimized for LC/MS/MS analysis to provide sensitivity and selectivity.

MS CONDITIONS

HPLC CONDITIONS

LC/MS/MS Compound MW MRM Column: **Atlantis™ dC**₁₈ 4.6 x 150 mm 5 μm Mobile Phase: A: 0.2% NFPA* in water flumequin 262→244 B: Methanol enoxacin 320 321→303 Gradient: Linear 10% B to 80% B in 10 min 210 200 \200

FLUOROQUINOLONE STRUCTURE



- Amphoteric behavior •pKa of acid functionality » 5 •pKa of basic functionality » 8-9
- Can be extracted using ion-exchange sorbents
 - •At pH below 6 using cation-exchange •at pH above 6 using anion-exchange

LC SEPARATIONS

Collsion

Energy

(eV)

20

20

Cone

(V)

50

50



Flow: 0.8 mL/min
Temperature: 30°C
Injection: 30-80 µL

LC-UV/Fluorescence/Single Quad MS Column: Atlantis[™] dC₁₈ 4.6 x 150 mm 5 µm Mobile Phase: 75% NFPA* (0.2%) in water 3% Methanol 22% Acetonitrile Flow: 1.2 mL/min

Temperature: 30°C

* CF₃-CF₂-CF₂-COOH nonafluoropentanoic acid (NFPA)

SPE CONDITIONS

Initial Tissue Extraction: Extract 2 g of kidney with 30 mL of 50 mM phosphate buffer (pH 7.4) centrifuge at 10 000 rpm for 10 minutes.



norfloxacin	319	320→302	50	23
Sarafloxacin	385	386→368	50	25
ofloxacin	361	362→344	50	20
Enrofloxacin	359	360→342	50	20
danofloxacin	357	358→340	50	25
lomefloxacin	351	352→334	50	20
ciprofloxacin	331	332→214	50	20

APCI+ or ESI+ can be used for FLQ analysis in the absence of matrix, ESI+ may be somewhat more sensitive. APCI+ was chosen to minimize matrix suppression effects and the suppression effect of the ion-pairing reagent. Matrix effect was < 20% for all compounds.

Considerations for Initial Tissue Extraction

Fluoroquinolone antibiotics can be extracted from tissue at low, high or neutral pH.

- At low pH, extraction into ethanol is preferred; higher matrix interference results from aqueous buffer extraction
- At high pH extraction is difficult using either ethanol or aqueous buffer; high amounts of matrix interference
- At neutral pH, extraction into aqueous buffer is preferred; higher matrix interference results from ethanol extraction

Neutral pH extraction into aqueous buffer was chosen for this study.

Figure 1. Isocratic separation of nine fluoroquinlones using LC/UV.



9.32 flumequine 262.2 > 244.1 7.67 321.2 > 303.1 enoxacin 7.87 norfloxacin 320.2 > 302.1 8.46 sarafloxacin 386.2 > 368.2 7.65 ofloxacin 362.2 > 344.1 8.06 enrofloxacin 360.2 > 342.1 8.02 danofloxacin 358.2 > 340.1 8.12 lomefloxacin 352.2 > 334.1 7.98 ciprofloxacin 332.2 > 314.1 11

Figure 2. Gradient separation of nine fluoroquinlones using LC/MS/MS.

RESULTS

Table 1. Recovery data for nine fluoroquinolones using LC/MS/MS.

	5 ppb 10 ppb				25 ppb			50 ppb			75 ppb			100 ppb				
Analyte	Recovery	STD	RSD	Recovery	STD	RSD	Recovery	STD	RSD	Recovery	STD	RSD	Recovery	STD	RSD	Recovery	STD	RSD
	(%)	Dev	(%)	(%)	Dev	(%)	(%)	Dev	(%)	(%)	Dev	(%)	(%)	Dev	(%)	(%)	Dev	(%)
Cipro	81.3	11.5	14.1	70.2	8.0	11.3	62.7	3.6	5.7	64.0	5.0	7.8	72.0	5.4	7.4	80.5	3.8	4.7
Lome	68.3	17.3	25.4	76.9	5.0	6.5	81.6	6.0	7.3	70.9	6.4	9.0	84.3	5.3	6.3	87.8	4.7	5.3
Dano	86.2	16.8	16.5	75.5	5.8	7.7	75.9	4.2	5.5	66.6	5.0	7.5	76.4	7.0	9.1	77.5	5.6	7.3
Enro	79.2	11.6	14.6	73.3	15.6	21.3	82.4	3.0	3.6	70.6	5.4	7.7	85.6	5.7	6.7	108.4	6.1	5.6
Oflo	77.0	7.6	9.8	85.3	5.3	6.2	77.1	5.8	7.5	72.7	5.1	7.0	81.5	5.7	6.9	91.9	7.1	7.7
Flum	70.9	12.0	17.0	75.0	10.1	13.5	75.2	7.2	9.6	65.9	7.4	11.2	67.4	4.9	7.3	68.0	4.6	6.8
Sara	87.0	11.6	13.3	71.9	7.8	10.9	91.1	4.0	4.3	69.6	4.7	6.7	82.8	4.3	5.2	82.0	4.5	5.5
Nor	70.1	12.0	17.1	70.2	5.8	8.2	61.8	4.6	7.5	62.9	4.9	7.9	67.2	5.6	8.3	75.7	4.5	5.9
Eno	63.9	7.6	12.0	69.4	7.7	11.0	59.3	6.0	10.1	63.0	6.2	9.9	67.5	6.9	10.3	78.7	6.0	7.6

FUTURE CONSIDERATIONS

Optimization of SPE for LC/Fluorescence



Neutralize with formic acid and

bring to volume with mobile phase buffer.

A dual cartridge SPE clean-up improves sample cleanliness and lowers detection limits. The following two chromatograms show the improved clean-up.



Figure 3. Top figure: The chromatogram resulting from a dual retention SPE sample clean-up using fluorescence detection and an

CONCLUSIONS

- Improved SPE protocol for bovine kidney
 - •Reproducible 70% recovery
 - •LOQ below 10 µg/kg
 - •Fast, straightforward method
- Improved HPLC separation •NFPA allows for near baseline separation for LC or LC/UV
 - •NFPA gives good separation from matrix for LC/MS/MS
- Detection independent method •Allows for MS/MS confirmation Ideal for fluorescence detection





