

Waters

Mass Spectrometry Systems



ISPPP 2003

November 12, 2003

The Impact of Column Peak Capacity on the Multi-dimensional Chromatography of Complex Peptide Mixtures

Martin Gilar

Life Sciences Chemistry, Waters Corporation




- Highly complex peptide samples
- Column peak capacity definition
- 1D RP-HPLC - where are the limits?
- Peak capacity prediction
- Peak capacity in 2D and MD-HPLC
- Productivity of MD separation
- Future development

- Proteomic samples = peptides
- 10-50 thousand proteins
- 100-500 thousand peptides
- Current MD-HPLC systems: peak capacity > 10 000
 - orthogonal selectivity - 1D x 2D x ...MD
 - SCX and RP-HPLC
- MS = additional separation dimension

Column peak capacity: Maximum number of peaks that can be separated on a column within a given gradient time.

$$P = 1 + \frac{t_g}{w}$$




1D HPLC peak capacity

Column type Lx ID (mm)	Particle size (mm)	Void volume (ml)	Gradient time (s=0.0123)	Gradient time (s=0.0246)	Gradient time (s=0.0492)
Symmetry300 C18					
50 x 4.6	5	0.66	37.5	17.8	8.9
Symmetry300 C18					
150 x 4.6	5	1.86	100	50	25
Symmetry300 C18					
150 x 4.6	3.5	1.90	100	50	25
Symmetry300 C18					
150 x 4.6	7	1.80	100	50	25
Symmetry300 C18					
300 x 4.6	5	3.46	186	93.1	46.5
			Gradient time	Gradient time	Gradient time
PolySULFOETHYL aspartamide 50 x 4.6	5	0.7	20	40	80

1D HPLC peak capacity

Symmetry300 C18, 150 x 4.6 mm, 5 μ m






Strategies for increasing the column peak capacity

- Long gradient time (shallow gradients)
- Extend the column length + gradient time
- Use small particle size sorbent

$$s = \Delta c \cdot \frac{t_0}{t_g} = 0.5 \cdot \frac{1}{50} = 0.01$$

Peak capacity prediction


$$P = 1 + \frac{\sqrt{N}}{4} \cdot \frac{B \cdot \Delta c}{B \cdot \Delta c \cdot (t_0/t_g) + 1}$$




Peak capacity prediction

$$P = 1 + \frac{\sqrt{N}}{4} \cdot \frac{B \cdot \Delta c}{B \cdot \Delta c \cdot (t_0/t_g) + 1}$$


Column efficiency N , $D_m = 3 \times 10^{-10} \text{ m}^2/\text{s}$




$$P = 1 + \frac{\sqrt{N}}{4} \cdot \frac{B \cdot \Delta c}{B \cdot \Delta c \cdot (t_0/t_g) + 1}$$



Five columns, 3 gradient slopes


 $s = 0.0123$ $s = 0.0246$ $s = 0.0492$ $B=55$ $D_m = 3 \times 10^{-10} \text{ m}^2/\text{s}$

Peak capacity prediction



150 x 4.6 mm column,
decreasing d_p

- Orthogonal selectivity - P multiplication
- SCX + RP-HPLC



- SCX peak capacity ~ 50 % of RP-HPLC values or less (polysulfoethyl aspartamide, PolyLC Inc.)
- SCX step gradient elution - peak capacity = no. of fractions

Column L x ID (SCX)	Gradient time (minutes)	Peak capacity
50 x 4.6	20	65
50 x 4.6	40	83
50 x 4.6	80	113

- SCX step gradient elution - extreme demand on 2D
- Fast 1D, very long 2D
- low peak/min productivity
- extensively long total analysis time
- complex sample enters MS
- limited number of identified peptides / dynamic range
- great demand on data interpretation (software)

2D-HPLC peak capacity

1D fractions SCX	RP column length (mm)	2D grad time RP (minutes)	Peak capacity	Time (hours)
10	150	60	2700	10
10	250	120	3800	20
10	500	240	5300	40
40	150	60	10800	40



* 5 μ m sorbent

Productivity?

Productivity of separation

Productivity =
peaks / min

Best for short,
efficient columns
using fast gradient



2D-HPLC productivity

1D fractions SCX	RP column length (mm)	2D grad time RP (minutes)	RP-HPLC peak capacity	Total peak capacity	Time (hours)	Productivity Peaks/min.
20	250 (3.5 μm)	240	657	13000	80	2.7
40	150 (3.5 μm)	60	354	14000	40	5.9
80	50 (3.5 μm)	15	175	14000	20	11.6
80	20 (1.8 μm)	6	178	14000	8	29.6

Short + efficient columns = best productivity

- Uses assumption of porous sorbents
- D_m , and $B = ?$ Changes with peptide MW.
- Model overestimates P for small d_p and fast gradients
- Extra-column peak broadening?
- Developed for 4.6 mm I.D. (lower P for nano-columns)
- Peak capacity \neq no. of separated peaks
 - 100 components injected \rightarrow 37 observed peaks (Giddings)
 - 1D & 2D will resolve $\sim 37\%$ components (assuming that no. of components injected = peak capacity of the system)

- RP-HPLC peak capacity ~ 200-400
- SCX-HPLC peak capacity ~50-100
- 2D-HPLC peak capacity
 - ▶ on-line ~5000-10,000
 - ▶ off-line ~10,000-40,000
- limited productivity of current setups
- long analysis time (second dimension)
- What is the expected sample complexity?
- What is the best pre-fractionation?

- Protein pre-fractionation prior to 2D-HPLC (peptides)
- Abundant protein depletion (serum proteomics)
- ~ 1 μm particles
- ultra-HPLC
- monolithic columns
- New orthogonal modes of 2D-HPLC
- Progress in MS/MS and data acquisition and handling
- Fast 2D-HPLC systems ~10,000 peak capacity



Amy E. Daly

Mariana Kele

Uwe D. Neue

Tad Dourdeville

Kenneth J. Fountain

Ying-Qing Yu

John C. Gebler