

*Diane Diehl, Kim Van Tran, Eric S. Grumbach, Jon Belanger, Bob Brennick, Michael Savaria, Uwe Neue, Jeffrey Mazzeo*  
Waters Corporation, 34 Maple St., Milford, MA 01757 USA

## OVERVIEW-

Pharmaceutical companies have been driven to cut costs by increasing sample throughput. Because of the mandate to run faster separations, researchers have tried to move away from hour-long separations on 15 or 25 cm columns. We show that by using column hardware in 20 mm lengths, packed with a variety of spherical particle packing materials, one can increase sample throughput. Columns packed with spherical particles can be used to run separations in under 5 minutes with flow rates of 3 mL/min or less without any compromise in chromatographic performance. We also present how an application run on a long column can be transferred to these shorter columns resulting in optimum chromatographic performance in a much shorter run time, as well as resulting in cost savings. The short columns can be used for quick method scouting runs, as well as for LC/MS applications without flow splitting. Hundreds of injections can also be made on the 20 mm length columns.

## INTRODUCTION-

Using current HPLC technology, we know that by optimizing the column length, particle size and the flow rate/temperature of the separation, we can achieve faster separations. One means of achieving high-throughput separations was to use guard cartridges — i.e. very short columns. However, these cartridges were never intended to be used for chromatography and therefore can suffer from poor peak shape. Therefore, we designed columns in 20 mm lengths in 2.1, 3.0, 3.9 and 4.6 mm inner diameters for optimum chromatographic performance and fast separations. These columns are available in the chemistries listed below. We have used these columns to scale down separations, to develop new separations for both LC/UV and LC/MS applications. We will show examples of a scaled-down separations and the cost savings that you can achieve.


The Intelligent Speed (*IS*<sup>TM</sup>) columns are available in the following chemistries:

XTerra<sup>®</sup> MS C<sub>18</sub>  
XTerra<sup>®</sup> RP<sub>18</sub>  
XTerra<sup>®</sup> MS C<sub>8</sub>  
XTerra<sup>®</sup> RP<sub>8</sub>  
Atlantis<sup>™</sup> dC<sub>18</sub>  
Symmetry<sup>®</sup> C<sub>18</sub>  
SymmetryShield<sup>™</sup> RP<sub>18</sub>  
Symmetry<sup>®</sup> C<sub>8</sub>  
SymmetryShield<sup>™</sup> RP<sub>8</sub>



## MAXIMIZING PEAK CAPACITY-

Peak capacity is a measure of the separation power of a gradient on a particular column.



Peak capacity is affected by:

- Gradient duration ( $t_g$ )
- Flow rates ( $F$ )
- Column length ( $L$ )
- Particle size ( $d_p$ )

Start with the simple equation for peak capacity



$$P = 1 + \frac{t_g}{w}$$

Make a few substitutions

$$P = 1 + \frac{\sqrt{\frac{L}{a \cdot d_p + b \cdot D_M \cdot \frac{t_0}{L} + \frac{d_p^2}{c \cdot D_M} \cdot \frac{L}{t_0}}}}{4} \cdot \frac{B \cdot \Delta c}{B \cdot \Delta c \cdot \frac{t_0}{t_g} + 1}$$


We can now generate a 3-dimensional plot to examine how the gradient time and flow rate effect the peak capacity of a separation.

Neue, U. D., Mazzeo, J. R. J. Sep. Sci. 2001, 24, 921-929.  
Cheng, Y-F., Lu, Z., Neue, U. Rapid Commun. Mass Spectrom. 2001, 15, 141-151.



We see above that for a 32 minute gradient run on a 4.6 x 20 mm  $IS^M$  column, a flow rate of 1.1 mL/min is predicted to maximize the peak capacity. But, for a 4 minutes gradient, a fast flow rate of 2.6 mL/min is predicted for maximum peak capacity.

We can now run experiments to test the theory. We observe that a flow rate of 3 mL/min offers the best solution as far as peak capacity, solvent consumption and fast separations.



Recommended starting flow rates, to be adjusted according to the separation.

| Column Dimensions | Flow Rate   |
|-------------------|-------------|
| 4.6 x 20 mm       | 3.0 mL/min  |
| 3.9 x 20 mm       | 2.16 mL/min |
| 3.0 x 20 mm       | 1.28 mL/min |
| 2.1 x 20 mm       | 0.63 mL/min |

## SCALING SEPARATIONS

To successfully scale separations run on long columns to the 20 mm length columns, we use the following equation to scale the gradient:

■ To scale a gradient

$$\frac{L_2}{L_1} \times t_{g1} = t_{g2}$$

$L_1$  = Long column length

$L_2$  = Short column length

$t_{g1}$  = Gradient time on long column

$t_{g2}$  = Gradient time on short column

## Conditions

Column: XTerra® MS C<sub>18</sub>, 4.6 x 150 mm, 5 µm

Mobile Phase A: Water

Mobile Phase B: Acetonitrile

Mobile Phase C: 100 mM NH<sub>4</sub>HCO<sub>3</sub>, pH 10

Flow Rate: 1.4 mL/min

## Gradient:

| Time<br>(min) | %A | %B | %C |
|---------------|----|----|----|
| 0.0           | 80 | 10 | 10 |
| 20.0          | 50 | 40 | 10 |
| 21.0          | 80 | 10 | 10 |
| 25.0          | 80 | 10 | 10 |

Injection Volume: 10 µL

Sample concentration: 20 µg/mL

Temperature: 30 °C

Detection: UV @ 254 nm


Sampling Rate: 5 pts/sec

Filter: 0 (no filter)

Instrument: Alliance® 2695 with 996 PDA

## Analytes:


1. Caffeine
2. Aniline
3. N-Methylaniline
4. 2-Ethylaniline
5. 4-Nitroanisole
6. N-N-Dimethylaniline



To scale the separation to a 20 mm column, we first scale the gradient using the previous equation:

$$\frac{20 \text{ mm}}{150 \text{ mm}} \times 20 \text{ min} = 2.7 \text{ min}$$

We ran a 2.7 minute gradient on a 4.6 x 20 mm IS<sup>TM</sup>, 3.5 µm column at 3 mL/min, for a total cycle time of 4 minutes – down from 25 minutes! We further optimized the method to a cycle time of 3 min.



## COST SAVINGS: CALCULATIONS

Assume that an HPLC is running about 67% of the year, or 4000 hr. We used the previous separation as an example:

|                           | 4.6 x 150 mm | 4.6 x 20 mm <i>IS</i> <sup>TM</sup> |
|---------------------------|--------------|-------------------------------------|
| Cycle time (min)          | 25           | 3                                   |
| # of Samples Run per Year | <b>9600</b>  | <b>80,000</b>                       |

This is 8.3x as many samples run in one year, on ONE HPLC by using an intelligently designed column!

We can now calculate the solvent cost savings assuming that 5000 samples will be analyzed for the study.


|   | 4.6 x 150 mm          | 4.6 x 20 mm <i>IS</i> <sup>TM</sup> |
|---|-----------------------|-------------------------------------|
| Cycle time (min)                          | 25                    | 3                                   |
| Total time for 5000 samples (hours)       | <b>2083 (87 days)</b> | <b>250 (10 days)</b>                |
| Flow rate (mL/min)                        | 1.4                   | 3                                   |
| Total solvent consumption (L)             | <b>175</b>            | <b>45</b>                           |
| Amount ACN Consumed (L)<br>(~ 25% is ACN) | 43.75                 | 11.25                               |
| Cost for ACN (\$42.50/L)                  | \$1860                | \$480                               |
| Cost for waste disposal (~\$2.50/L)       | \$438                 | \$113                               |
| Total solvent costs                       | <b>\$2298</b>         | <b>\$593</b>                        |

For even further cost savings, we can scale the inner diameter from a 4.6 to a 2.1 mm and calculate the savings:

■ To scale a flow rate for different internal diameters

$$\frac{(d_2)^2}{(d_1)^2} \times F_1 = F_2$$

$d_1$  = Diameter of original column  
 $d_2$  = Diameter of second column  
 $F_1$  = Flow rate on original column  
 $F_2$  = Flow rate on second column




2.1 x 20 mm *IS*<sup>TM</sup>, 3.5  $\mu$ m ■ Scale from a 4.6 mm i.d. to a 2.1 mm i.d. column.

■ Flow rate scales from 3 mL/min to **0.6 mL/min.**

|   | 4.6 x 150 mm  | 2.1 x 20 mm <i>IS</i> <sup>TM</sup> |
|---|---------------|-------------------------------------|
| Flow rate (mL/min)                        | 1.4           | 0.6                                 |
| Total solvent consumption (L)             | <b>175</b>    | <b>9</b>                            |
| Amount ACN Consumed (L)<br>(~ 25% is ACN) | 43.75         | 2.25                                |
| Cost for ACN (\$42.50/L)                  | \$1860        | \$96                                |
| Cost for waste disposal (~\$2.50/L)       | \$438         | \$23                                |
| Total solvent costs                       | <b>\$2298</b> | <b>\$119</b>                        |

Cost savings of \$2179!

## LC/MS SEPARATIONS-

ConditionsColumn: Atlantis<sup>TM</sup> dC<sub>18</sub>, 2.1x 20 mm / IS<sup>TM</sup>, 3  $\mu$ mMobile Phase A: H<sub>2</sub>O

Mobile Phase B: MeOH

Mobile Phase C: 1% HCOOH in H<sub>2</sub>O

Flow Rate: 0.4 mL/min

| Gradient | Time (min) | %A | %B | %C |
|----------|------------|----|----|----|
|          | 0.0        | 50 | 40 | 10 |
|          | 1.0        | 30 | 60 | 10 |

Injection Volume: 2  $\mu$ LSample Concentration: 10  $\mu$ g/mLTemperature: 30  $^{\circ}$ CInstruments: Alliance<sup>®</sup> 2795 and Waters ZQ<sup>TM</sup>

| Analytes       | MW    |
|----------------|-------|
| Cinoxacin      | 262.2 |
| Oxolinic Acid  | 261.2 |
| Nalidixic Acid | 232.2 |

With a 2.1 mm i.d. column,  
lower flow rates enable DIRECT  
flow into the MS and the  
separation is only 1 minute long!

No flow splitting!

## CONCLUSIONS-

- Shorter analysis times can be achieved with shorter particulate columns with excellent results
  - Significant cost savings
  - Good column lifetimes
- New 20 mm length Intelligent Speed (IS<sup>TM</sup>) columns are available
  - Scale-down longer methods
  - Perform Method Development
  - Simplify LC/MS
- Waters can assist with a complete solution:
  - IS<sup>TM</sup> columns
- Choose from a wide range of chemistries and configurations
  - Method development technical support
  - Waters Instrumentation and Service

WATERS CORPORATION  
34 Maple St.  
Milford, MA 01757 U.S.A.  
T: 508 478 2000  
F: 508 872 1990  
www.waters.com

Waters, Micromass, Alliance, Symmetry, SymmetryShield, XTerra and Atlantis are trademarks of Waters Corporation.  
All other trademarks are the property of their respective owners.  
©2003 Waters Corporation

