

Marianna Kele, Uwe Neue, Kevin Wyndham, Pam Iraneta, Jeffrey Mazzeo and Tom Walter
Waters Corporation, Milford, MA 01757, USA




INTRODUCTION




During the last few years instrument and column manufacturers introduced several new approaches to improve performance in HPLC practice. Instrumentation for high temperature isothermal and thermal gradient LC and recently an extended pressure limit (15000 PSI, 1000 atm), low dispersion Ultra Performance Liquid Chromatography (UPLC) instrument became available. A new generation of small particle packed columns (<2.0 µm) are now available from several manufacturers to accommodate the need for shorter analysis time and better resolution.

This presentation gives an overview of different liquid chromatography techniques. It compares performance that can be achieved under traditional HPLC, High Temperature LC and high pressure, Ultra Performance Liquid Chromatography (UPLC) conditions. Advantages and challenges of these techniques are compared. The physico-chemical properties of small particles and the chromatographic performance of small particle packed columns are presented.

COMPARISON OF POROUS SMALL PARTICLES

Physico-chemical Properties


| | Particle Size [µm] | Particle Size Distribution (90/10 ratio) | Specific Surface Area [m²/g] | Average Pore Diameter [Å] | Specific Pore Volume [cm³/g] | Particle Shape |
|---------|--------------------|--|------------------------------|---------------------------|------------------------------|--|
| Brand 1 | | | | | | |
| A | 1.44 | 1.93 | 174 | 109 | 0.53 |  |
| Brand 2 | | | | | | |
| B | 1.75 | 1.58 | 185 | 130 | 0.68 |  |
| Brand 3 | | | | | | |
| C | 1.79 | 1.68 | 196 | 67 | 0.36 |  |

| | Particle Size [µm] | Particle Size Distribution (90/10 ratio) | Specific Surface Area [m²/g] | Average Pore Diameter [Å] | Specific Pore Volume [cm³/g] | Particle Shape |
|---------|--------------------|--|------------------------------|---------------------------|------------------------------|---|
| Brand 4 | | | | | | |
| D | 1.44 | 2.13 | 188 | 101 | 0.52 |  |
| Brand 5 | | | | | | |
| E | 1.66 | 2.71 | 201 | 85 | 0.51 |  |
| Brand 6 | | | | | | |
| F | 1.70 | 2.70 | 182 | 135 | 0.66 |  |




COMPARISON OF LC TECHNIQUES

| | Current Limits | Particle Size [µm] | Optimum Reduced Linear Velocity (chromatographic) | Interstitial Porosity | Specific Permeability [m ²] | Optimum Separation Impedance | Number of Theoretical Plates | Number of Theoretical Plates Produced in Unit Time at k'=4 N/t | Advantages | Disadvantages |
|--|---|---------------------|---|-----------------------|--|------------------------------|------------------------------|--|---|---|
| | | | | | B | E | N | | | |
| HPLC Porous Particle Packed Columns | 6000 PSI 3000 PSI* 10 ml/min 10 Hz | 1.8-10 | 3-5 | 0.38-0.42 | 4x10 ⁻¹⁵ -8x10 ⁻¹⁴ | 3000-4500 | 25000 | 800 | Wide selection of packing materials with different chemistry Multiple vendors of instruments | Slow |
| HPLC Porous Monolith Columns* | | 1-2 (skeleton size) | | 0.6-0.8 | 7x10 ⁻¹⁴ -1.2x10 ⁻¹³ | 500-1000 | 25000 | 3000 | Fast analysis Minimum sample clean up | Limited availability of stationary phases 4 µm particle packed column performance; Solvent consumption |
| High Temperature (60-200 °C) HPLC Porous Particle Packed Columns | | 3-10 | | 0.38-0.42 | 7x10 ⁻¹⁵ -1.2x10 ⁻¹³ | 3000-4500 | 25000 | | High efficiency separations | Limited availability of stable stationary phases Analyte degradation |
| High Pressure LC (UPLC) Porous Particle Packed Columns | 15000 PSI 2 ml/min 40 Hz | 1.5-2.0 | | 0.36-0.40 | 2-3.5x10 ⁻¹⁵ | 3500-5500 | 25000 | 9000 | High efficiency fast separations | Limited availability of stationary phases and packed columns Single vendor of instrument |


Chromatography Performance of 1.7, 3.5 and 5 µm Particle Packed Columns and Monolith



Columns: 2.1x50 mm ACQUITY UPLC™ BEH C18; 2.1x50 mm 3.5 µm XTerra MS C18; 2.1x50 mm 5.0 µm Symmetry C18; 4.6x50 mm Chromolith C18 (Merck) Solvent: Acetonitrile-Water 70:30 (v/v) Temperature: 25.0 °C; Solute: Hexylbenzene




Stability of 1.7 µm Particle Packed Columns Under Fast Gradient Conditions




Column: 2.1x100 mm ACQUITY UPLC™ BEH C18; Gradient: from 5% AcN to 95% AcN in 0.5 min, at 95% AcN for 1.2 min, to 5% AcN in 0.1 min; Maximum pressure: 11000 PSI Temperature: 30 °C; Solutes: alkylbenzenes; Instrument: Acquity UPLC (Waters)

Sensitivity to Extra-column Contribution




Chromatograms Obtained on 1.7, 3.5 and 5.0 µm Particle Packed and on a Monolith Column at Optimum Flow Rate (F_{opt})



Columns: 2.1x50 mm ACQUITY UPLC™ BEH C18; 2.1x50 mm 3.5 µm XTerra MS C18; 2.1x50 mm 5.0 µm Symmetry C18; 4.6x50 mm Chromolith C18 (Merck)

Chromatographic Performance



Column size: 2.1x50 mm (except Manufacturer A, Column2: 2.1x100 mm; Solvent: Acetonitrile-Water 70:30 (v/v); Temperature: 30 °C; Solute: Propylbenzene

SUMMARY

Experimental data demonstrate the advantages of using small particles (<2.0 µm) in packed columns. The efficiency and peak capacity in combination with speed of separation can be improved under UPLC conditions compared to traditional HPLC. Reduced plate heights (2.0-2.5) typical for 3-10 µm particle packed columns could be routinely achieved in high pressure applications. A limited selection of stationary phases with < 2.0 µm particle size is available at this point on the market. The performance of small particle packed columns varies from vendor to vendor.

The authors would like to acknowledge the contribution of Susan Serpa and Patty David