Waters HPLC 2004

COMPARISON OF DIFFERENT REVERSED-PHASE PACKING MATERIALS BASED ON BRIDGED HYBRID PARTICLES

Nicole Lawrence, Kevin Wyndham, John O'Gara, Ken Glose, Pamela Iraneta, Bonnie Alden, Cheryl Boissel, and Thomas Walter Waters Corporation, 34 Maple Street, Milford, MA 01757-3100, USA Contact: Nicole_Lawrence@waters.com

OVERVIEW

Since the launch of our first generation Methyl Hybrid particle (XTerra®)¹ we have continued to explore the use of hybrid organic/inorganic particles as reversed-phase HPLC packing materials. We recently developed the use of ethyl bridged alkoxysilanes as particle precursors.² By employing a 4:1 molar ratio of inorganic (SiO₂) to organic substituent groups (SiO_{1.5}CH₂CH₂SiO_{1.5}), novel silsesquioxane-based materials have been shown to be excellent base particles for the preparation of efficient and resilient reversed-phase (C_{18} , C_{8}) packing materials.



In this report we present the synthesis and characterization of different reversed-phase bondings on Bridged Hybrid particles. We will also evaluate the reversed-phase chromatographic performance of these new packing materials, and compare these results with Methyl Hybrid and conventional monofunctional C₁₈ and sterically protected C₁₈-bonded silica columns.

BRIDGED HYBRID BONDING CHEMISTRIES



- Surface derivatization of 5 µm materials by reaction with alkylchlorosilanes and endcapping. (SA= 189 m²/g; APD = 148 Å)
- Bridged Hybrid particles contain 6.4 %C before bonding
- Bridged Hybrid materials have 2.9-3.3 µmol/m² ligand coverage

CHROMATOGRAPHIC EVALUATION



- 23.4 ± 0.1°C.
- standard set of test analytes.³
- C₁₈ columns.
- (USP Tailing Factors = 0.9-1.3).

COLUMN STABILITY UNDER ACIDIC CONDITIONS

- bonded phase hydrolysis.
- different columns.²

Test	Proto
Challenge:	1%
Wash Step:	1%
Test:	1%
Repeat:	unt



Figure 1: Chromatographic comparison of different bondings on 5 µm Bridged Hybrid particles. Column: 4.6 x 150mm. Mobile Phase: 65 / 35 methanol / 20 mM KH2PO4 / K2HPO4, pH 7.00, 1.4 mL /min,

• Chromatographic comparison tests performed using standard conditions and a

• Bridged Hybrid columns have high resolution and unique selectivities.

• Bridged Hybrid columns show similar selectivity as XTerra[®] MS C₁₈ and Symmetry[®]

• Bridged Hybrid columns have excellent peak shape for all analytes

• Generally accepted failure mechanism for silica columns at low pH is acid-catalyzed

• Accelerated low pH aging test was developed to compare chemical stability of

- ocols (3.0 x 50 mm columns, 80°C):
- 5 TFA in water (0.6 mL/min, 240 min)
- TFA in Acetonitrile (1.3 mL/min, 90 min)
- TFA in water
- til a 50% loss in original ethyl paraben retention



Figure 2: Accelerated low pH aging of different 5 µm Bridged Hybrid columns and conventional monofunctional Silica C₁₈ columns.

- Acid stability increases with chain length and trifunctional bondings.
- Bridged Hybrid C_{18} bonded particles show < 10% loss in k at 113 h, which is comparable to the sterically protected Silica $Z-C_{18}$.

COLUMN STABILITY AT PH 10, 50°C

- Generally accepted failure mechanism for silica columns at high pH is basecatalyzed particle dissolution.
- · Accelerated alkaline pH aging test was developed to compare chemical stability of different columns.¹





Figure 3: Loss of original efficiency (5 sigma method) for acenaphthene on different 5 µm Hybrid columns and conventional monofunctional Silica C₁₈ columns. For unbonded Bridged Hybrid the efficiency loss of the void marker, uracil, was reported.

- Silica C_{18} -columns fail within 15-17 hours.
- Bridged Hybrid columns are still good after 135 hours (Unbonded included).
- Drastic increase in high pH stability attributed to chemical hydrolytic stability of the base hybrid particle.

COLUMN STABILITY AT pH 12.3 (0.02 N NaOH), 50°C



Figure 4: Loss of original efficiency for different 5 µm Bridged Hybrid columns at pH 12.3.

- Bridged Hybrid C₁₈ bonded particles last for 60 hours before 50% loss in efficiency.
- Bridged unbonded particles last for 12 hours under these vigorous conditions.

CONCLUSIONS

- Chromatographic evaluation of Bridged Hybrid columns show excellent peak shape, unique selectivity, and great efficiency.
- Bridged Hybrid columns show great performance at both high and low pH:
 - C_{18} bonded phase shows great retention stability with 1% TFA (80°C).
 - Bonded and unbonded Bridged Hybrid columns display excellent chemical stability at pH 10-12.3, 50°C.

REFERENCES

- 1) Am. Lab. 1999, 31, 36. 2) Anal. Chem. 2003, 24, 6781
- 3) LC-GC 1994, 12, 468.