

Claude R Mallet\*, Ziling Lu, Jeff R Mazzeo, Applied Technology, Waters Corporation, 34 Maple St, Milford, MA 01757-3696

## Abstract

Today, LC/MS/MS is without question the preferred platform for the analysis of organic molecules from complex matrices. Various fields such as pharmaceuticals, toxicology, environmental and clinical, have made the transition from the traditional one-dimensional system (i.e. LC/UV ... etc) to hyphenated technologies (LC/MS, LC/MS/MS ... etc). Most applications utilize the electrospray ionization (ESI) interface (about 80 % of published papers) versus atmospheric pressure chemical ionization (APCI). Several reasons can be presented to explain this trend. ESI is an interface that is relatively easy to use, exhibits low solvent consumption and can be used for large analytes (up to 100 KDa), wide polarity range and applicability to thermally labile compounds.

Unfortunately, ESI is prone to a phenomenon called "ion suppression"<sup>1,2,3</sup> Ion suppression can lead to adverse effects during quantitation (i.e. poor accuracy and precision). Previous papers have reported that ion suppression is a direct result of endogenous material present in biological samples. However, the measured result is the combination of several operating conditions and parameters of an SPE/LC/MS/MS system. Little has been done to effectively monitor and/or choose optimized conditions for the extraction, clean up, separation and analysis. This work presents a simple setup for quantification of ion suppression/enhancement. Several mobile phase additives, ion pairing agents and SPE extracts were measured against a standard reference. The results demonstrated that a mixed mode ion exchange protocol leads to minimal ion suppression/enhancement (< 10 %) for compounds that were investigated.

1- Miller-Stein C, Bonfiglio R, Olah TV, King RC, Am. Pharm. Rev., 2000. 54

2- Matuszewski BK, Constanzer ML, Chavez-Eng CM, Anal. Chem., 2003, 75: 675

3- Mallet, CR, Lu, Z, Mazzeo, JR, Neue, UD, Rapid Commun. Mass Spectrometry, 2002, 16: 805

\* Corresponding author: claude\_mallet@waters.com





Ion suppression or enhancement is a known phenomena that occurs with a mass spectrometer equipped with an electrospray interface (ESI). Several papers in the literature explain in detail the formation of ions with this type of source, but several parameters, effects or observations are still unclear. For example, is ESI concentration or mass flow dependent? Some papers suggest that ESI is concentration/mass flow combination. Most applications use the ESI source in combination with LC (mostly reversed-phase column). The problem of suppression or enhancement occurs when additives in the mobile phase can either increase or decrease the intensity of the target analyte. Sample extracts also produce the same effect.



# Can ion suppression be quantitated?

### SPE extracts



Ammonium bicarbonate

-30.8

<sup>3</sup> -43 €<sup>-9,8</sup> -12.7 € -17

-26.7 -32.6 C



• = 50 mM

# Salt & Ion pairing additives\*



monium biphosphate

A = Propranolol B = Trimethoprim C = Pipenzolate D = Resperidone E = Terfenadine F = Methoxy-verapamil G = Benextramine H = Reserpine 1 = Fumaric Acid 2 = Malic acid 3 = Etidronic Acid 4 = Clodronic acid 5 = Niflumic acid 6 = Canrenoic acid 7 = Cholic acid 8 = Raffinose







Rat plasma (un-centrifuged & unfiltered)

umes are not optimized for LC/



A = Propranolol B = Trimethoprime C = Pipenzolate D = Resperidone E = Terfenadine F = Methoxy-verapamil G = Benextramine  $\hat{H} = Reserptine$ 

# **Conclusions**

On average, up to 95 % of the ESI signal was suppressed when rat plasma samples were prepared using a protein precipitation protocol for both solutions of both basic and acidic test analytes. The same conclusion can be drawn for human plasma. However, when a mixed mode SPE method is used to prepare plasma extracts, the results shows an average of only 10 % suppression and even an enhancement effect for some compounds. This is certainly not a novel observation, and has been reported in the literature since the introduction of LC/MS/MS for such analytical applications. It certainly shows that in order to gain lower LOQ and better accuracy, an optimized extraction and clean up protocol becomes crucial.

Dilute with 0.5 mL H<sub>2</sub>O