

Amy E. Daly, Martin Gilar, Bonnie A. Alden, Ying-Qing Yu, John C. Gebler
Waters Corporation, 34 Maple Street, Milford, MA 01757, USA

Overview

- Proteolysis is thought to be a well understood process, however, many exceptions to the basic cleavage rule (cleave after K or R, but not when followed by P) have been reported (Table 1).
- In this study the selected proteins were digested and analyzed by LC-MS, LC-MS/MS, and MALDI-MS in order to evaluate the fidelity of trypsin (Figures 1-4).
- The examples demonstrating the variability of proteolytic cleavage are given (Figure 1). This variability makes it difficult to use protein digests as HPLC and MS standards and potentially complicates the proteome analysis of global protein digests.
- Ionic surfactant (RapiGest™ SF) was found to enhance the speed and completeness of the digestion (Figure 2) (Yu et al. *Anal Chem.* 2003 Nov 1, 75(21): 6023-8, Yu et al. *Rapid Commun Mass Spectrom.* 2004; 18(6): 711-5, Suder et al. *Rapid Commun Mass Spectrom.* 2004; 18(7): 822-4.)
- A set of standards was prepared by digesting selected proteins with trypsin using RapiGest™ SF for improved tryptic digestion.
- RapiGest™ SF decomposes to MS (LC/MS and MALDI-TOF) non-interfering by-products upon acidification
- Digestion using RapiGest™ SF results in fewer incompletely cleaved tryptic peptides. However, certain motifs in the primary amino acid sequence are still resistant to cleavage (Figure 3). Some of those motifs differ from those reported in the literature (Table 1).
- Standards were prepared to be free of salts, undigested protein, trypsin, and alkylation and reduction agents.
- Standards are useful for MS/HPLC instrument testing, calibration, MALDI-MS peptide mass fingerprinting and 2D-HPLC method development.

Experimental

Yeast alcohol dehydrogenase bovine serum albumin, bovine hemoglobin, and rabbit phosphorylase b were reduced and alkylated (DTT, iodoacetamide, Sigma) prior to digestion (sequencing grade trypsin, Promega), yeast enolase and cytochrome c were digested without reduction/alkylation. The acid-labile surfactant RapiGest™ SF was added to the digestion buffer (50 mM ammonium bicarbonate, pH 7.9). Progress of the digestion was monitored by LC-MS. The complete digestion of protein was achieved only in the presence of surfactant. Digested proteins were separated using BioSuite™ C₁₈ 3 µm PA-A, 150 x 2.1 mm columns with on-line MS detection (Waters Micromass® ZQ™ 4000 single quadrupole instrument). Mobile phases were A: 0.02% TFA in water and B: 0.016% TFA in acetonitrile. Gradient conditions were 0.8% B/minute with a flow rate of 0.2 ml/min. Samples were also analyzed by MALDI-MS using alpha-cyano-4-hydroxy-cinnamic acid as the matrix

HPLC System: Waters Alliance® 2795 Separations Module (Milford, MA) with a 996 PDA detector and online MS detection using a Waters Micromass® ZQ™ 4000

MALDI TOF MS: Waters MALDI® LT TOF MS, Reflectron mode

Figure 1: LC-MS analysis of cytochrome c tryptic digest. A—incomplete digestion, B—advanced digestion (extended digestion time and/or addition of trypsin). Some peptides partially resistant to cleavage (T9-10, T12-13, and T13-14) were digested in case B, however, the T4-5 peptide was highly resistant to tryptic cleavage.




Figure 2: Comparison of an enolase digest with and without RapiGest™ SF. Enolase was digested with sequencing grade trypsin at a ratio of 50:1 (w/w) protein to trypsin with and without 0.05% RapiGest™ SF present in the sample. It is clear that enolase is more completely digested in the presence of RapiGest™ SF.




Figure 3: LC-MS analysis of MassPREP™ protein digestion standards. Number of tryptic peptides detected and identified by their molecular weight is indicated. A complete and reproducible digestion of proteins was achieved with the aid of RapiGest™ SF acid labile surfactant. The most resistant missed cleavages are indicated with an asterisk.





Table 1: Motifs that are partially resistant to trypsin cleavage. The column on the left lists the partially resistant motifs that are listed in the Swiss-Prot database. The column on the right lists the resistant motifs identified in this study (Figure 1 and Figure 3 peaks labeled with an asterisk). Unique motifs are listed in red.

Swiss-Prot Partially Resistant Motifs	Resistant Motifs Found in this Study
CKY	EKD
DKD	GKE
CKH	RKV
CKD	DKD
KKR	DKL
RRH	DKG
RRR	VKV
CRK	QKC
DRD	
RRF	
KRR	

Figure 4: LS-MS/MS analysis of MassPREP™ enolase digestion standard. Peaks labeled with an asterisk result from secondary chymotryptic activity of Promega trypsin. Interestingly, only certain available sites were chymotryptically cleaved.



Conclusions

- The LC-MS monitoring of protein tryptic digestion suggests that proteolysis rarely proceeds to absolute completion. Therefore, it is difficult to reproducibly prepare protein digestion standards.
- The use of RapiGest™ SF acid-labile surfactant as an additive to the digestion buffer increased the speed and completeness of the digest.
- RapiGest™ SF improves digestion by denaturing/unfolding the proteins, making them more amenable to proteolysis and generating fewer incompletely cleaved peptides.
- The surfactant had a positive impact on the completeness of the digestion, however, highly resistant motifs were not cleaved.
- Several new resistant motifs were found in this work (Table 1).
- Our results suggest that RapiGest™ SF is useful for preparing more complete and more reproducible global protein digests.
- The five protein digests were prepared and qualified by LC-MS and MALDI-MS. They were found to be useful as samples for 2D-HPLC method development.