

Fang Xia*, Jie Y. Cavanaugh, Martin Gilar, Uwe D. Neue, Jeffrey R. Mazzeo, Diane M. Diehl
Waters Corporation, 34 Maple St., Milford, MA 01757-3696 (* Send Correspondence to Fang_Xia@waters.com)

OVERVIEW

Atlantis™ dC₁₈ columns are a fully LC/MS compatible line of reversed-phase (RP) columns designed for retaining and separating both polar and non-polar compounds. This new silica-based packing material was designed with the best combination of pore size, ligand density, ligand type, and is fully compatible with 100% aqueous mobile phases.

The overloading profiles of synthetic peptides on analytical columns are shown first. We also show the separation of a crude synthetic peptide mixture on both analytical and semi-preparative columns. Finally, a preparative load of peptide is achieved successfully on an analytical dimension column. The fraction collector was triggered by the mass-to-charge ratio of the target peptide to obtain a high purity (>99%).

These results indicate that Atlantis™ dC₁₈ RPLC columns provide good efficiency, high mass loading, and ease of scale-up for synthetic peptide crude samples.

INTRODUCTION

Synthetic peptides are too polar to retain on traditional RPLC columns, because traditional columns dewet under 100% aqueous conditions necessary for retaining polar analytes. The Atlantis™ columns are specially designed to meet the needs of chromatographers looking for polar compound retention.

The Waters AutoPurification™ system is capable of running preparative chromatography in an automated fashion. FractionLynx™ software installed with MassLynx™ triggers, tracks, controls and documents runs so the user can carry out purifications in an automated fashion and an unattended mode. This scalable MS- and UV-based system provides a simple but powerful automated purification process for successful compound discovery and development.

MASS OVERLOADING




Figure 1: Mass overloading of synthetic peptide crude samples on an analytical Atlantis™ dC₁₈ columns. Mobile phase A: 0.1% trifluoroacetic acid in water; Mobile phase B: 0.1% trifluoroacetic acid in acetonitrile.

- Atlantis™ dC₁₈ columns are capable of retaining **synthetic peptide crude** under highly aqueous conditions.
- Displacement profile** were observed when overloading the column with synthetic peptide, because the adjacent hydrophilic impurity were pushed forward by the enrichment of the target peptide.
- Displacement helps load even more samples on column without losing resolution.

SEPARATION OF A SYNTHETIC PEPTIDE




Figure 2: Separation of synthetic peptides on analytical and semi-preparative Atlantis™ dC₁₈ columns. Mobile phase A: 0.2% formic acid; in water; Mobile phase B: 0.2% formic acid in acetonitrile. Target peptide sequence information: NH₂-DRNFLRF-COOH.

Waters ZQ™

ESI⁺

Capillary (kV)	3.5	LM Resolution	15
Cone (V)	25	HM Resolution	15
Extractor	3.0	Ion Energy	0.3
RF Lens	0.3	Multiplier (V)	650
Source Temp (°C)	100	Scan Range (m/z)	320-1920
Desolvation Temp (°C)	250	Continuum	
Cone Gas Flow (L/Hr)	50	Scan Time (sec)	2.2
Desolvation Gas Flow (L/Hr)	260	Delay Time (sec)	0.1

- Atlantis™ dC₁₈ columns are capable of retaining **synthetic peptides** under highly aqueous conditions.
- Achieve linear scale-up from analytical to preparative chromatography.

MASS DIRECTED FRACTION COLLECTION




Figure 3: Mass triggered fraction collection and the re-injection result.

- Achieve **preparative** peptide loading on an **analytical** column.
- MS directed fraction collection results in high purity of the target samples.

CONCLUSIONS

- Atlantis™ dC₁₈ RPLC columns are useful tools for the purification of **synthetic peptide crude** under highly aqueous mobile phase conditions.
- Displacement profile is observed when overloading peptide, which helps load preparative amount of peptide on an analytical column.
- The mass-directed fraction collection purification systems ensure high purity of target peptide samples.
- Atlantis™ dC₁₈ preparative columns are available in various dimensions for ease of scale-up, and high mass loading.

