


Ziling Lu, Diane M. Diehl and Jeffrey R. Mazzeo, Waters Corporation, Milford, MA 01757, USA

OVERVIEW

A new Oasis® solid-phase extraction (SPE) sorbent has been developed. Oasis® WCX is a mixed-mode sorbent that combines reversed-phase and Weak Cation eXchange mechanisms of retention. This material is especially useful for quaternary amines and strong bases, as well as bases that are unstable at high pH levels. This sorbent can also be used for retaining polar and hydrophobic bases, a clear advantage over silica-based WCX-type sorbents. We have developed robust SPE protocols that are highly selective and sensitive for the clean-up of bases from biological matrices such as urine and rat plasma. We can achieve excellent SPE recoveries for quaternary amines, polar and hydrophobic bases.

INTRODUCTION

To address the need for new SPE sorbents in the pharmaceutical, environmental and life science fields, a new weak ion-exchange SPE material was developed. This sorbent is a new addition to the Oasis® family of polymeric SPE products. Oasis® MCX material is a mixed-mode, *strong* cation exchange SPE sorbent that is too retentive for strong bases such as quaternary amines. In these applications, both the sorbent and analyte remain charged and cannot be easily released from the sorbent. The new Oasis® WCX material is also a mixed-mode sorbent, but contains a *weak* cation exchange functionality whose charge can easily be controlled by the pH of the solution. This sorbent is ideal for the clean-up of quaternary amines from samples, as well as for all other types of basic analytes.



HPLC Conditions

Column: XTerra® MS C₁₈ 2.1 x 20 mm *IS*TM, 3.5 µmMobile Phase A: 10 mM NH₄HCO₃, pH 10Mobile Phase B: MeOH with 10 mM NH₄HCO₃

Flow Rate: 0.4 mL/min

Gradient: Time (min)	%A	%B
0.0	95	5
3.0	5	95
4.0	5	95
4.1	95	5
5.0	95	5

Injection Volume: 10 µL

Instrument: Waters 2777 Sample Manager and Waters 1525µ Binary

HPLC Pump

MS/MS Conditions

Waters Micromass® Quattro UltimaTM ESI+

Source Temp: 150 °C

Desolvation Temp: 350 °C


Cone Gas Flow: 50 L/Hr

Desolvation Gas Flow: 550 L/Hr

Collision Cell: 2.2e⁻³ bar (Ar gas)

MRM Transitions:	Cone (V)	CID (eV)
Valethamate m/z 306.1 → 218.9	35	20
Protriptyline m/z 264.0 → 191.1	60	25
Atenolol m/z 266.9 → 144.9	45	25

ANALYTES



SPE METHODS

Oasis® WCX µElution Plate

Condition: 200 µL MeOH


Equilibrate: 200 µL H₂OLoad: 150 µL urine or 1:1 diluted rat plasma, spiked with 10 pg/µL each analyte. To disrupt protein binding, add 2% H₃PO₄ (total sample volume) to protriptyline onlyWash 1: 200 µL 25 mM phosphate buffer in H₂O, pH 7

Wash 2: 200 µL MeOH


Elute: 50 µL (25 µL x 2) 2% FA in MeOH

Dilute: 100 µL 2% NH₄OH in H₂O

Inject: 10 µL



Representative LC/MS/MS data. Recoveries are the SPE recovery for rat plasma samples on the Oasis® WCX µElution plate. Excellent recoveries are obtained for all three classes of bases on this material.



Oasis® MCX µElution Plate

Condition: 200 µL MeOH


Equilibrate: 200 µL H₂OLoad: 150 µL 1:1 diluted rat plasma, spiked with 10 pg/µL each analyte. To disrupt protein binding, add 2% H₃PO₄ (total sample volume)

Wash 1: 200 µL 0.1 N HCl

Wash 2: 200 µL MeOH

Elute: 50 µL (25 µL x 2) 5% NH₄OH in MeOHDilute: 100 µL H₂O

Inject: 10 µL



Representative LC/MS/MS data. Recoveries are the SPE recovery for rat plasma samples on the Oasis® MCX µElution plate. Excellent recoveries are obtained for both atenolol (polar base) and protriptyline (hydrophobic base). However, valethamate (quaternary amine) is retained on the sorbent.

Utilizing the flow chart can help to determine which Oasis® sorbent to use in SPE method development.

Oasis® WCX versus a Silica-Based Weak Cation Exchanger 10-mg 96-well Plates

Condition: 500 µL MeOH

Equilibrate: 500 µL H₂O

Load: 0.25, 0.5, 1.0, 2.0 and 2.5 mL saline, spiked with 20 pg/µL of each analyte

Wash 1: 500 µL 25 mM phosphate buffer in H₂O, pH 7


Wash 2: 500 µL MeOH

Elute: 250 µL (125 µL x 2) 2% FA in MeOH

Dilute: (1) 250 µL 2% NH₄OH in H₂O for protriptyline and valethamate(2) 750 µL 2% NH₄OH in H₂O for atenolol

Inject: 10 µL

COMPARISON VS. SILICA-BASED MATERIAL



A comparison of the results for saline spiked with the three bases after SPE clean up with the Oasis® WCX and a commercially available silica-based WCX material is shown above. Increasing amounts of spiked saline were loaded onto the sorbents. Excellent recoveries were seen for all three analytes under all loading conditions on the Oasis® WCX. However, on the silica-based WCX, the polar analyte was not retained during the load step and 80% or less recoveries were observed for the quaternary amine and hydrophobic analytes.

CONCLUSIONS

The mixed-mode Oasis® WCX SPE sorbent provides a means of selective, fast and robust sample preparation with high recoveries for quaternary amine, polar, and hydrophobic basic analytes. Unlike silica-based SPE materials, there is no breakthrough for polar analytes. This material provides superior recoveries to eliminate sample breakthrough as seen on silica-based SPE materials.