Waters

Retention Mechanisms of Mixed-Mode Cation-Exchange SPE Sorbents

Michael S. Young, Uwe D. Neue and Kevin M. Jenkins, Waters Corporation, 34 Maple Street, Milford, MA 01757 USA michael_s_young@waters.com

Introduction

Mixed-mode SPE sorbents combine an effective reversed-phase chemistry with ionexchange sites on a single particle. Mixed-mode cation-exchange sorbents are often employed for the extraction of organic bases from complex matrices. For example, these sorbents are utilized for the isolation of basic drugs from urine or other biological fluids or for the extraction of basic pesticides from plant or animal tissue This type of analysis demonstrates one of the great advantages obtained using mixed mode cation-exchange; the ability to effectively recover basic drugs from sample matrices that have been acidified to release protein-bound residues. However, the benefits of mixed mode SPE extend well beyond the ability to recover drug residues under acidic conditions. For example, the highly polar aminotriazole cannot be effectively extracted using reversed-phase sorbents at any pH, but are effectively recovered using mixed-mode cation-exchange. Another compelling advantage of the mixed-mode sorbents is the ability to effectively retain both acids and bases in the same experiment.

In this presentation, we will then look at some applications that demonstrate each of three analytical categories that mixed-mode cation-exchange provides a clear advantage; 1. extraction of organic bases from acidified matrices, 2. extraction of highly polar bases from aqueous samples, and 3. multi-residue extraction of acids and bases from environmental samples. Finally, we will demonstrate benefit of mixedmode sorbent chemistries containing weak cation-exchange sites for reterntion of strong bases and quaternary ammonium compounds.

Retention Behavior



Figure 2. Mixed-Mode Cation-Exchange



Reversed-phase sorbents are the most commonly used SPE materials. Among the sorbents used for reversed-phase SPE are the alkyl bonded silicas (C8, C18 etc.) and polymeric materials such as polystyrene/divinylbenzene. Oasis HLB is a reversedphase sorbent produced by co-polymerization of hydrophobic (divinylbenzene) and hydrophilic (vinyl pyrrolidone) monomers. SPE devices are typically utilized for extraction of drugs, pesticides, pollutants, metabolites and other analytes of interest, from various matrices. Figure 1 (above left) shows theoretical retention of acids, bases and neutrals on of a reversed-phase sorbent. At low pH, acids and neutrals are retained, but the bases in cationic form, are not. At high pH, bases and neutrals are retained, but the acids in anionic form, are not. Now, consider Figure 2 (above right) that shows the theoretical retention of the compounds on a mixed-mode cationexchange sorbent. This diagram indicates that good retention of all classes of compounds is possible at low pH.

Example 1. Extraction of Organic Bases From Acidified Matrices

This type of analysis is appropriate for drug or drug metabolite residue analysis in biological fluids such as plasma, urine or saliva. Most often the compounds of interest are bases and the sample is acidified to release the residues from protein binding. Figure 1 shows that there is minimal reversed-phase retention of basic compounds at low pH. However, Figure 2 shows that mixed-mode retention mechanism is *favored* under such conditions. This example shows the determination of MDMA (ecstasy) and metabolites in urine. This example and examples 2 and 3 utilize Oasis MCX, a mixed-mode strong cationexchange sorbent.

SPE Protocol Oasis MCX Cartridae 1 cc, 30 ma

Condition 1mL methanol/ 1 mL water	
Load 0.5 mL sample	mm (
Wash #1 1 mL 0.1 N HCI	locks drugs or resin
Wash #2 1 mLmethanol	 removes acid and neutral interferences
Elute 1 mLmethanol (10%NH₄OH)	 releases drug as free bases

The LC-MS chromatogram below is from a urine sample obtained from an individual suspected of driving under the influence of illicit drugs (HMMA and MDA are metabolites of MDMA).



LC/MS System Waters[™] Alliance 2695 Separations Module Waters/Micromass ZQ[™] mass spectometer



Example 2. Extraction of Highly Polar Bases

Figure 3 shows that the retention mechanism of highly polar bases on mixed-mode cation-exchange sorbents is mostly cation-exchange (assume pKa of 8-9 for each This type of analyte is poorly retained on base). reversed-phase sorbents, but is retained, under acidic conditions, using the mixed-mode sorbent.



Amitrole (aminotriazole) is a common herbicide that is a highly polar base. This example shows the determination of amitrole at sub ppb levels in river water of high humic content.





Example 3. Multi-Residue Extraction of Acids and Bases

To effectively extract acidic and basic compounds using standard liquid extraction methods or traditional SPE methods, two separate extractions must be performed on the same sample. Consider the retention diagram presented in Figure 2. At low pH, mixed-mode cationexchange sorbents will retain acidic and neutral compounds by the reversed-phase mechanism and basic compounds by the mixed-mode mechanism. Therefore, the mixed-mode sorbent can be used to simultaneously extract acids, bases and neutrals in one SPE experiment. We have applied this technology to the GC-MS analysis of a test mixture containing 33 compounds (acids, bases and neutrals) listed in EPA method 8270C. Structures for a representative acid, base and neutral compound from this list are shown below.



RTX 5-MS, 30 m x 0.25 mm (ID) Helium carrier gas @ 1 mL/min 35°C for 4 minutes, 8°/min, to 300°



1. pyridine	14. 2-
2. picoline	15. trie
3. aniline	16. 2-
4. phenol	17. 3-
5. benzyl alcohol	18. dit
6. o-cresol	19. 1-
7,8. <i>m,p-</i> cresol	20. 2-
9. o-toluidine	21. tet
10. phentermine	22. 2-
11. chloroaniline	23. 4-
12. dichlorophenol	24. dip
13. phenylenediamine	25. ph

Mass Spectrometer: Waters/Micromass GCT™ with Aailent 6890 GC

26. aminobiphenyl

29. methapyrilene

28. nitroquinoline oxide

31. dimethylbenzidine

32. acetamidofluorene

33. dichlorobenzidine

30. dimethylaminoazobenzene

acids bases neutrals

27. dinoseb

- nethylnaphthalene :hlorophenol hitroaniline hitroaniline penzofuran -aminonaphthalene -aminonaphthalene trachlorophenol methyl-5-nitroanilin nitroaniline
- phenylamine nenacetin

A Novel Mixed-Mode SPE Sorbent for Quaternary Ammonium Compounds and Strong Bases

These types of analytes are well retained by the mixed-mode mechanism on a strong cation-exchange based sorbent (Oasis® MCX). However, elution must be accomplished at high ionic strength using highly concentrated salt solutions or strong acids. Such an eluant cannot be evaporated and reconstituted in mobile phase and the high ionic strength is not compatible with LC-MS analysis. However, consider the retention and elution behavior shown in Figure 4. Cationic species retained on the weak ionexchanger, Oasis[®] WCX, are easily eluted using a mildly acidic eluent.

Figure 4. Retention and Elution of Quats on Mixed-Mode Sorbents



An SPE protocol and the structure of the sorbent are shown below. Also shown are LC-MS analysis of river water samples for paraquat and diquat using the Oasis® WCX Protocol.



©2004 Waters Corporation