# Waters

Ziling Lu, Diane M. Diehl, Claude Mallet and Jeffrey R. Mazzeo Waters Corporation, Milford, MA01757

## ABSTRACT

of choice. However, many chemists feel that SPE method development is too time consuming or too difficult. To address these issues, we took a systematic approach to SPE method development and outlined all the steps including sample pretreatment, sample solubility, SPE sorbents, wash steps and solvents, elution solvents, calculation of recovery and method reproducibility, using acids, bases and neutrals. In this presentation, we outline all the important steps for SPE method development, highlighting the pitfalls and how to avoid them. Based on these data, we then outline a straightforward approach to SPE method development.

### INTRODUCTION

The Oasis<sup>®</sup> family of copolymeric SPE sorbents offers five chemistries to satisfy all sample clean up needs. Oasis® HLB is based on a hydrophilic-lipophilic balanced copolymer that offers superior retention of both polar and non-polar analytes. In addition, four mixed-mode chemistries offer additional retention and clean-up benefits when analytes are ionizable. Oasis<sup>®</sup> MCX material is a mixed-mode, strong cation exchange SPE sorbent. It is too retentive for strong bases and quaternary amines, both the sorbent and analyte remain charged and cannot be easily released from the sorbent. The Oasis® WCX material is also a mixed-mode sorbent, but contains a *weak* cation exchange functionality whose charge can easily be controlled by the pH of the solution. Oasis<sup>®</sup> MAX material is a mixed-mode, strong anion exchange SPE sorbent. It is too retentive for strong acids that contain highly ionic species, such as phosphate groups. The Oasis® WAX material is also a mixed-mode sorbent, but contains a weak anion exchange functionality whose charge can be controlled by the pH of the solution.

## CHOOSING THE SORBENT

These are the recommended starting methods for the various sorbents. Solid phase extraction (SPE) has become the sample preparation tool The key point in understanding what sorbent to use is determining the For all sorbents: nature of the analytes. The following flow chart outlines the sorbent Prepare sample selection process. Condition/Equilibrate with MeOH and Water Load Sample Solution



Whenever the analytes are ionizable (i.e. acids or bases), it is an procedure. Prior to SPE method development, we determined that valethamate If the analyte is bound to proteins in the matrix, the analyte may pass advantage to select a mixed-mode ion exchange sorbent. This type of Recovery of the Extraction Procedure (RE)<sup>1</sup> is as follows: does not bind to plasma proteins. This is done by spiking a known through the sorbent with the proteins. sorbent offers the best selectivity and the cleanest extracts. In the figure amount of valethamate into rat plasma, treating aliquots with or shown below, we can see the differences in signal-to-noise for rat without acid, and taking the samples through the SPE method. The The flow rate during this step is also crucial – make sure to observe plasma spiked with amitriptyline at 0.1 ng/mL taken through protein same results were obtained for both sets of samples. Therefore, we the drops and adjust the vacuum so discrete droplets (i.e. not a stream) precipitation, reversed-phase and mixed-mode ion exchange clean-up Additionally, in LC/MS/MS analyses, the matrix can cause the are flowing from the device. did not use acid in the pretreatment step. We utilized the suggested steps. Clearly, the samples prepared using Oasis<sup>®</sup> MCX are the starting SPE protocol for the Oasis<sup>®</sup> WCX µElution Plate: suppression or enhancement of the signal. Therefore, it is important to cleanest and offer the best sensitivity. Condition: 200 µL MeOH <u>Wash Steps</u> measure this effect.



# **Developing SPE Methods: Outlining a Systematic Approach**

Oasis<sup>®</sup> WCX:

Oasis<sup>®</sup> WAX:

Wash 1: 25 mM Phosphate

Wash 2: 100% MeOH

Elute: 2% FA in MeOH

Wash 1: 2% Formic Acid

Elute: 2% NH₄OH in MeOH

Wash 2: 100% MeOH

Buffer, pH 7

### STARTING SPE METHODS

Oasis<sup>®</sup> HLB: Wash: 5% MeOH in water Elute: 100% MeOH

Oasis<sup>®</sup> MCX: Wash 1: 2% Formic Acid Wash 2: 100% MeOH Elute: 2% NH₄OH in MeOH

Oasis<sup>®</sup> MAX: Wash 1: 2% NH₄OH Wash 2: 100% MeOH Elute: 2% Formic Acid in MeOH

# TIPS AND TECHNIQUES FOR TROUBLESHOOTING

Definition of Recovery

Because the term "recovery" can be confusing, it is important to state a simple definition for the overall recovery of the SPE extraction

$$\% RE = \frac{Re \, sponse_{Extracted \, Sample}}{Re \, sponse_{Post-Extracted \, Spiked \, Sample}} \times 100$$

Matrix Effect (ME)<sup>1</sup> is as follows:

% Matrix Effect = 
$$\frac{\text{Response}_{\text{Post-Extracted Spiked Sample}}}{\text{Response}_{\text{Non-extracted Sample}}} \times 100$$

### Plasma Sample Pre-Treatment

To optimize the contact time with the sorbent and to prevent the plugging of wells, we recommend the following steps:

- 1. Dilute 1:1 with water
- 2. Disrupt protein binding by adding  $H_3PO_4$  to a final concentration of 2% (use base if compound is acid labile)
  - Note: Steps 1 & 2 can be combined into one.
- 3. Add the internal standard the amount of organic solvent should be less than 10% of the final volume – this will help to prevent protein precipitation and well plugging
- 4. If possible, mix the solution before loading onto the sorbent.

Many times analytes will bind strongly to proteins in biological samples. If this binding is not disrupted, the analyte will pass through the SPE device and result in low SPE extraction recoveries. Therefore, it is important to determine whether or not binding is occurring. In this example, protriptyline binds strongly to rat plasma proteins. Only  $H_3PO_4$  is strong enough to disrupt this binding. Note that the addition of NH<sub>4</sub>OH does not disrupt the binding.

2ndSET_PlasmaDilute_HLB_neutral_final50%MeOH_4
Oasis <sup>®</sup> HLB, no additive
2ndSET_PlasmaDilute_HLB_2%NH4OH_final50%MeOH
Oasis <sup>®</sup> HLB, 2% NH <sub>4</sub> OH
PlasmaDilute_HLB_2%H3PO4_final50%MeOH_3
$Oasis^{\mathbb{R}} HLB$ , 2% $H_3PO_4$
2ndSET_PlacmaDilute_HLB_Poct_final50%McOH_2

Breakthrough in the Loading Step If a sample is in too high an organic solvent, for example, the supernatant from an extraction step, the sorbent may not retain the analyte and breakthrough may occur.

Equilibrate: 200 µL H<sub>2</sub>O The wash steps in a method help to remove the interfering sample Load: 100 µL Spiked rat plasma (diluted 1:1 with H<sub>2</sub>O, 1 ng/µL valethamate) components and are crucial in reducing matrix effects. A good Wash 1: 200 µL 25 mM Phosphate Buffer, pH 7 methodology for determining the best wash and elution steps is to run Wash 2: 200 µL MeOH a 2D Optimization. Let's walk through the Oasis<sup>®</sup> HLB 2D Elute: 50 µL 2% Formic Acid in MeOH Optimization. Solutions of base and acid are made containing 0 to SPE Extraction Recovery: 102% 90% methanol. The analyte is spiked onto the SPE device and then eluted with each of the solutions. These eluents are analyzed and the data plotted as in the next figure. CONCLUSIONS









From these plots, an optimized method can be developed:

Wash 1: Base with 5% MeOH (remove proteins to prevent clogging of wells; bases are in neutral state for more retention)

Wash 2: Base with 40% MeOH (removes hydrophilic bases and neutrals and all acids)

Wash 3: 100% Water (removes residual ammonium hydroxide) Elute: Acid with 70% MeOH (101% recovery from rat plasma)

# METHOD DEVELOPMENT EXAMPLE - OASIS® WCX

The analyte is valethamate, a quaternary amine. Following the sorbent selection flow chart, the Oasis® WCX material should be used.



- SPE Method development can be simplified by first correctly selecting the sorbent
- Once the sorbent is selected, following the outlined tips for the various SPE steps will result in a robust, reproducible SPE method.

## **REFERENCES**

<sup>1</sup>Matuszewski, B.K., Constanzer, M.L., Chavez-Eng, C.M. Anal. Chem. **2003**, *75*, 3019-3030.