# Waters

# **Deglycosylation and Sample Cleanup Method for Mass Spectrometry Analysis of N-linked Glycans**

Ying Qing Yu, Martin Gilar, Jennifer Kaska, and John C. Gebler

Life Sciences R&D, Waters Corporation

A rapid method for analysis of glycans is presented. This method describes deglycosylation of glycoproteins aided by an enzyme-friendly surfactant (RapiGest<sup>™</sup> SF), sample cleanup using a HILIC chromatography performed in a 96-well microtiter plate (MassPREP™ HILIC µElution Plate), and MALDI MS analysis of the resulting glycans.

lycosylation is one of the most important types of posttranslational modification (PTM) in proteins. Due to the high degree of heterogeneity, the characterization of glycans is a challenging task. Mass spectrometry (MS) is a primary tool for biopolymer analysis; however, the characterization of (native) glycans is complicated by the time-consuming sample preparation and their poor MS ionization efficiency. A typical sample preparation method for MS involves a chemical or enzymatic cleavage of glycans, followed by salts, surfactants, and protein residues re m oval. Purified native glycans can be directly analyzed by MALDI-TOF MS.

The efficient sample deglycosylation is a key requirement for a successful and sensitive glycan analysis. Nevertheless, the quantitative glycan release (e.g., using enzymes) is rarely achieved, since the glycosylated sites of the proteins are often obstructed by the protein secondaryand tert i a ry structure.

The goal of this work was to develop a rapid and efficient deglycosylation of N-linked glyc op roteins with a glycosidase (PNGase F) aided with the enzyme-friendly surfactant, RapiGest<sup>™</sup> SF. This was followed with a novel micro-scale hydrophilic-interaction chromatography (HILIC) solid-phase extraction (SPE) plate (Waters® MassPREPTM HILIC µElution Plate) for a rapid sample cleanup prior to MALDI MS analysis using highly purified MALDI matrix (Waters<sup>®</sup> MassPREP<sup>TM</sup> MALDI Matrix, DHB).

# **Experimental Conditions**

#### **Deglycosylation of N-linked proteins**

The glyc oproteins were solubilized in 0.1% (w/v) RapiGest<sup>™</sup> SF solution prepared in 50 mM NH4HCO3 buffer, pH 7.9. Protein samples (e.g., ovalbumin) were reduced with 10 mM DTT for 45 min at 56 °C and alkylated with 20 mM iodoacetamide in the dark for 1 h at room temperature. The enzyme PNGase F (2.5-5 units) was added, and the protein solutions we reincubated for 2 h at 37 °C.

# LC-MS analysis of the protein deglycosylation

The RP HPLC instrument (CapLC<sup>®</sup> XE, Waters) was equipped with a microbore RP-HPLC column (Waters<sup>®</sup> Atlantis<sup>®</sup> dC<sub>18</sub> column, 3.5  $\mu$ m, 1.0  $\times$  100 mm). The LC separation was hyphenated with a Waters<sup>®</sup> Micromass<sup>®</sup> Q-Tof micro<sup>TM</sup>. Mobile phase A was made of 0.1% formic acid in Milli - Q<sup>®</sup> water (Mllipore Corp., Billerica. Massachusetts). Mobile phase B was made of 0.1% formic acid in 100% acetonitrile. A linear gradient was run from 0 to 60% B in 30 min (2% B per min). Separation was carried out with 35  $\mu$ L/min flow rate; the column temperature was set at 40 °C.

# Glycan cleanup using a 96-well microelution HILIC SPE plate

The N-linked glycans released from glycoproteins we re extracted using the 96-well, MassPREPTM HILIC µElution Plate attached to a vacuum manifold. Use of this SPE device involves an initial wash and instruction for oligosacequilibration of the sample well(s), sample charides sample cleanloading, sample well washing to remove up using the MassPREP™ undesired products, and final elution of the HILIC  $\mu$ Elution plate is isolated glycans. Figure 1 shows the opti-

Wash 200 بلا of Milli-Q water
Condition 200 µL of 90% acetonitrile
Constitute sample 80-90% acetonitrile The final volume: 100–750 µL
Load Constituted sample
Wash 200 µL of 90% acetonitrile
Elute 25–50 µL of 25 mM ammonium citrate in 25% acetonitrile

Figure 1: A general illustrated.

mized MassPREPTM HILIC µElution Plate SPE protocol for both neutral and the sialylated glycans. The entire process requires less than 20 min. The HILIC plate performance was evaluated with maltoheptaose standard. Load, wash, and elution SPE fractions we re quantitatively analyzed by an HPLC system with evaporative light scattering detection (ELSD). The mass balance revealed no breakthrough in the load fraction. Most of the material eluted in the first 25 µL elution. Total mass balance was 90%. Recovery was estimated to be a p p roximately 70%.

#### MALDI-QTOF MSexperiments

Ultra pure MassPREPTM MALDI matrix, DHB (2,5-Dihydroxybenzoic acid) was used for MALDI-TOF analysis. The matrix was reconstituted in 500 µL of pure ethanol to a final concentration of 20 mg/mL. Purified glycan solutions we re mixed with DHB matrix in one to one ratio; one µL was placed onto a stainless steel MALDI target. Waters<sup>®</sup> Q Tof Ultima MALDI was used to determine the molecular weight of the released glycans and perform MS-MS experiments to characterize the structure of the glycans. The typical collision energy used here was 70 to 120 V.

#### Results

#### Protein deglycosylation under denaturing conditions

In earlier reports, we described the use of a mild and enzyme-friendly surfactant, RapiGest<sup>TM</sup> SF, for denaturation of proteins prior to proteolytic enzymatic digestions (1). It was found that this surfactant improves the speed and completeness of enzymatic proteolysis, most noticeably for globular and membrane proteins (2). Therefore, we investigated the use of RapiGest<sup>™</sup> SF in conjunction with PNGase F for the enzymatic release of N-linked glycans.

Figure 2 shows the extent of the deglycosylation reaction of



Figure 2: LC-MS spectra of deglycosylated ovalbumin are shown. (a) Ovalbumin was solubilzed without the use of denaturant and was not deglycosylated. (b) Ovalbumin was denatured using 0.1% OG and deglycosylated. (c) Ovalbumin was denatrued in 0.1% RapiGest<sup>™</sup> SF solution and deglycosylated. The MS scans were deconvoluted to the MW of the protein. Complete deglycosylated was observed after 2 h deglycosylation for the RapiGest<sup>™</sup> SF solubilized ovalbumin.

chicken ovalbumin solubilized in 0.1% RapiGest<sup>TM</sup> SF (Figure2c) digested with PNGase F for 2 h in 50 mM ammonium bicarbonate solution. The prog ress of deglycosylation is apparent in comparison to a control ovalbumin sample (Figure2a) with no enzyme added. The deglycosylation was also carried out with the addition of 0.1% of non-ionic surfactant, n-octyl- $\beta$ -gly $\infty$ pyranoside (OG) (Figure2b).

The LC–MS analysis of samples produced the ESI spectra featuring the multiply charged protein states, which we redeconvoluted using Waters® Micromass® Masslynx<sup>TM</sup> MaxEnt<sup>TM</sup> 1 software (deconvduted MS spectra are shown in the right panel in Figure 2). As expected, no signal corresponding to the MW of deglycosylated p rotein was found in the control sample (Figure 2a). Interestingly, no distinguishable deglycosylation was also observed in the OG-mediated deglycosylation (Figure 2b). Multiple peaks between 44–45 kDa re p resent the various N-linked glycoforms of ovalbumin. The reaction in the presence of RapiGest<sup>TM</sup> SF shows nearly complete deglycosylation; the protein mass shifted and a prominent peak was detected at approximately 43 Kda, which is consistent with the MW of the unmodified protein.

## Glycan sample cleanup prior to MALDI-MS analysis

The glycans we re extracted using the MassPREP<sup>TM</sup> HILIC µEl ution plate. In a HILIC mode, the hyd rophilic glycans are retained due to a partitioning separation mechanism between the organic mobile phase and a layer of water adsorbed on the surface of sorbent. Since the high concentration of organic solvent is necessary to ensure good retention of glycans, the samples we re first diluted with ACN to a final concentration of 80–90%. Some precipitation of glycans may







Figure 4: MALDI MS spectra of 5 pmol BSA tryptic peptides. (a) Control sample, SDS was not removed by any SPE method. (b) MassPREP™ HILIC µElution Plate was used to remove the SDS prior to MALDI TOF analysis.

occur if they are present at high concentrations. It is not recommended to centrifuge samples prior to loading to the HILIC µe l ution plate. After plate conditioning (sample cleanup section in experimental), glycan samples we re loaded by gravity (Figure 1).

#### MALDI-Q Tof MS-MS of glycans released from ovalbumin

The MALDI-Q TOF MS spectra of underivatized N-linked glycans released from 10 pmol Ovalbumin were obtained (Figure 3). MS–MS fragmentations of selected ions we re performed to validate the glycan structures. For example, collision induced dissociation of the complex glycan ion of mass to charge ratio of 1948.734 (M + Na) was shown (Figure 3). This ion is observed in the MS mode with low ion intensity, howe ver, enough fragmentation ions were produced in the MS–MS mode to determine its structure(GlycoSuite database, Proteome Systems, Ltd.).

## Surfactant removal using the HILIC µElution plate

The MassPREP<sup>TM</sup> HILIC µElution Plate facilitates the re m oval of impurities including the surfactants, such as RapiGest<sup>TM</sup> SF from the sample. It can be used for surfactant re m oval in general, for example SDS from peptides/glycopeptides. Fig u re4 shows the MALDI MS analysis of the bovine serum albumin (BSA) tryptic digest. No signal was observed for the sample contaminated with 0.1% SDS, while BSA tryptic peptide signals were observed in high abundance without any ion suppression caused by the presence SDS.

#### Conclusions

We have developed a method suitable for fast and robust analysis of glycans released from glyc o p roteins. The method utilizes an enzyme friendly surfactant (RapiGest<sup>TM</sup> SF) that was shown to greatly accelerate a deglycosylation reaction via glyc o p rotein denaturation, which makes the glycans more accessible to enzymatic cleavage. A complete deglycosylation of proteins was achieved after 2 h incubation with PNGase F. The MassPREP<sup>TM</sup> HILIC µElution plate was utilized to extract and desalt the glycans prior to their MS analysis using MassPREP<sup>TM</sup> MALDI Matrix, DHB. The SPE method is fast and requires minimum sample manipulation.

#### References

- (1) Y.Q. Yu, M. Gilar, P.J. Lee, E.S Bouvier, J.C. Gebler, Anal Chem. 75 (21), 6023-6028, 2003.
- (2) Y.Q. Yu, M. Glar, J.C. Gebler, Rapid Commun. Mass Spectrom. 18(6), 711-715, 2004.

Waters Corporation 34 Maple Street, Milford, MA 01757 tel. (508) 478-2000, fax (508) 478-1990 www.waters.com