

Waters

XBridge™ Prep Columns: Scalability and Loadability for Preparative Separations

Fang Xia, Jie Y. Cavanaugh, and Diane M. Diehl, Waters Corporation

BEH Technology™, the second generation of patented organic-inorganic hybrid particle technology (HPT), is the new benchmark for HPLC columns. Waters XBridge™ Prep columns reach a new level of maximum loadability and direct scalability.

XBridge™ columns were designed to be the most pH-stable phases commercially available, while still providing maximum efficiency, peak shape, and robustness. For method development consideration, we offer C₁₈, C₈, phenyl and RP₁₈ chemistries, available 2.5, 3.5, and 5 μm particle sizes and dimensions from ana-

lytical to prep. XBridge™ Prep columns are manufactured with the patent pending Optimum Bed Density (OBD™) design, which helps us to achieve direct scale-up from analytical to preparative columns, with the same efficiency and excellent column lifetimes.

Experimental Conditions

Scalability

Columns: XBridge™ C₁₈ 5 μm 4.6 × 100 mm; XBridge™ Prep C₁₈ 5 μm 19 × 100 mm

Mobile Phase A: 10 mM ammonium bicarbonate buffer at pH 10

Mobile Phase B: Acetonitrile/100 mM ammonium bicarbonate buffer, pH 10 (90/10)




Figure 1: Scale-up of a critical pair of antifungal drugs from analytical to preparative XBridge™ columns. (A) XBridge™ C₁₈, 5 μm 4.6 × 100 mm. (B) XBridge™ Prep C₁₈ 5 μm 19 × 100 mm. Analytes: (1) econazole, (2) miconazole.




Figure 2: Separation of five basic drugs on XBridge™ Prep column in high-pH mobile phase. Analytes in order of elution: labetalol, quinine, diltiazem, verapamil and amitriptyline.

Flow Rate: 1.06 mL/min (analytical); 18 mL/min (preparative)
 Gradient: 10-min linear from 5% to 95% B
 Injection Volume: 30 μ L (analytical); 510 μ L (preparative)
 Sample: Econazole and miconazole in DMSO (100 mg/mL each)
 Instrument: Waters AutoPurification™ System

Loadability

Columns: XBridge™ Prep C₁₈ 5 μ m 19 \times 50 mm
 Mobile Phase A: 0.1% diethylamine in water
 Mobile Phase B: 0.1% diethylamine in acetonitrile
 Flow Rate: 23.9 mL/min
 Gradient: 8-min linear from 5% to 95% B
 Injection Volume: 660 μ L
 Sample: Labetolol (50 mg/mL), quinine (50 mg/mL), diltiazem (50 mg/mL), verapamil (100 mg/mL) and amitriptyline (50 mg/mL) in DMSO
 Instrument: Waters AutoPurification™ System

Results

The retention and separation of two antifungal drugs on the analytical XBridge™ C₁₈ column is shown in Figure 1A. Under the total load of 6 mg, we observe very symmetric peaks. The mass load was proportionally scaled-up and run on the preparative XBridge™ Prep C₁₈ column, as shown in Figure 1B. Note the direct scale up, excellent peak shapes, and total mass load of 102 mg.

The separation and loadability of five basic analytes on XBridge™

Prep C₁₈ column under high pH mobile phase conditions is shown in Figure 2. We successfully loaded 198 mg of bases on a 19 \times 50 mm column without sacrificing peak shape.

Conclusions

XBridge™ Prep columns provide highly efficient separations, direct scale-up, and maximum loadability, crucial for isolation of critical mixture components.

© 2005 Waters Corporation. Waters, Beh Technology, OBD, and Autopurification are trademarks of Waters Corporation.

Waters Corporation

34 Maple Street, Milford, MA 01757
 tel. (508) 478-2000, fax (508) 478-1990
www.waters.com