IDENTIFICATION OF N-LINKED GLYCOSYLATION SITES USING MALDI TIME-OF-FLIGHT MASS SPECTROMETRY

Waters

INTRODUCTION

Glycopeptides are typically generated using proteolytic enzymes such as trypsin, to cleave at specific sites. The resulting glycopeptides tend to give weak ion signals upon mass spectrometry analysis (ESI or MALDI) due to their large molecular weight. The glycosylation site(s) identification using tandem mass spectrometry is hindered by the low ion intensity and complicated by the fragment ions generated from both the peptide back bone and the glycosidic cleavages from the sugar moiety.

We explored a method that uses unspecific enzyme, pronase, to generate smaller glycopeptides; the majority of resulting glycopeptides has an amino acid sequence between two to six. The glycopeptides were enriched and desalted using a micro scale hydrophilic interaction solid phase extraction device (HILIC SPE) prior to MALDI-Q-Tof MS analysis. The protonated glycopeptide ions were enhanced dramatically using ammonium citrate doped MALDI Matrix, 2, 5-Dihydroxybenzoic Acid (Fig. 1). Sodiated ions were also observed as minor ions. Collision induced dissociation was performed on both the protonated and sodiated ions. MS/MS fragmentation spectra revealed that the proton tends to retain on the peptide moiety while the sodium ion tends to retain on the sugar moiety (Fig. 2, 3). Characteristic fragment ion pairs were observed and used to identify glycosylation sites for both the protonated and the sodiated ions. Model proteins, horseradish peroxidase (HRP) and α 1- acid glycol proteins (AGP) were used as examples to illustrate how N-linked glycosylation sites can be identified with minimum sample manipulations and data interpretation.

A) contro B) 10 mM Ammonium Citrate Figure 1. MALDI-QTof MS analysis of glycopeptides generated from HRP. A) The DHB matrix was solubilized in 100% ethanol without any additives. B) DHB matrix contains 10mM ammonium tris citrate.

HILIC SPE Protocol

The glycopeptides were desalted using HILIC SPE (MassPREPTM, Waters) before MALDI analysis.



METHODS

Pronase Digestion

Glycoprotein (50 – 100 mg) was solubilized in the surfactant, RapiGest[™] SF (0.1%, w/v) solution with 50 mM ammonium bicarbonate. The sample was reduced using 10mM dithiothretol and alkylated using 15 mM iodoacetamide. The pronase to protein concentration is 1 to 50. Enzymatic digestion was performed at 37°C overnight. In the case of AGP, the sialic acid residues were removed by adding 1 unit of $\alpha_2(3,6,8,9)$ -neuraminidase after the pronase digestion and incubated for another 12 hours.

MALDI MS/MS

Micromass MALDI-QTof Ultima[™] (Waters) was equipped with a fixed nitrogen laser source. Argon gas was used as the collision gas. The typical collision energy used was in the range of 70 to 120 V depends on the size of the analyte ions. MALDI Matrix used was 2,5-dihydroxy benzoic acid, DHB (MassPREPTM Waters).

FindPept tool in the SwissProt database was used to identify the possible peptide sequences by submitting the peptide masses and searching against a protein with a known primary sequence.



RESULTS

MS/MS fragmentation of the **HRP** glycopeptides (protonated and sodiated) reveals characteristic ion pair



Figure 2. MS/MS fragmentation spectra of the same HRP glycopeptides. A) The precursor ion is protonated. B) The precursor ion is sodiated. Characteristic ion pairs were labeled.

Proposed fragmentation pathway for the protonated glycopeptides (A) and the sodiated glycopeptides (B)



Ying Qing Yu, Jennifer Fournier, Martin Gilar and John C. Gebler Life Sciences R&D, Waters Corporation, Milford, MA USA 01757



MS/MS fragmentation of the **AGP** glycopeptides (protonated and sodiated)



Figure 3. MALDI-Q-Tof MS/MS fragmentation on A) a protonated glycopeptide, B) a sodiated glycopeptide. The characteristic ion-pairs were labeled along with their identity. FindPept tool identifies the glycosylation site at amino acid residue number 72.