# Waters

## INTRODUCTION

Detecting and measuring modifications in protein structure that affect biological activity is required in protein chemistry research. For protein biopharmaceuticals, such changes must be monitored throughout process development, stability testing, and quality assurance because structural changes can affect efficacy and safety. While peptide mapping can be used, it is time-consuming and labor-intensive. It is attractive to consider using molecular weight measurement of the intact protein to determine the degree of modification. Such an approach would require measurement of exact mass in combination with spectral interpretation that preserves quantitative information. For this approach to be used, high quality, exact mass spectra are required. Of equal importance, the spectra must be analyzed using an algorithm that does not require prior knowledge of the sample composition and that preserves relative quantitative information about the protein components of the sample. A Tof mass spectrometer (LCT Premier) in conjunction with the MaxEnt1algorithm meets these basic criteria.

## **METHODS**

Instruments

Waters 2796 Alliance Liquid Chromatography System

Waters Micromass LCT Premier ESI-orthogonal Tof mass spectrometer

Waters Micrmass MassLynx 4.0 software for control and data acquisition

### Column:

MassPREP<sup>™</sup> on-line desalting cartridge, 2.0 x 10 mm

## Separation

Flow rate 250 µL/min 0.1% formic acid in water (Buffer A) 0.1% formic acid in acetonitrile (Buffer B) Gradient 10-60% Buffer B, 2 minutes Temperature 40 °C

## Samples:

Standard proteins, bovine RNase B and yeast enclase, were purchased from Sigma Chemical Company, St. Louis, MO

Samples were dissolved and diluted to 10 pmol/µL each in aqueous 0.1% formic acid.

 $10\mu$ L injections = 100 pmol of protein on column

## **Tuning Tips**

- Tune for maximum declustering. On the LCT Premier one adjusts the source pressure to 2.0-3.0 eO. At a set source pressure increase the voltage of Ion Tunnel 1 and Aperture 1 to decrease m/z peak widths caused by noncovalent clusters. Do not increase Ion Tunnel 1 and Aperture 1 voltages past the point when charge stripping begins. (Do not allow the charge-state envelope to distort to high m/z).
- Tune with a protein which has a MW as close a possible to the MW of the proteins to be analyzed.
- Acquire data over the widest possible m/z range to determine the chargestate distributions for the proteins of interest







charge states. Deconvolution applies only to data in the spectrum view window

## TO DOWNLOAD A COPY OF THIS POSTER VISIT WWW.WATERS.COM/POSTERS

## RAPID MEASUREMENT OF PROTEIN MODIFICATIONS BY COMBINING EXACT MASS MEASUREMENT AND ADVANCED SPECTRAL DECONVOLUTION

Figure 5. MaxEnt 1<sup>™</sup> parameter window.

Beth L. Gillece-Castro, Ziling Lu, Thomas E. Wheat, and Jeffrey R. Mazzeo Waters Corporation; 34 Maple Street; Milford, MA01757

- forms of a protein.