DETECTOR CONSIDERATIONS IN DEVELOPING PERFORMANCE QUALIFICATION (PQ) CRITERIA FOR HPLC SYSTEMS

Vaters

INTRODUCTION

Peptide mapping is a particularly important tool for the characterization of therapeutic proteins. When this technique is used in a validated analytical method, it is necessary to characterize the system and to apply appropriate qualification procedures. A performance qualification test based on peptide analysis is required for predicting system performance with real samples. Such qualification protocols, however, only apply to a single instrument configuration. It is common practice, however, to use different detectors for specific purposes at different stages of development. We have, therefore, established performance standards for reversed-phase HPLC in conjunction with both UV and MS detection for systems dedicated to running peptide maps. The qualified optical detectors include tunable UV and photodiode array. The mass spectrometers include single quadrupole, oa-Tof, and tandem quadrupole-oa-Tof. The criteria and specifications include retention time reproducibility and chromatographic selectivity. Mass accuracy was included for MS configurations. Since injector reproducibility and injector linearity can only be judged in conjunction with the selected detector, separate specification criteria were defined were each of the five detectors. The combinations of optical and mass detectors were also tested This qualification method is based on the separation of a well-characterized peptide analyte mixture on a tested HPLC column specific for peptide analysis. The consistent separation chemistry and test materials ensure that the test may be efficiently adapted to the selected detector. The developed procedure could also be used as a routine system suitability test.

MATERIALS AND METHODS

Test Samples

Waters MassPREP[™] Peptide Standard Mixture—dissolved in H₂O with 0.02% TFA and 10% ACN to a concentration of 6.7 pmoles/µL for use with ZQ and UV or to a concentration of 2 pmoles/µL for use with LCT, Q-Tof *micro™*, and UV Waters BioSuite[™] Peptide Standard–dissolved in H₂O with 0.02% TFA to a concentration of 10 pmoles/µL for use with ZQ and UV or to a concentration of 3 pmoles/ μ L for use with LCT, Q Tof μ , and UV.

Instrumentation

Waters AQUITY UPLC[™] System

Detectors

Waters ACQUITY UPLC[™] TUVe Waters ZQ Single Quadrupole Mass Spectrometer Waters LCT Premier oa-Tof Mass Spectrometer Waters Q-Tof *micro™* Hybrid Tandem Mass Spectrometer

Column

Waters AQUITY UPLC[™] BEH130 C18 1.7µm, 2.1 x 50mm Peptide Separation Technology

Mobile Phase

A: 0.02% TFA in water

B: 0.018% TFA in Acetonitrile

Gradient Methods

LC Method for MassPREP Pentide Mixture LC Method for BioSuite Pentide Standard

	Time	Flow Rate (mL/min)	% A	% B	Time	Flow Rate (mL/min)	% A	% B
	0.00	0.200	100.0	0.0	0.00	0.200	100.0	0.0
	10.00	0.200	50.0	50.0	3.00	0.200	80.0	20.0
	11.00	0.200	20.0	80.0	3.50	0.200	20.0	80.0
	11.10	0.200	100.0	0.0	4.00	0.200	100.0	0.0
	17.00	0.200	100.0	0.0	8.00	0.200	100.0	0.0





Table 1a and b: Absolute and Relative Retention

	Peak A	Peak B	Peak C	Peak D	
	(4.3+/- 0.5)	(6.6+/- 1.0)	(8.0+/- 1.0)	(10.3+/- 1.0)	
lnj #1	4.26	6.60	7.97	10.27	lnj #1
ln j#2	4.26	6.60	7.97	10.28	ln j#2
Inj #3	4.29	6.60	7.97	10.28	Inj #3
lnj #4	4.24	6.60	7.97	10.29	lnj #4
lnj #5	4.26	6.61	7.97	10.30	lnj #5
Inj #6	4.22	6.61	7.97	10.30	Inj #6
Mean	4.26	6.60	7.97	10.29	Mean
Stdev	0.0235	0.0052	0.0000	0.0121	
% RSD	0.551	0.078	0.000	0.118	

Figure 2: Chromatograms with Alternate Detectors



Qualification of Chromatographic Elements: The chromatogram in Figure 1 shows the separation of the MassPrep[™] Peptide Standard Mixture with the defined method. The qualification criteria include the absolute retention times and the repeatability of those retention times over six replicate injections. The specification range shown in the table could be narrower. A second sensitive test of all the factors that affect chromatography is the relative retentions shown in Table 1b. The same test can be used for any of the detectors or combination of detectors that might be used in the system, as shown in Figure

TO DOWNLOAD A COPY OF THIS POSTER VISIT WWW.WATERS.COM/POSTERS



<u>RTB RTC RTD</u> RTA RTB RTC 1.21 1.21 1.21 1.29 1.21 1.21 1.21

ACQUITY TUVe 214nm



LCT Premier TIC





Qualification of System Quantitative Performance: The observed peak areas in a chromatogram are affected by chromatographic performance, by the response of the detector, and by the integration algorithms. Successful tests of the reproducibility of peak areas at constant injection and linearity of response are only possible if all system components are performing properly. Here the same test is used for MS and for UV detection.



Qualification of MS and MS/MS Functions — The use of a mass spectrometer as the information rich detector in a peptide mapping system expands the range of required qualification tests. As shown above, the same tests of chromatographic and quantitative elements can be used with both UV detectors and mass spectrometers. Further tests are required to ensure that the MS system produces qualitative data suitable for the analytical problem. The tests of mass resolution and mass accuracy are obvious. The data is collected in LC_MS mode, rather than by infusion, to most closely reflect the performance of the detector in a peptide mapping experiment. The same test sample and UPLC[™] method can be used for all the detectors. However, different dilutions are required to reflect the different instrument sensitivities. The target specifications are also different as required by the different mass analyzers. When the system includes a Q-Tof *micro™*, MS/MS function must be tested. The same peptide test can be used here, specifying the known mass and observing the product ions expected for that sequence. The process is simplified because the same sample and chromatographic method can be used for all the mass analyzers. Only differences in expected performance need be defined.

303.16 371.17

728.33

710.32

775.35

Ziling Lu, Thomas E. Wheat, Beth L. Gillece-Castro, and Jeffrey R. Mazzeo Waters Corporation, Milford, MA01757 USA

QUALIFICATION REQUIREMENTS

Retention Time Reproducibility – Qualifies chromatographic elements, i.e., Solvent Delivery System, Temperature Control, Column, and Mobile Phase.

Quantitative Reproducibility—Qualifies combination of chromatographic elements, injector, and detector for precision.

Quantitative Linearity—Qualifies combination of chromatographic elements, injector, and detector for usable ranges.

MS Functions—Qualify mass accuracy and resolution for impact on peak identification, as well as, effect on quantitative properties.

MS/MS Function—Qualify selection of peptide in an LC/MS chromatogram and the generation of the correct product ions in the collision cell.

PQ Tests for System Configurations					
System Configuration	PQ Tests to Run				
UV-based system	Injector Reproducibility and Linearity Tests				
	Preliminary Peptide Mapping Run				
LC_UV	System Reproducibility Test				
UV/ZQ-based system	Injector Reproducibility and Linearity Tests				
or ZQ-based system	Preliminary Peptide Mapping Run				
	System Reproducibility Test				
LC_UV_MS or LC_MS	Mass Accuracy				
UV/LCT-based system	Injector Reproducibility and Linearity Tests				
or LCT-based system	Preliminary Peptide Mapping Run				
	System Reproducibility Test				
LC_UV_MS or LC_MS	Mass Accuracy				
	Mass Resolution				
	Mass Sensitivity				
UV/QTof-based system	Injector Reproducibility and Linearity Tests				
or QTof-based system	Preliminary Peptide Mapping Run				
	System Reproducibility Test				
LC_UV_MS/MS or LC_MS/MS	Peptide Identification				
	Mass Accuracy				
	Mass Resolution				
	Mass Sensitivity				

SUMMARY

- Procedures have been developed for System-level performance qualification of ACQUITY UPLC[™] Peptide Mapping Systems
- The qualification protocols are based on the analysis of peptides under conditions similar to those used for routine peptide mapping
- The same test samples and chromatographic methods can be used for the ACQUITY UPLC[™] system with any detector or combination of detectors that would be routinely used in a peptide mapping system
- The only required operating adjustments are in sample concentration for detectors of widely disparate sensitivity
- Suitable specifications can be readily defined for all the instrument combinations using actual peptide data