# Vaters

## INTRODUCTION

Peptide maps that are used in the characterization of biopharmaceuticals must completely resolve all possible peptides derived from the sample, including those representing a variety of minor chemical modifications. Development of such separations is a time-consuming and often labor-intensive procedure. It is greatly influenced by the experience of the scientists involved. In making this procedure more efficient, various columns are evaluated. Those columns are selected because they have properties that are expected to interact with peptides in useful ways. These variables include column dimensions and particle size, but the greatest significance is often ascribed to bonded phase chain length, pore size, and base material. The relative importance of these properties is difficult to estimate since there are few examples that compare columns that differ in only one parameter. Such comparisons are developed in these experiments. Where necessary, small amounts of prototype packing materials were synthesized to permit evaluation of each relevant property.

### **METHODS**

| LC System:             | Waters 2796 Separation Module                                  |  |  |  |
|------------------------|--|--|--|--|
| UV Detection:          | Waters 2487 Dual Wavelength Absorbance                         |  |  |  |
|                        | Detector. Wavelength 214 nm                                    |  |  |  |
| MS System:             | Waters Micromass ZQ™ Mass Spectrometer                         |  |  |  |
|                        | Electrospray Ionization (+)                                    |  |  |  |
| Mobile Phase:          | With "TFA" modifier:   |  |  |  |
|                        | A = 0.02% Trifluoroacetic Acid in Water                        |  |  |  |
|                        | B = 0.018% Trifluoroacetic Acid in Acetonitrile                |  |  |  |
|                        | With "FA" modifier   |  |  |  |
|                        | A = 0.1% Formic Acid in Water                                  |  |  |  |
|                        | B = 0.1% Formic Acid in Acetonitrile                           |  |  |  |
| Flow Rate:             | 0.3 mL/min   |  |  |  |
| Injection Volume       | : 20 μL  |  |  |  |
| Columns: Wate          | rs   |  |  |  |
| <b>A</b> - BioSuite™ C | 18 PA-A, 2.1 x 150 mm 3.0 µm particles, 120 Å pores            |  |  |  |
| <b>B</b> - BioSuite™ C | <sup>18</sup> PA-B, 2.1 x 150 mm 3.5 µm particles, 300 Å pores |  |  |  |

- **C**-Symmetry<sup>TM</sup>  $C_{18}$ , 2.1 x 150 mm 5 µm particles, 300 A pores
- C- Prototype BEH<sup>™</sup> (Bridged-Ethyl-Hybrid)C<sub>18</sub>, 2.1 x 150 mm, 4.5 µm, 300 Å
- D- Prototype BEH™ (Bridged-Ethyl-Hybrid)C<sub>18</sub>, 2.1 x 150 mm 3.5 µm, 130 Å

#### RESULTS

### Effect of Peptide Properties on Retention

In an effort to better define the basis for column selection for peptide mapping, the MassPREP™ Enclase Digestion Standard was used as a model system. The peptides in this mixture were categorized based on those chemical properties most likely to affect retention and selectivity. These include relative hydrophobicity and hydrophilicity based on common models, calculated isoelectric point, pl, and size. These categories are tabulated below. Note that a peptide may appear in more than one category, sometimes with unexpected results. For example, the largest peptide is the most hydrophobic, but the second largest is among the most acidic.

| Size<br>(mass)                        | Hydrophobic | Hydrophilic | Basic<br>(pl = 10.1) | Acidic<br>(pl = 3.6 to 3.9) |
|---------------------------------------|-------------|-------------|----------------------|-----------------------------|
| T21 (3737)                            | T21         | T3          | T5                   | T27                         |
| T27 (3257)                            | T35         | T5          | T16                  | T45                         |
| T35 (1872)                            | T16         |             | T50                  | T14                         |
| , , , , , , , , , , , , , , , , , , , |             |             |                      | T37                         |





drophilic peptides near T50 show a useful difference.





Figure 3: Effect of Gradient Slope on Selectivity: More shallow gradients are usually used to increase resolution. In this example of the map of the tryptic digest of enolase, rearrangements in selectivity often occur with such changes Two peaks that are separated in the steep gradient merge in the shallow while two other peaks reverse elution order.

#### TO DOWNLOAD A COPY OF THIS POSTER VISIT WWW.WATERS.COM/POSTERS

## SYSTEMATIC DEVELOPMENT OF PEPTIDE MAPS FOR PROTEIN CHARACTERIZATION

terns observed with TFA as a modifier. The other three columns do differ from one another in ways that may be of benefit for different samples.

In these examples, the peptides with the greater change in retention are the species with the greater excess of ion pairing sites over ion suppression sites. The same mechanism can explain the highlighted change in selectivity between formic acid and TFA. Small changes in modifier can be used to improve the selectivity of the map and may be done in a planned way where the structure is known.

Ziling Lu, Thomas E. Wheat, Beth L. Gillece-Castro, and Jeffrey R. Mazzeo Waters Corporation; Milford, MA01757 USA

- Test both steeper and shallower slopes
  - Shallower than 0.25%/col. vol. generally degrades sensitivity
- Steeper than 2-3%/col.vol. noticeably compromises resolution
- Final adjustments
- Revisit temperature and modifier concentration as needed
- Consider dividing into multiple shorter methods for specific modifications