

OVERVIEW

Reversed-phase LC columns designed for polar compound retention often have low ligand densities and/or novel ligands in order to maintain pore wetting and enhance retention. Potential weaknesses of these types of columns include shortened lifetimes at low pH (e.g., TFA-containing mobile phases), poor peak shapes above pH 5 and MS bleed. To solve these issues, we developed an improved C18 column for polar compound retention (Atlantis® T3) using a trifunctional ligand and a propriety endcapping processes. We have seen dramatic improvements in low pH lifetimes (up to a 5-fold improvement) and high pH lifetimes (up to 2.5-fold improvement). Minimal dewetting is observed when the flow is stopped, then restarted, when using 100% aqueous mobile phases. MS response of column bleed was also measured and compared to other RPLC columns. In addition, these new columns exhibited similar selectivity to Atlantis® dC₁₈ columns for easy method transfer. Lastly, we show additional applications that are now possible with this new column chemistry.

LOW pH STABILITY TEST




Figure 1: Comparison of efficiency loss for C18 stationary phases during exposure to **0.5% TFA** mobile phases at 60°C.

Trammell, B. C. et al J. Chromatogr. A, 1060 (2004) 153-163.

HIGH pH STABILITY TEST




Figure 2: Comparison of stability of C18 stationary phases during exposure to **50 mM TEA pH 10** mobile phases at 50°C. (2a): plate count; (2b): capacity factor
Analyte: acenaphthene
Test mobile phase: 65:35 (MeOH:20 mM K₂HPO₄, pH 7) at 50°C
Wyndham, K.D. et al Anal. Chem. 2003 (75) 6781-6788.

DEWETTING




Figure 3: Results of the dewetting study on Atlantis® T3 and a conventional C18 column.

Walter, T. H. et al J. Chromatogr. A 1075 (2005) 177-183

MASS SPECTROMETRY

HPLC Conditions:
Mobile Phase A: 0.1% HCOOH in water
Mobile Phase B: 0.1% HCOOH in acetonitrile
Flow Rate: 1.5 mL/min with 250 µL/min flow split to mass spectrometer
Gradient: 10 min gradient from 5% to 100% B with 3 min initial hold

Column Temperature: 50°C
Column Dimension: 4.6 x 50 mm
Instrument: Waters Alliance® 2695

MS Conditions:
Ion Source: ES+
Capillary: 3.5 kV
Extractor: 3.0 V
Source Temp: 150 °C
Cone Gas: 60 L/hr
Ion Energy: 0.5
Cone: 33.0 V
RF Lens: 0.3 V
Desolvation Temp: 300 °C
Desolvation Gas: 509 L/hr
Instrument: Waters Micromass® ZQ™




Figure 4: Comparison of column bleed in mass spectrometry.

SELECTIVITY COMPARISON




Figure 5: Comparison of column selectivity.
Column: 4.6 x 150 mm, 5 µm
Flow Rate: 1.0 mL/min
Analytes: (1) thiourea, (2) 5-fluorocytosine, (3)adenine, (4) guanosine-5'-monophosphate, (5) thymine
Test mobile phase: 10 mM NH₄COOH, pH 3 (v/v) at 30°C

APPLICATIONS

Water Soluble Vitamins




Figure 6: Separation of water soluble vitamins

Column: Atlantis® T3 4.6 x 150 mm, 5 µm
Mobile Phase A: 0.1% CF₃COOH in water
Mobile Phase B: Acetonitrile
Flow Rate: 1.4 mL/min
Gradient: Time %A %B
0 100 0
4 97 3
6 85 15
15 80 20
24 80 20

Injection Volume: 10.0 µL
Sample: As indicated in the chromatogram
Detection: UV @ 260 nm

Catecholamines and Metabolites

Compound	USP Tailing	Sample Concentration
1. Norepinephrine	1.37	25 (µg/mL)
2. Epinephrine	1.34	25 (µg/mL)
3. Dopamine	1.29	10 (µg/mL)
4. 3,4- Dihydroxyphenylacetic acid	1.08	25 (µg/mL)
5. Serotonin (5-HT)	1.49	30 (µg/mL)
6. 5-Hydroxy-3-indoleacetic acid	1.18	25 (µg/mL)
7. 4-Hydroxy-3-methoxyphenyl acetic acid (HVA)	0.9	25 (µg/mL)




Figure 7: Separation of catecholamines and metabolites

Column: Atlantis® T3 4.6 x 150 mm, 5 µm
Mobile Phase A: 10 mM NH₄COOCH₃, pH 5.0
Mobile Phase B: Acetonitrile
Flow Rate: 1.0 mL/min
Isocratic: 98/2 (A/B)
Injection Volume: 10.0 µL
Sample: As indicated in the chromatogram above
Column Temp: 30 °C
Detection: UV @ 280 nm
Instrument: Waters Alliance® 2695 with 2996 PDA

CONCLUSIONS

- Atlantis® T3 columns provide:
- Enhanced retention for polar analytes
 - Superior stability under low and high pH conditions
 - Minimal dewetting when the flow is stopped, then restarted, when using 100% aqueous mobile phases
 - Low MS bleed
 - Comparable selectivity to Atlantis® dC₁₈ columns
 - Good peak shapes
 - Available in a variety of dimensions and particle sizes