DEVELOPING A SYSTEMATIC APPROACH TO METHOD TRANSFER FROM UPLC[™] TO PREPARATIVE HPLC

Waters

Fang Xia, Jie Y. Cavanaugh, Thomas E. Wheat, Diane M. Diehl and Jeffrey R. Mazzeo Waters Corporation, Milford, MA 01757, USA

OVERVIEW

Ultra Performance LC (UPLC[™]) has been widely accepted by chromatographers because of the improvements over HPLC in the sensitivity, resolution and speed of separations. As scientists begin to use UPLC[™] for impurity and metabolite profiling, they will need to transfer the methods to preparative LC to isolate and purify their compounds for further research. Therefore, it is necessary to systematically transfer UPLC[™] assays not only to HPLC, but, more importantly, to preparative chromatography. In this presentation, we discuss a systematic approach with key guidelines and observations to expedite the method transfer from UPLCTM to preparative LC. Parameters that need to be considered include column chemistry and dimension, injection volume, temperature, mobile phase composition, flow rate, gradient method, cycle time, detector sensitivity, injector design and system volume. Step by step calculations are included for a better understanding of evaluating method transfer validity. In addition, observations on reducing total run time in preparative LC will be covered.



GOALS

- Scale-up the separation from a 2.1×50 mm, 1.7 μm UPLC[™] column to a 19 × 150 mm, 5 μm preparative column
- Scale-down to a 19×50 mm, 5 µm column
- Maximize the loading on the 19×50 mm, 5 µm preparative column

SCALING FLOW RATES

ACQUITY UPLC™: AutoPurification[™]:

1.0 mL/min30.6 mL/min

- 1.0 mL/min flow rate in UPLC[™] was used for the fastest scouting separation. Geometrically scaled flow rate would be 81.9 mL/min on the prep column! But: $F_{\text{opt}} \propto \frac{1}{dp} \quad \text{Therefore:} \quad \frac{F_{\text{opt Prep}}}{F_{\text{opt UPLC}}} \propto \frac{dp_{\text{UPLC}}}{dp_{\text{Prep}}} = 27.8 \text{ mL/min}$
- Flow rates of 27.8 ±10% mL/min were run and 30.6 mL/min resulted in the best separation.

SYSTEM VOLUME COMPENSATION

2.60 mL

- ACQUITY UPLC™: 0.105 mL AutoPurification[™]:
- Differences in system volumes can result in differences in retention times, selectivity and gradient formation.
- An initial hold time may be added to the prep run to compensate for this difference.

SCALING FROM UPLC[™] TO PREP



Injection volume is scaled according to the column volumes

Target injection volume =

Target Column Volume Original injection volume X **Original Column Volume**



Ease of Migration from HPLC to UPLC[™]

ACQUITY UPLC[™] columns and XBridge[™] columns

Same column chemistry ensures the ease of scale-up

The initial hold time is calculated using the following equation:



GRADIENT SCALING

To maintain the same selectivity, segment durations in column volumes (c.v.) need to remain the same when scaling the gradient table from UPLC[™] to Prep LC.

Table 1: Gradient Tables for Various Dimension **Preparative Columns**

UPLC 2.1×50 mm, 1.7 μm	Steps	Time (min)	Flow Rate (mL/min)	%A	%В	Segment Duration Time (min)	Segment Duration (c.v.)
	Initial hold	0	0.6	95	5	0	0
		2	0.6	5	95	2	11.55
		3	0.6	5	95	1	5.78
		3.1	0.6	95	5	0.1	0.58
		5.0	0.6	95	5	1.9	13.87
UPLC 2.1x50 mm, 1.7 μm	Initial hold	0 2 3 3.1 5.0	0.6 0.6 0.6 0.6 0.6	95 5 5 95 95	5 95 95 5 5	0 2 1 0.1 1.9	0 11.55 5.78 0.58 13.87

Prep LC 19×150 mm, 5 μm	Steps	Time (min)	Flow Rate (mL/min)	%A	%В	Segment Duration Time (min)	Segment Duration (c.v.)
	Initial hold	0.78	30.6	95	5	0	0
		16.80	30.6	5	95	16.05	11.55
		24.83	30.6	5	95	8.03	5.78
		25.64	30.6	95	5	0.80	0.58
		44.90	30.6	95	5	19.26	13.87

Prep LC 19x50 mm, 5 μm	Steps	Time (min)	Flow Rate (mL/min)	%A	%В	Segment Duration Time (min)	Segment Duration (c.v.)
	Initial hold	0.20	30.6	95	5	0	0
		5.55	30.6	5	95	5.35	11.55
		8.22	30.6	5	95	2.68	5.78
		8.49	30.6	95	5	0.27	0.58
		14.91	30.6	95	5	6.42	13.87

Excellent scaling is observed from UPLC[™] to Prep

MAXIMIZE LOADING

XBridge[™] Prep OBD[™] C₁₈ 19 x 50 mm, 5 µm



Sample Concentration: 190 mg/mL total in DMSO Injection Volume: 1 mL Temperature: Ambient

Detection: UV @ 245 nm

Waters AutoPurification[™] system Instrument:

CONCLUSIONS

Guidelines for successful method transfer

- Scale and optimize the flow rates
- Calculate initial hold time for Prep based on the two system volumes
- Maintain the same segment duration (in column volumes) for geometric gradient scaling
- Scale the injection volume according to the column volumes
- Maximize loading on the prep column
- By following these simple guidelines, scaling from UPLC[™] to Prep LC is readily achieved

TO DOWNLOAD A COPY OF THIS POSTER VISIT WWW.WATERS.COM/POSTERS

UPLCTM SEPARATION

are both based on BEH Technology™.

from UPLC[™] to preparative LC.

