

IMPROVING GLYCOPEPTIDE FINGERPRINTS WITH HIGH RESOLUTION CHROMATOGRAPHY

Beth L. Gillece-Casto, Thomas E. Wheat, Ziling Lu, and Jeffrey R. Mazzeo
Waters Corp., Milford MA

INTRODUCTION

Characterization of glycoproteins entails analysis of the primary structure of the protein and analysis of the attached oligosaccharides. Oligosaccharides are routinely characterized by glycan mapping after release from the protein. Peptide maps can provide information about the site-specific micro heterogeneity of the glycans if the glycopeptides are resolved. The enhanced chromatographic resolution associated with Ultra Performance Liquid Chromatography™ has been demonstrated for peptide mapping.¹ The technique improves resolution by a factor of three or more. The Peptide Separation Technology Columns for UPLC™ include both 130Å and 300Å pore size materials. Coupling UPLC™ to oa-ToF mass spectrometry allows high sensitivity identification of glycopeptides through exact mass measurement and deconvolution of complex mass spectra. Optimized chromatographic resolution of individual glycoforms of glycopeptides yields high quality mass spectra. This spectral information can be used to develop UPLC™ methods for quantitative glycopeptide mapping with UV detection.

METHODS

UPLC Conditions

Injection Volume: 10 µL of 3-5 pmol/µL digest solution
Columns: Peptide Separation Technology
ACQUITY UPLC™ BEH 130 C₁₈ 1.7 µm, 2.1 x 100mm
ACQUITY UPLC™ BEH 300 C₁₈ 1.7 µm, 2.1 x 100mm
X-Bridge™ BEH C₁₈ 3.5µm, 2.1 x 150mm
Vydac 238MS™ 5µm, 2.1 x 250 mm
Temperature: 40°C (unless otherwise noted)
Flow Rate: 100 µL/min (unless otherwise noted)
Solvent A: 0.1% TFA in water
Solvent B: 0.08% TFA in acetonitrile
Mixer: High Sensitivity Peptide Analysis mixer
[P/N 205000403]
Gradient Table:

Time (min)	%A	%B
Init	100	0
2	100	0
118	50	50
120	25	75

MS Conditions

LCT Premier™ oa-ToF mass spectrometer
Capillary: 3000V
Cone: 35V
Source Temperature: 100°C
Source Pressure: 1.4 x e⁰ Torr
Scan Range: 300-3000 m/z
Scan Rate: 2 scans/sec

RESULTS

Glycopeptides with Neutral N-linked Glycans




Figure 1. Small portion of mouse IgG tryptic digest UPLC™ chromatograms. The glycoforms have been significantly resolved.




Figure 2. Mouse IgG glycoforms were also detected with MS, and extracted ion chromatograms of the triply-charged glycopeptides are shown. Four peaks were observed for the molecular weights indicating




Figure 3




Figure 3. Mouse IgG glycoforms were detected with UV at 214 nm. The peptide sequence is shown above. Area was integrated for the four glycoforms shown.

Figure 4




Figure 4. Mouse IgG glycoforms were detected by UV (4a) Using HPLC. Resolution of glycopeptides on the X-Bridge

Glycopeptides with Sialylated N- and O-linked Glycans

Figure 5




Figure 5. Tryptic digest of bovine fetuin analyzed with UPLC on a 300Å pore size UPLC™ column with TFA mobile-phase modifier. Glycopeptide peak shapes are narrow, symmetrical and yield useful spectra.

Figure 6




Figure 6. The intensity of the total ion chromatograms of the fetuin trypsin digest is comparable whether the source is tuned for molecular ions or fragments. Numerous sialylated glycopeptides were observed by fragment analysis. The flow rate in this case was 200 µL/min. Labeled peptides are described in Table 1

Table 1

Peak	Peptide	Major Glycoforms	Resolution
A	T10	O- 1 sialo-hex-hexNAc " +1 sialic acid	Overlapped
B	T11-12	N- Trisialo-Triantennary " +1 sialic acid	Resolved
C	T12	N- Trisialo-Triantennary " +1 sialic acid	Partially Resolved
D	T7	N- Trisialo-Triantennary " +1 sialic acid	Overlapped
E	T20	O- 3 sialo-hex-hexNAc 3 " +1 sialic acid	Partially Resolved
F	T13	N- 2 sialo-hex-hexNAc Trisialo-Triantennary " +1 sialic acid	Resolved Partially Resolved

Figure 7




Figure 7. Deconvoluted mass spectra of the 3 sialic acid-containing region of Peak B. The spectra were collected in the same run at low (4a) or high (4b) aperture1 potentials. Under high energy conditions fragments are observed. With good chromatographic resolution useful glycan structural information is obtained in MS only mode.

Figure 8




Figure 8. The largest glycopeptides observed in these experiments represents missed cleavages in the fetuin digest. The glycoforms have molecular weights approximately 10 kDa, and were separated on the 130Å pore size, 1.7 µm particle packing material.

Figure 9




Figure 9. The fetuin peptide T11-12 has two major glycoforms. Both are complex triantennary structures with either 3 or 4 sialic acid termini. Elution order is reversed when formic acid is replaced with TFA.³

MS Conditions for Glycan Detection

Glycopeptides were detected using the classic method of source fragmentation to create diagnostic glycan ions.² MS scans were alternated between normal source conditions (Aperture 1 = 15V) and fragmentation conditions (Aperture 1 = 50V). The fragment ion, m/z 657, identifies the presence of N- and O-linked sialylated glycans. Sialylated glycan termini, both N- and O-linked, produce m/z 292 (NeuAc) and m/z 657 (NeuAc-Gal-GlcNAc).

Figure 10




Figure 10. Comparison of the separation of tryptic peptides of bovine fetuin. Both the 300Å and 130Å pore size Peptide Separation Technology™ 1.7µm columns provide excellent resolving power for peptides.

CONCLUSION

- Chromatographic resolution provides excellence in mass spectra of glycopeptides.
- High resolution separation of glycopeptides has been achieved on 300Å and 130Å pore size 1.7 µm particles.
- Good peak shape is obtained for sialylated and non-sialylated glycopeptides including positional isomers.
- Changing the gradient slope, temperature or mobile phase modifier can optimize separations.
- Data obtained in MS only mode produced spectra characterizing complex sialylated glycan structures.

References

- Mazzeo, J.R.; Wheat, T.A., Gillece-Castro, B.L. and Lu, Ziling *BioPharm International*, **2006**, 19:1, 1-9
- Huddleston; M.J., Bean; M.F. and Carr, S.A. *Anal. Chem.* **1993**, 65 877-884.
- Medzhiradzky, K.F.; Besman, M.J. and Burlingame, A.L. *Anal Chem.* **1997**, 69, 3986-3994.

TO DOWNLOAD A COPY OF THIS POSTER VISIT WWW.WATERS.COM/POSTERS