

Waters® Alliance™ LC/MS System



LC/MS API Methods Development:

Optimizing existing HPLC methods
for acidic analytes in negative mode

Michael P Balogh, Waters Corporation, Milford, MA

Detection of hypericin
and related compounds

Quantitation

UV Profiles

Analytical Conditions

References

Key words:

Botanicals, St. John's
Wort, Method
Conversion for LC/MS,
Natural Products

Natural products analysis: hypericin and related compounds from plant matrix

Sometimes existing chromatographic methods are less than ideal for adaptation to LC/MS. In this case, a method proposed for the analysis of hypericin, (the suspected active ingredient in the natural antidepressant St. John's Wort prepared from the plant *Hypericum perforatum*, L.) was determined to be unsuitable. The gradient only separated two of the four currently known dianthrone and inhibited detection and quantitation through the use of a strong acid with the weakly anionic hypericin. The method cited here solves a number of problems and provides:

- reproducible quantitation
- resolution of all related dianthrone in the sample
- detection of the unique wavelengths associated with dianthrone
- improved limits of detection for the mass spectrometer

An improved method¹ was derived from published work² for the somewhat anionic hypericin molecule using TEA (triethyl ammonium acetate) resulting in:

- an improved UV profile for the unique dianthrone lambda max at 588 nm ([Figure 1](#)) and,
- highly accurate quantitation of the separated dianthrone ([Figure 2](#)).

Many current studies have adopted 300 mg therapeutic doses containing 0.3% hypericin (or 0.9 mg). Determining adulteration and substandard preparations is of primary importance. In addition, this method provides useful separation of the four currently known biosynthetic and dehydrated forms ([Figure 3](#)).

Chromatographic conditions:

Column: Waters Symmetry[®] C₈, 2.1 mm x 150 mm

Flow Rate: 0.3 mL/min

Mobile phase A: 100-mM triethylamine acetate, pH 7.0 (TEA acetate)

B: methanol

C: acetonitrile

Initial conditions: 30:39:31 (A:B:C) A 15-minute linear gradient to final conditions of 10:50:40 with a 10 minute reequilibration to initial conditions.

¹Balogh, M.P. and Li, J.B., *Analysis of Hypericin with HPLC-PDA-MS: Advantages of Multispectral Techniques*, LCGC, in print spring 1999.

²Piperopolous, G., Lotz, R., Wixforth, A., Schmierer, T., Zeller, K.P., *J.Chrom. B*, 695, 309-316 (1997).

Waters, Alliance and Symmetry are trademarks of Waters Corporation




Figure 1 UV Profile of hypericin with 0.5% TFA mobile phase and 100 mM TEA mobile phase (lambda max 588 nm)




Figure 2 Calibration - 2 ng to 40 ng on-column




Figure 3 Negative mode electrospray acquisition for two samples of St. John's Wort Preparation:

- pseudoprotohypericin (521 m/z; biosynthetic precursor)
- protohypericin (519 m/z)
- pseudohypericin (505 m/z; biosynthetic precursor)
- hypericin (503m/z)

The PDA extracted wavelength (588 nm) is shown in the lower trace for reference.

