# Waters Alliance<sup>®</sup> LC/MS System



Key Words Algaecide Biocide Copper Pyrithione Zinc Pyrithione Electrospray ionization, ESI Toxic Compounds

# Analysis of Toxic Organometallic Compounds: Zn / Cu Pyrithione

Makiko Ichiki and Jun Yonekubo, Nihon Waters, Japan

## Background

Zinc and copper pyrithione (ZnPT and CuPT, respectively) are widely-used biocides used in a variety of applications, such as marine paint for boats. Because the Japanese Environmental Agency (JEA) will add ZnPT to its list of toxic environmental compounds, a suitable analytical method by LC/MS was investigated and developed.

## **Analytical Conditions**

Mass spectrometry employed the Waters ZMD, using positive electrospray ionization (ESI<sup>+</sup>). Data was acquired in both full-scan and Selected Ion Recording (SIR) mode. HPLC utilized a Waters Alliance<sup>®</sup> HT Separations Module and a Waters Symmetry<sup>®</sup> C<sub>18</sub> column (2.1x30 mm, 3.5  $\mu$ m) at a flow rate of 0.3 mL/min. The following solvent program was used for separation:

<u>Solvent A</u> =  $2mM AcONH_4$  <u>Solvent B</u>= (MeOH)

Waters Alliance®HT 2790 HPLC Gradient Timetable

<u>Time</u>	<u>_A%</u>	<u>_B%</u>	<u>_Flow</u> (mL/min)	<u>Curve</u>
0.00	80.0	20.0	0.30	1
10.00	0.0	100.0	0.30	6



Compound 1 name: Zn - Cu PT Coefficient of Determination: 0.998317 Calibration curve: 278.068 \* x + -2234.83 Response type: External Std, Area Curve type: Linear, Origin: Include, Weighting: Null, Axis trans: None



## Figure 1: Mass Spectra of CuPT and ZnPT

Mass spectra of CuPT (top) and ZnPT (bottom) are shown. Both compounds give strong [M+H]<sup>+</sup> ions which differ by one mass unit. Also, note that the isotope pattern for each compound depends upon which metal is chelated.

#### Table 1: Mass Intensity vs. Concentration of CuPT

<u>Name</u>	Conc	<u>Rt</u> (min)	<u>Response</u>	<u>Cal.</u>	<u>%Dev</u>
CuPT Blank	<u>(pbp)</u>	<u>(min)</u>		<u>- conc</u>	
CuPT 5b	5.0	8.51	632	10.3	106.2
CuPT_10b	10.0	8.53	1332	12.6	26.5
CuPT_50b	50.0	8.56	10223	40.4	-19.2
CuPT_100b	100.0	8.56	23341	88.3	-11.7
CuPT_250b	250.0	8.56	68441	254.2	1.7
CuPT_500b	500.0	8.56	137106	501.5	0.2

## Figure 2: CuPT Calibration Curve

Analysis of a series of CuPT standards under SIR conditions (m/z 316) gave excellent linearity (left) and sensitivity over a concentration range from 5-500 ppb. The table above lists the results of this analysis.

## Figure 3: Actual ZnPT Spectra

Zinc PT gave unusual results when the protonated molecular ion (m/z 317) was monitored. Additional experiments showed that ZnPT was being converted to the copper chelate when injected onto the HPLC column. This conversion was not observed when the zinc compound was flow-injected directly into the mass spectrometer.



#### Quantify ZnPT (as CuPT) Summary Report (5.0 - 100 ppb)

<u>Name</u>	Туре	<u>Std.Conc</u> (ppb)	<u>Rt</u> (min)	<u>Response</u>	<u>conc</u>	<u>%Dev</u>
Blank_AcN_0	Blank				-	
ZnPT_5b_1	Standard	5.0	8.48	832	7.4	47
ZnPT_10b_1	Standard	10.0	8.61	1013	9.1	-9
ZnPT_25b_1	Standard	25.0	8.61	2415	22.2	-11
ZnPT_50b_1	Standard	50.0	8.61	5516	51.3	3
ZnPT_100b_1	Standard	100.0	8.61	10704	100.0	0

Compound 1 name: Zn - Cu PT

Correlation coefficient: r = 0.998681,  $r^{2} = 0.997364$ Calibration curve: 106.573 \* x + 46.0665Response type: External Std, Area Curve type: Linear, Origin: Exclude, Weighting: Null, Axis trans: None



#### Figure 4: ZnPT Calibration Curve

When ZnPT was analyzed as the copper analogue and monitored at m/z 316, excellent linearity and sensitivity were obtained over the concentration range from 2-100 ppb. It is felt that this lower dynamic range is related to the conversion from the zinc to copper form.



#### 25 ppb ZnPT x 7 masschromatograms (316 amu)



**Figure 5: Replicate ZnPT Injections** Replicate injections of 25 ppb ZnPT Standard. Peak area reproducibility = 2.4% RSD.

## **Figure 6: Signal to Noise** The signal to noise ratio at 1 ppb ZnPT standard is 5.5:1. Similar data were obtained for CuPT.

# Conclusions

LC/MS offers sufficient sensitivity and linearity to detect and quantitate low ppb levels of both zinc and copper pyrithione. In addition to the quantitative power, LC/MS data clearly demonstrated the conversion of zinc pyrithione to its copper form. It is doubtful if this transformation would have been detected and elucidated with any other type of detector.