

# Waters column

*A Reprint from*

*Spring 1995*

*If you would like a subscription,  
or more information,  
please contact your nearest  
Waters Office.*

*Check the Waters website  
for up-to-date information:  
<http://www.waters.com>*

Published by Waters Corporation  
34 Maple Street  
Milford, MA 01757 USA  
Tel: 508-478-2000  
Toll-free: 1-800-252-4752  
Fax: 508-872-1990

Editor: Uwe D. Neue, Ph.D.  
ISSN #1084-0540


# HPLC Purification of Phosphorothioated DNA for Antisense Therapeutic Investigations

## Objective:

The objective of this application note is to demonstrate the utility of reversed-phase HPLC for large scale purifications of DMT protected synthetic oligonucleotides.

## Details:

Synthetic oligonucleotides are successfully used in a wide variety of applications ranging from use as hybridization probes to utilization as primers for DNA sequencing and the polymerase chain reaction\* (PCR). Most recently, studies have shown the effectiveness of using standard (i.e. phosphodiester, P=O) or phospho backbone modified (i.e. phosphorothioated, P=S) oligonucleotides to selectively inhibit the production of targeted proteins within cells. As such, "antisense and related technology investigations" require milligram to multigram amounts of purified synthetic DNA. Reversed-phase HPLC provides a rapid and cost effective technique for the isolation of either phosphodiester or phosphorothioated full length oligonucleotide products for these applications.



*Milligram quantities of DMT protected synthetic oligonucleotides can be rapidly purified from undesired synthesis reaction by-products using Bondapak HC<sub>18</sub> HA reversed-phase packing in an 8 x 100mm Radial-Pak cartridge. This same packing is available in a variety of column dimensions and in bulk packaging for gram scale purifications as required for antisense therapeutic investigations.*

## System:

The HPLC system used consisted of a Waters W600E Multisolvant Delivery System, Waters 486 U.V. Detector, and Waters 860 ExpertEase<sup>™</sup> Data System (signal via SIM/LAC/E<sup>™</sup>).

## References:

Agrawal, S., J. Goodchild, M.P. Civeira, A.H. Thornton, P.S. Sarin and P.C. Zamecnik. 1988. "Oligodeoxynucleoside Phosphoramidates and Phosphorothioates as Inhibitors of Human Immunodeficiency Virus", Proceeding of the National Academy of Sciences USA. Vol 81. Pages 3297-3301.

Van der Krol, A.R., J.N.M. Mol and A.R. Sluitje. 1988. "Modulation of Eukaryotic Gene Expression by Complementary RNA or DNA Sequences", Bio Techniques. Vol 6. Pages 958 - 973.

Moffat, A.S. 1988 (Nov/Dec). "Researchers Pursue Anti-Sense Technology In Quest for Novel Drugs and Agriproducts", Genetic Engineering News. Vol 8. Pages 1 and 7.

Klauser, A. 1990. "Antisense Start-Ups Surveyed", Bio/Technology. Vol 8. Pages 303 - 304.

\* See U.S. Patent No. 4683202 to Cetus Corporation