

Waters column

A Reprint from

Spring 1995

*If you would like a subscription,
or more information,
please contact your nearest
Waters Office.*

*Check the Waters website
for up-to-date information:
<http://www.waters.com>*

Published by Waters Corporation
34 Maple Street
Milford, MA 01757 USA
Tel: 508-478-2000
Toll-free: 1-800-252-4752
Fax: 508-872-1990

Editor: Uwe D. Neue, Ph.D.
ISSN #1084-0540

Rapid Polyolefin Additive Separation Using a Waters Nova-Pak[®] C₁₈ Column

Objective:


This separation demonstrates that flow programming, along with a gradient is a powerful technique for separating a large number of polymer additives in a short amount of time. This gradient does not hydrolyze the phosphite ester into multiple peaks, as evidenced by the single, symmetrical Irgafos 168 peak.

Details

See reference for microwave and ultrasonic extraction methods.

System

The system used consisted of a Waters Powerline 600 Solvent Delivery System, with the Waters 490 UV Detector, Waters 712 WISP Autosampler, and the Waters 860 ExpertEase[™] Data System.



Gradient HPLC with flow programming can be used to separate these 12 polyolefin additives in less than 17 minutes. The Irgafos 168 phosphite ester appears as a single peak, with no evidence of degradation.

References:

Nielson, R., "Extraction and Quantitation of Polyolefin Additives", Journal of Liquid Chromatography, Vol. 14, #3, 1991

