

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

¹ Waters Corporation, Milford, MA 01757; ² Millipore Corporation, Bedford, MA, USA

METHODS

A chimeric anti-phosphotyrosine IgG1 mAb (PTG1) was expressed in two CHO cell lines (CHO-S and DG44) and purified by Protein A chromatography using two different protocols. Null cell cultures, (expressing no PTG1 protein), were grown under identical conditions and Protein A purified. Five protein standards were spiked in PTG1-RT and the resulting protein mixture was denatured with 0.1% RapiGest (15 min at 60 °C), reduced with 10 mM dithiothreitol (DTT) for 30 min at 60°C, alkylated with 20 mM iodoacetamide (IAM) for 30 min at 37°C, and enzymatically digested with Promega trypsin (1:20 w/w ratio) at 37°C overnight. Following digestion, the protein amounts loaded on-column on the 2D-LC system were: **4,000 fmoles** ALB (bovine alpha-lactoglobulin), **800 fmoles** PHO (rabbit glycogen phosphorylase b), **320 fmoles** ADH (yeast alcohol dehydrogenase), **80 fmoles** BSA (bovine serum albumin), and **16 fmoles** ENL (yeast enolase).

A nanoACQUITY™ UPLC® system with 2D technology was used for peptide separations. A reversed-phase/reversed-phase (RP/RP) method was developed that uses the pH of the mobile phases to change the selectivity of a peptide separation in two separate dimensions [1,2].

- **First Dimension (1D) pH=10:** 1.0 mm x 50 mm XBridge C₁₈ column (5 µm particles); flow 10 µL/min. Mobile phase 20 mM ammonium formate in water (Solvent A) and ACN (Solvent B).
- **Online dilution** (1:10) of the eluent from 1D before analyte trapping onto the 2D column.
- **Trap column:** A 5-µm Symmetry C₁₈ trap (2.7 µL volume) was used to trap peptides between the two LC dimensions.
- **Second Dimension (2D) pH=2.4:** 0.3 mm x 150 mm analytical column BEH C₁₈ 1.7 µm, kept at 65 °C and operated at 12 µL/min.

Fractions were eluted in ten steps (1: 10.8% Eluent B, 2: 12.4%, 3: 14.0%, 4: 15.4%, 5: 16.7%, 6: 18.6%, 7: 20.4%, 8: 25.0%, 9: 30.0%, 10: 50.0% B. Each step was mixed in a 1:10 ratio with 0.1% TFA in water (pH=2.1) before trapping. Low pH separations in the second chromatographic dimension used a 30 min gradient from 7 to 35% acetonitrile (0.1% FA-formic acid). See Figure 1 for the diagram of the 2D-LC system.

A) Discovery and Label-free Quantitation: Data independent, alternate scanning LC/MS^E experiments were performed with a SYNAPT HDMSTM mass spectrometer.

- Acquisition time was 0.5 sec
- m/z range: 50-1990
- Fixed CE at 5 eV for MS; CE ramp 15 - 35 eV for high-energy fragmentation
- Database search—ProteinLynx Global Server (PLGS) 2.4 with Identity^E Informatics [3].

B) Targeted Quantitation: MRM analyses were performed on a Waters XEVOTM TQ MS tandem quadrupole mass spectrometer coupled with the ACQUITY UPLC system.

- 2.1 x 150 mm BEH130 C₁₈ column (1.7 μm particles) at 35 °C.
- Flow rate of 300 μL/min; 3-40% B over 10 min (B: acetonitrile with 0.1% FA; A: water with 0.1% FA)
- Other parameters: Potential 3.5 kV, CV 37 V, source temp 90 °C, MS1/MS2 isolation window 0.75 Da (FWHM), adjustable collision energy depending of precursors/fragment (15-30 eV) and 10 ms dwell time.

INTRODUCTION

- Residual host cell proteins (HCPs) from recombinant production are a major component of biopharmaceutical process-related impurities. HCPs can elicit an unpredictable immune response in patients.
- The composition of HCPs is extremely heterogeneous and changes with production and purification procedures.
- Current analytical methods for measuring HCPs, (typically ELISA, gels, blots), are expensive, subjective, time-consuming to develop, and require prior knowledge about the identities of the contaminant proteins.

- Spike a mixture of 5 proteins in the concentration range of 10-1,000 ppm (ng/mg) in the protein biopharmaceutical
- Enzymatic digestion of the protein sample to peptides
- 2D-LC/MS^E with IDENTITY^E to discover HCPs
- Use the top 3 peptides from each protein for absolute label-free quantitation based on the known concentration of the spiked proteins
- Use VERIFY^E for mining the MS^E dataset to extract appropriate MRMs for monitoring HCPs by Xevo TQ MS
- Analyze the sample digests by tandem quadrupole MS for targeted, high-throughput quantitation (e.g. using isotopically labeled peptides)

Chromatographic Performance of the 2D-LC System




Figure 3. Reproducibility of 2D chromatography: extracted mass chromatograms of the ENL T43 peptide (VNOIGTLESSEK, $[M+2H]^{2+}=644.86$) from four consecutive 2D-LC experiments. All second dimension chromatograms were performed at 12 $\mu\text{L}/\text{min}$ using a 30 min gradient (3–35% ACN, 0.1% FA). T43 peptide eluted only in Fraction 3/5 (16.7% B₀) in the 1st dimension chromatography. The amount of ENL digest loaded on column in each 2D experiment was 20 fmoles.




Figure 4. Chromatographic performance (e.g. RT reproducibility and peak width) is maintained during 1D fractionation: mass chromatograms of EN1-T43 peptide obtained under four fractionation conditions:

- (A) "simulated" 1D run using a single elution step (from 10.8 to 50% ACN)
- (B) fraction 2 out of 3 (from 10.8 to 18.6% ACN)
- (C) fraction 4 out of 5 (from 15.4 to 18.6% ACN)
- (D) fraction 5 out of 10 (from 15.4 to 16.7% ACN)




Figure 5. Performance of high-pH fractionation: mass chromatogram of ENH T43 peptide recorded for 3 successive fractions during a 10-step 2D fractionation experiment: fraction 4/10 corresponds to a step elution from 14.0 to 15.4% ACN, fraction 5/10 was eluted from 15.4 to 16.7% ACN and fraction 6/10 was eluted from 16.7 to 18.6% ACN.

- **LC-MS^E Data:** a combined MS^E dataset collected from 10-step 2D RP/RP separations performed in triplicate.

- **Protein Identification Informatics:** PLGS 2.4 featuring Identity^E Bioinformatics.
- **Database:** protein database was compiled from 12,943 mouse Swiss Prot entries, 654 Golden hamster proteins, the sequences of LA, ADH, PHO, BSA, ENL, porcine trypsin, protein A (*S aureus*) and the heavy and light chain sequence of PTG1. The final randomized database contained 27,212 entries.

Pro	Protein Description	CHO S Cells				DG-44 cells						
		I	II	III	IV	A1	B1	A2	B2	C	D	
1	Nascent polypeptide associated complex subunits Mus musculus		50			1739	36	2189		3725		
2	Nucleolin Mesocricetus auratus					1511	35	3588		2115		
3	Heterogeneous nuclear ribonucleoprotein isoform 1 Mus musculus					1365	24	3681		2115		
4	Elongation factor isoform Mus musculus					1241	18	1697	30	2472		64
5	Procollagen C endopeptidase enhancer 1 Mus musculus		142	665		2043	18	1697	30	2472		
6	Laminin subunits Mus musculus					1260	12	161	45			
7	Actin isoform Mesocricetus auratus					1125	1780	1033	1134	1558		115
8	Calnexin Mus musculus		537	1339		1365	24	3681		2115		25
9	Nidogen 1 Mus musculus					71	4	826	47	141		
10	Glyceraldehyde 3 phosphate dehydrogenase (muscle isoform)					1558	59	75	73	646	88	92
11	Sialic acid transferase 1 Mus musculus					71	4	826	47	141		
12	Glycogen phosphorylase b rabbit - PHO (BIO 8000)		778	778	778	778	778	778	778	778	778	778
13	Plasminogen activator 1 H19 Mus musculus					1125	1780	1033	1134	1558		
14	Lipoprotein lipase Mesocricetus auratus						884	547	1149	707		
15	75 kDa glucose regulated protein Mesocricetus auratus						210	628	546	251	244	56
16	Peptide arginyl tyrosine phosphatase Mus musculus					53	2	342	120	222		
17	T complex protein 1 subunits Mus musculus					51	3	342	120	222		
18	Serine protease HTRA1 Mus musculus					71	3	342	120	222		
19	High mobility group protein isoform Mus musculus		74	78		71	3	342	120	222		
20	80S ribosomal protein S3 Mus musculus					595	155	312	149	224		
21	Alpha-thalassaemia 1 Alpha Box 1 (2000 fmoles)		200	130	502	281	319	312	149	224	447	358
22	Lysosomal alpha glucosidase Mus musculus					708	711	312	149	224		
23	Hemostatic derived growth factor Mus musculus		73			229	304	304	304	235		
24	Nucleolin sensitive inhibitor of protein kinase C Mus musculus					58	29	304	304	235		
25	Pyruvate kinase isoform M2 Mus musculus					24	612	22	300	207		39
26	Activated RNA polymerase I transcriptional coactivator p15					24	216	12	95	194		
27	Heat shock cognate HSP 90 beta Mus musculus					125	117	117	117	138		
28	Nucleophosmin Mus musculus					125	117	117	117	138		
29	Insulin like growth factor binding protein 4 Mus musculus		64			50	10	103	95	121		
30	Complement C1 protein necrosis factor related protein 4					68	93	123	67	237		
31	Eukaryotic translation initiation factor isoform Mus musculus					21	38	118	33	103		
32	Alcohol dehydrogenase Y-Ado (200 fmoles)		101	98	111	98	109	99	110	118	97	
33	Interleukin enhancer binding factor 2 Mus musculus							96	90	103		
34	Guanine nucleotide binding protein subunits Mus musculus							104	92	65		
35	Vimentin Mesocricetus auratus							85	79	91		69
36	Endothelial differentiation related factor 1 - BSA							50	71			
37	Serum albumin precursor Allergen Bos 1-100 (60 fmoles)		61	57	67	56	57	68	57	69	55	28
38	Collin Mus musculus							60	58	39		31
39	Penicillinase 1 Mus musculus							60	58	39		24
40	Helicase isoform Mus musculus					509		599		599		28
41	Penicillinase B Mus musculus							95	40			24
42	Peptidyl prolyl cis trans isomerase A Mus musculus							92	44			39
43	Nonreducing chromosomal protein HMG 17.1 Mus musculus		95	94								
44	60S ribosomal protein L23A Mus musculus						17		35			
45	DNA binding protein A Mus musculus							16	28			
46	60S ribosomal protein L30 Mus musculus								23	181		
47	Ubiquitin associated protein 2 Mus musculus		30	74		503						
48	Tubulin isoform Mus musculus		158	75								
49	Thrombospondin 1 Mus musculus											
50	Polyadenylation binding protein 1 Mus musculus					373						19
51	Heat shock protein HSP 90 alpha Mus musculus								89			
52	Heat shock cognate 71 kDa protein Mus musculus								152			68
53	Enolase yeast - ENL (16 fmoles)											
54	Enolase Mus musculus		6	10								8
55	Enolase Mesocricetus auratus											

Table 1. HCP concentrations (ng/mL) measured in PTG1 mAb expressed in CHO (samples I-IV) and DG-44 cells (samples A-D). Two protein A purification protocols were compared across ten samples: I vs II - CHO cells grown under regular conditions; III vs IV - CHO cells non-expressing the PTG1 target; A1/2 vs B1/2 biological replicates for DG-44 cells grown under regular conditions; C vs D - non-expressing DG-44 cells. High-abundance proteins are highlighted in RED (>1,000 ng/mL), medium-abundance are either in YELLOW (500-1000 ng/mL) or GREEN (100-500 ng/mL) and low-abundance HCPs are GREY (< 100 ng/mL).

Sample ID	Cell line	PTG1 expressed?	No of HCPs Identified	Common HCPs	HCP concentration (ng/mL)	
					Mass Spec	ELISA
I II	CHO	Y	5	4	4,870	7,013
	CHO	Y	14	4	9,752	12,228
III IV	CHO	N	0	0	0	30
	CHO	N	31	0	9,093	4,138
A1 B1	DG-44	Y	34	21	61,138	3,684
	DG-44	Y	50	21	109,437	8,042
A2 B2	DG-44	Y	22	15	35,031	4,282
	DG-44	Y	32	15	105,123	8,041
C D	DG-44	N	0	0	0	6
	DG-44	N	28	0	1,454	1,318

Table II. Comparison of total HCP concentrations (ng/mL) measured by the mass spectrometry and ELISA (using the Cygnus Technologies CHO HCP mammalian kit).

The MS^E data was automatically mined using VERIFY^E software for the best proteotypic peptides from 20 HCPs. The MRM assay, implemented on a triple quadrupole instrument, used a 15-min gradient for monitoring 58 MRMs from 29 peptides (2 transitions per peptide), belonging to twenty HCPs.




Figure 6. Correlation between the MS^E and MRM data: (A) concentration of 78 kDa glucose regulated protein measured across 10 different PTG1 preparations (samples I-IV and A-D) measured using the MS^E data; (B) peak areas produced by two MRM transitions of peptide VEI-IANDQGNR from the same protein.




Figure 7. Histogram plot showing the peak area RSDs for 58 MRM transitions measured in 5 replicates for sample B1. The average RSD was 13.8%.

CONCLUSIONS

- The LC/MS^E-based assay allows the identification and quantification of low-abundance HCP contaminants in biopharmaceuticals over four orders of magnitude in concentration.
- Protein A purification of mAbs using different purification protocols produces different HCP patterns. Cell lines used for mAb purification can affect the HCPs identity and concentration significantly.
- ELISA assay seems to underestimate total HCP concentration.

References

1. Gilar M, Olivova P, Daly AE, Gebler JC *J. Sep. Sci.*, **2005**, 1694.
2. Stapels MD, Fadgen K Current Trends in Mass Spectrometry, *Spectroscopy* supplement, March **2009**.
3. Li GZ, Vissers JP, Silva JC, Golick D, Gorenstein MV, Geromanos SJ *Proteomics*, **2009**, 1696.