

Tabisam Khan¹, Gordon Kearney¹ Amir Farook¹ and Sheng Long-Sheng²¹ Waters Corporation, Manchester, UK. ² Chinese Pharmaceutical University, 24 Tong Jia Xiang, Nanjing, China.

Presented at ASMS 2002, Orlando, Florida, USA, 3rd - 6th June 2002

OVERVIEW

- *Ginkgo biloba* is a tree that has, since ancient times, been used as a component of traditional Chinese herbal medicines.
- Ginkgolides, the main components of preparations of Ginkgo leaf and root, are now being developed as therapeutic agents.
- In order to quantify known components of *Ginkgo biloba* leaf, fast gradient HPLC separations were coupled to MRM analyses on a Quattro Micro triple quadrupole mass spectrometer (Micromass UK Ltd., Manchester, UK).
- Quantification methods show good sensitivity, dynamic range and precision.
- In order to ascertain the identity of other components in the leaf preparation a slow HPLC gradient was coupled to high-resolution, accurate-mass analysis using a Q-ToF Ultima hybrid quadrupole/oaToF mass spectrometer (Micromass UK Ltd., Manchester, UK).
- A range of compounds and structural isomers were identified in the mixture.

INTRODUCTION

Medical treatment has long been connected with the natural products and actives extracted from plants. This has lead to many advances in medicines and has long been the pivotal force in fighting ailments. It has been known and appreciated for many years that certain plants can exhibit healing effects when prepared as part of a mixture, but when taken in isolation can cause detrimental effects. This knowledge is now being used more and more in sophisticated biological and bioanalytical studies of natural products.

Despite advances both in utilising synthetic approaches to drug design, and in sophisticated structure-activity studies, there is still a great need for compounds with a unique mechanism of action. Major breakthroughs have resulted primarily from the study of natural products. Some of the most important drugs have been isolated from plant sources; for example most antibiotics and anti-cancer drugs.

The maidenhair tree, *Ginkgo biloba*, is an ancient Chinese plant that has been cultivated for its health-promoting properties. Ginkgolides, the main active ingredients of *Ginkgo biloba*, have not only helped to explain the pharmacological basis of several traditional medicines, but have also provided a valuable new class of therapeutic agents. Research on the biochemical effects of *Ginkgo biloba* extracts is still at a very early stage. Although the terpene fraction of *Ginkgo biloba*, which contains the ginkgolides, may contribute to the neuroprotective properties of the *Ginkgo biloba* leaf, it is also likely that the flavonoid fraction, containing free radical scavengers, is important in this respect. The structures of the ginkgolides are shown (Figure 1).




Figure 1.

EXPERIMENTAL**Extraction Procedure**

The extraction procedure used by the China Pharmaceutical Company is described as follows.

- 50g Ginkgo Biloba leaf powder was refluxed with 200ml 60% acetone for 1hr, then filtered.
- Extraction solvent was recovered and remaining residue was extracted with ethyl acetate.
- Ethyl acetate was recovered and remaining residue was extracted with water.
- Water was extracted with diethyl ether.
- Acetone, ethyl acetate and diethyl ether fractions were combined and evaporated.
- Samples were reconstituted in 60% acetone.

For the commercial samples obtained from the UK the method above was not necessary as the samples were obtained in tablet form. They were simply dissolved in 10mL of mobile phase and centrifuged to remove all the stabilisers. The resulting solvent was decanted and filtered. The samples were then ready for injection.

Quattro Micro*HPLC Pump Conditions**Mobile Phases*

A	Acetonitrile
B	Water

The gradient Timetable contains 6 entries, which are:

Time	A%	B%	Flow (ml/min)
0.00	95.0	5.0	0.300
1.00	50.0	50.0	0.300
2.50	50.0	50.0	0.300
3.00	5.0	95.0	0.300
7.00	5.0	95.0	0.300
7.50	95.0	5.0	0.300

Injection Volume (µl)	5.0
Stop Time (mins)	12.00

Symmetry 100 x 2.1mm C₁₈ 3.5µ HPLC column (Waters Corp., Milford, MA, USA).

Mass Spectrometer Conditions

Polarity	ES-
Capillary (kV)	3.00
Cone (V)	20.00
Extractor (V)	3.00
RF Lens (V)	0.5
Source Temperature (°C)	150
Desolvation Temperature (°C)	350
Cone Gas Flow (L/Hr)	64
Desolvation Gas Flow (L/Hr)	671

Q-ToF Ultima*HPLC Pump Conditions**Mobile Phases*

A	Acetonitrile
B	Water

The gradient Timetable contains 6 entries, which are:

Time	A%	B%	Flow (ml/min)
0.00	95.0	5.0	0.300
2.00	95.0	5.0	0.300
60.0	40.0	60.0	0.300
61.0	5.0	95.0	0.300
65.0	5.0	95.0	0.300
66.0	95.0	5.0	0.300

Injection Volume (µl)	5.0
Stop Time (mins)	72.00

Symmetry 100 x 2.1mm C₁₈ 3.5µ HPLC column (Waters Corp., Milford, MA, USA).

Mass Spectrometer Conditions

The instrument was operated using a Lockspray ionisation source and W mode ToF ion optics.

Polarity	ES-
Capillary (kV)	3.00
Cone (V)	60.00
Source Temperature (°C)	120
Desolvation Temperature (°C)	350
Cone Gas Flow (L/Hr)	50
Desolvation Gas Flow (L/Hr)	600
Resolution (FWHM)	20,000
Scan Range (m/z)	100 to 900
Acquisition Mode	Centroid
Reference Mass	Leucine Enkephalin at m/z 554.261

RESULTS

Quattro Micro

Calibration lines were constructed over the range 0.5pg/uL to 5000pg/uL for a mixture of the three standards, with each individual calibration point being injected in duplicate. The resulting coefficients of correlation are shown below:

Ginkgo_A	$R^2=0.999370$
Ginkgo_B	$R^2=0.997926$
Bilobalide	$R^2=0.999375$
Ginkgo_A_isomer	$R^2=0.999726$

Each calibration line showed a good linear response with a high degree of precision for the method. The calibration graphs are shown (Figures 2, 3, 4 and 5).




Figure 2. Bilobalide




Figure 3. Ginkgolide A




Figure 4. Ginkgolide B




Figure 5. Ginkgolide A (isomer)

The samples were then injected in duplicate and the mean amount of active ingredients are shown (Table 1)

Samples	Ginkgo A	Ginkgo B	Bilobalide	Ginkgo A isomer
Commercial sample1	1722.1	884.6	889.3	31.5
Commercial sample2	328.0	185.9	172.4	5.8
Nanjing sample 1	1366.2	1122.9	N/A	58.1
Nanjing sample 2.	1323.3	1138.0	N/A	56.2

Table 1.

Q-ToF Ultima

The base peak intensity (BPI) chromatogram for the analysis of an extract of ginkgo leaf is shown (Figure 6).




Figure 6.

A variety of components can be seen and their accurately determined masses, together with matched elemental compositions are given (Table 2). This allows many of the peaks to be identified.

Table 2.

Retention Time/min	Measured Mass	Elemental Composition	ppm deviation	Identity
6.36	445.21	C ₂₁ H ₃₃ O ₁₀	6	
6.53	455.1202	C ₂₀ H ₂₃ O ₁₂	2.7	Gingkolide D?
7.63	325.0921	C ₁₅ H ₁₇ O ₈	-0.6	Isomer of Bilobalide
9.15	429.2134	C ₂₁ H ₃₃ O ₉	2.3	
9.53	441.1395	C ₂₀ H ₂₅ O ₁₁	-0.4	
9.97	441.1808	C ₁₄ H ₃₃ O ₁₅	-2.6	
10.62	465.2339	C ₂₁ H ₃₇ O ₁₁	0.7	
11.19	425.1465	C ₂₀ H ₂₅ O ₁₀	4	
11.79	451.255	C ₂₁ H ₃₉ O ₁₀	1.6	
12.56	325.0918	C ₁₅ H ₁₇ O ₈	-1.8	Isomer of Bilobalide
13.75	425.1454	C ₂₀ H ₂₅ O ₁₀	1.4	
14.1	325.0918	C ₁₅ H ₁₇ O ₈	-1.6	Bilobalide
14.5	609.1453	C ₂₇ H ₂₉ O ₁₆	-0.5	
15.95	423.129	C ₂₀ H ₂₃ O ₁₀	-0.2	Isomer of Gingkolide B (Gingkolide J or M?)
16.5	439.1239	C ₂₀ H ₂₃ O ₁₁	-0.3	Gingkolide C
17.18	447.0936	C ₂₁ H ₁₉ O ₁₁	2	
18.56	593.1474	C ₃₄ H ₂₅ O ₁₀	4.4	
19.65	431.0984	C ₂₁ H ₁₉ O ₁₀	1.4	
21.05	377.0859	C ₁₈ H ₁₇ O ₉	-3.5	
21.82	411.1651	C ₂₀ H ₂₇ O ₉	-1	
22.15	423.1296	C ₂₀ H ₂₃ O ₁₀	1.1	Isomer of Gingkolide B (Gingkolide J or M?)
23.03	423.1295	C ₂₀ H ₂₃ O ₁₀	0.9	Gingkolide B
23.07	407.1331	C ₂₀ H ₂₃ O ₉	-1.6	Gingkolide A
26.76	551.3066	C ₂₆ H ₄₇ O ₁₂	-0.3	
32.12	405.1201	C ₂₀ H ₂₁ O ₉	3.8	

Table 2.

CONCLUSION

Ginkgo Biloba's main components of interest, Gingkolide A (including the isomer), B and the Bilobalide were quantified using the methodology described. The results showed that commercial samples and the extracted Nanjing samples can be quantified using this method using the Quattro Micro tandem mass spectrometer.

The data from the Q-ToF highlighted the accurate mass measurements use in identifying the main components of the Ginkgo Biloba. The active ingredients were accurately mass measured to within 5ppm and the unknown compounds were also accurately mass measured and a structure was proposed. Further analysis of the unknowns can result in viable structures being proposed. With intense data mining possible metabolism, breakdowns and transformations can be determined.

**Author to whom all correspondence
should be addressed:**

Amir Farook
Waters Corporation
(Micromass UK Limited)
Floats Road, Wythenshawe
Manchester, M23 9LZ
Tel: + 44 (0) 161 946 2400
Fax: + 44 (0) 161 946 2480
e-mail: amir.farook@micromass.co.uk

WATERS CORPORATION
34 Maple St.
Milford, MA 01757 U.S.A.
T: 508 478 2000
F: 508 872 1990
www.waters.com

Made in the United Kingdom

Waters
RIGHT ON TIME.

