

James I. Langridge¹, A. Millar¹, P. Young¹, R. O'Malley¹, N. Swainston¹, J. Skilling², J. Hoyes¹, and K. Richardson¹¹ Waters Corporation, Floats Road Wythenshawe, Manchester.² Maximum Entropy Data Consultants, Cambridge, UK.

Presented at ASMS 2002, Orlando, Florida, 3rd-6th June, 2002

Overview

- The aim of the following experiment was to demonstrate the utility of a *de novo* sequencing algorithm automatically applied to an entire LC-MS/MS dataset.
- A mixture of two proteins, present in equimolar amounts, was digested with trypsin and analysed by LC-MS/MS. The digest mixture was injected into, and separated by, a capillary liquid chromatography system prior to analysis with a Q-ToF *Ultima* API mass spectrometer.
- LC-MS/MS data were acquired using a nanoelectrospray source to provide exact mass measurement.
- All MS/MS data was acquired, processed and analysed in an automated manner.
- The MS/MS continuum data were filtered, combined, and transformed into singly charged, mass measured spectra.
- After conversion to XML, these spectra were *de novo* sequenced using a Bayesian probability based algorithm.
- These results were compared to the results of searching the same dataset against the Swiss-Prot database.

Introduction

Mass spectrometry has rapidly become the method of choice for the identification and characterisation of proteins. Large quantities of MS/MS spectra can be acquired and database-searched in an automated manner from LC-MS/MS experiments, allowing rapid identification of multiple proteins contained within a single sample.

A challenge remains, however, in that many of the acquired MS/MS spectra do not provide matches when searched against known protein sequences. The nature of database searching is such that only peptides whose molecular weight matches those within the databank will be identified. Consequently, many good quality spectra containing a single amino acid substitution remain unmatched. Furthermore, for organisms whose genomes are not completely sequenced electrospray ionisation tandem mass spectrometry (ESI-MS/MS) is an ideal technique since it can be used to provide high quality sequence data from individual peptides produced by the enzymatic digestion of protein. This high quality sequence information can either be used to probe existing protein or nucleotide databases using BLAST™ based homology searching or used for the generation of oligonucleotide primers (cDNA). The enhanced sensitivity, resolution and mass measurement accuracy in data produced by the Q-ToF *Ultima* API mass spectrometer considerably aids the process of inferring *de novo* amino acid sequence. Under favourable circumstances (e.g. tryptic peptides) it is possible to generate long unambiguous stretches of amino acid sequence.

In this work we describe a new fully automated program, that removes the need for manual intervention by automatically applying the *de novo* algorithm to an entire LC-MS/MS dataset. The resultant sequences are compared to the results of performing a database search using the same dataset. The entire process is controlled via a single Java interface.

Experimental

Data acquisition

All data were acquired using a Q-ToF *Ultima* API, hybrid quadrupole orthogonal acceleration time-of-flight mass spectrometer (www.micromass.co.uk).

The sample consisted of two proteins tryptically digested and mixed in equimolar amounts. 100 fmole aliquots of the samples were injected onto the LC-MS/MS set-up.

The analytical system used for the analysis consisted of a ten port valve (the stream select module) of a Micromass CapLC attached directly to the **Z** SPRAY™ source of the mass spectrometer. This was configured with a pre-concentration column in series with a nanoscale analytical column. The trapping column was packed with a C₁₈ stationary phase (300 µm ID x 5 mm), the analytical column used was a 150mmx 75µm column (www.lcpackings.com) packed with PepMap C₁₈ material.

This was set up to spray via a transfer capillary into the ESI source of the mass spectrometer through a nano-LC sprayer. Approximately 3000 volts were applied to the spray tip. A small amount of nebulising gas (approx. 5 psi) was also introduced around the spray tip to aid the electrospray process. A splitter gave a resultant flow through the analytical column of 200 nL/min with the pump programmed to deliver a flow of 2 µL/min.

The mass spectrometer was operated in a data dependant acquisition (DDA) mode whereby following the interrogation of MS data, ions were selected for MS/MS analysis based on their intensity and charge state. Collision energies were chosen automatically based on the m/z and charge state of the selected precursor ions.

All data were acquired with an internal reference ion in order to provide mass measurement accuracies of less than 5 parts per million (ppm) in both the MS and MS/MS mode of acquisition.

Data processing

ProteinLynx™ Global SERVER 2.0 was used to define a 'Workflow' template, comprising:

- the post-acquisition processing routines required to reduce the raw continuum MS/MS data set to a database-searchable form (**figure 1**);
- the *de novo* sequencing parameters (**figure 3**).
- the database to be searched and the associated query parameters (**figure 2**).




Figure 1. The post-processing parameter setup dialog




Figure 2. The database query setup dialog

100fmol pyruvate kinase

2.0fMol ES+ TIC 2.6fMol

3.0fMol ES+ TIC 2.0fMol

4.0fMol ES+ TIC 2.6fMol

1.1pfmol MBP 1.4fMol ES+ BP1 1.4fMol

Time

Figure 4. The de novo sequencing setup dialog

The chromatogram from the LC ESI- MS/MS analysis is shown in **Figure 4** where it can clearly be seen that there are many peptides eluting from the column in the 25 minute period, between 10 to 35 minutes.

The data were processed and searched against the SwissProt v. 39 (86,865 entries, FASTA format) database according to the parameters defined in the WorkFlow template:

- | M/z | SwissProt result | mass ppm | delta | de novo result | mass ppm | delta | matching database entry |
|---------|--------------------------|----------|-------|-----------------------|----------|-------|-------------------------|
| 420.768 | APIIAVTR | 2.971 | | APIIAVTR | 2.971 | | KPYI_RABIT |
| 434.743 | MOHLIAR | 5.189 | | MOHLIAR | 5.189 | | KPYI_RABIT |
| 439.241 | LFEEALR | -8.089 | | LFEEVOL | -8.089 | | KPYI_RABIT |
| 448.241 | ADQGRPEPQVIK | 3.378 | | ADQOLLCELLPI | -10.85 | | ALFA_RABIT |
| 469.234 | AAAEYEVK | 3.062 | | QEYEVNR | -23.784 | | ALFA_RABIT |
| 470.745 | ELSDIAHR | 3.219 | | ELSDLAHR | 3.219 | | ALFA_RABIT |
| 473.259 | VNLAMNVGK | 9.486 | | VNLAMNVGK | 9.486 | | KPYI_RABIT |
| 495.758 | GSQTAEVELK | 2.591 | | GSQTAEVELK | 2.591 | | KPYI_RABIT |
| 503.289 | KLFEEALR | 3.078 | | KLFEEALR | 3.078 | | KPYI_RABIT |
| 510.263 | GDYPLEAVR | -1.99 | | TAYPELAVR | 33.739 | | KPYI_RABIT |
| 559.807 | GSQTAEVELKK | -0.447 | | GSQTAEVELKK | -0.447 | | KPYI_RABIT |
| 571.301 | GDLGIEPAEK | 3.718 | | LDGLDELPAEK | 3.718 | | KPYI_RABIT |
| 577.292 | SSGTSYPDVLK | -1.813 | | SSGTSYPDVLL | -1.813 | | TRYI_BOVIN |
| 586.322 | LDIDSPATIR | -3.398 | | LDLDSALPTAR | -3.398 | | KPYI_RABIT |
| 588.997 | KGVNLPGAADVLP
PAVSEK | 5.968 | | RELLPGAADVLP
AVSEK | 5.968 | | KPYI_RABIT |
| 599.292 | ITLDNAYMEK | 5.324 | | ITLDNAYMEK | 5.324 | | KPYI_RABIT |
| 607.977 | RFDEILASDGIM
VAR | 0.843 | | MVNMWLEASD
GLQLL | -11.008 | | KPYI_RABIT |
| 634.621 | MNFSHGTHEYHA
ETIK | 6.923 | | MNFSHGTFNME
NMTLK | -6.545 | | KPYI_RABIT |
| 666.853 | GILAADESTGSIK | 2.11 | | LGLAADESTTEL
L | -6.323 | | ALFA_RABIT |
| 708.027 | QKGPDFLVTVEEN
GGFLGSK | 7.879 | | KRDEEKVTEVE
NYLEQL | 5.989 | | KPYI_RABIT |
| 724.904 | VYVDDGLISLVQK | -0.089 | | VYVDDGLLSLEL
L | -7.852 | | KPYI_RABIT |
| 745.853 | LQSIGTENTEENR | 7.008 | | LQSLGTEEDEEK
L | 16.339 | | ALFA_RABIT |
| 818.95 | GVNLPGAADVLP
AVSEK | -0.764 | | GVNLPQAVDLA
EGLLGA | -0.764 | | KPYI_RABIT |
| 821.892 | PDVQEAWAEDVD
LR | -3.041 | | PDVQEAWAEDV
DLR | -3.041 | | KPYI_RABIT |
| 833.41 | FDEILASDGIMV
AR | 2.179 | | DFELLEAEELLE
CL | -11.318 | | KPYI_RABIT |

Table 1. Comparison of the results obtained by searching the SwissProt databank and applying the *de novo* sequencing algorithm

Table 2. *de novo* sequencing results obtained for peptides not matched by searching the SwissProt database. blastp searching of the sequences produced the protein matches displayed

Results

Of the 30 ions searched against the database, 25 were matched to either Pyruvate Kinase or Fructose-bisphosphate Aldolase A. These assignments were validated by the post-search filter.

- For 11 ions, the *de novo* algorithm produced identical sequences to those returned by database searching; differing only by I or L (these residues being isobaric and therefore indistinguishable by low energy MS/MS). See example in **Figure 5**.
- In 14 cases, a high degree of homology was observed.
- For the 5 ions for which no database entry was found, the sequence produced by the *de novo* algorithm was searched using blastp (www.ncbi.nlm.nih.gov) against nrDB (www.ncbi.nlm.nih.gov), and each was matched to one of the two protein components of our mixture - Pyruvate Kinase or Fructose-bisphosphate Aldolase A.
- blastp matches the sequence (**Figure 6**) from the first of the unmatched ions from **Table 2** to Pyruvate Kinase. This match was not made during the initial database search as the cleavage is non-tryptic:

(L)IVLTESGR

- Similarly, the second of the unmatched ions (**Figure 6**) is also matched by blastp to Pyruvate Kinase. This match was not made during the initial database search as the cleavage is again non-tryptic:

(R)PVAVA LDTK

- Similarly for **521.751...**

(M)AMGSVEASYK




Figure 5. Electrospray MS/MS spectrum for the doubly charged precursor ion m/z 434.7, annotated against the peptide sequence MQHLLAR




Figure 6. Electrospray MS/MS spectrum for the doubly charged precursor ion m/z 437.7, annotated against the peptide sequence IVLTESGR




Figure 7. Electrospray MS/MS spectrum for the doubly charged precursor ion m/z 457.2, annotated against the peptide sequence PVAVALDTK

Conclusion

The full automation of a *de novo* sequencing algorithm has been shown.

Results from database searching and *de novo* sequencing have been compared and contrasted.

This approach can be used to identify those spectra initially unmatched by a database search, due to non-tryptic cleavages of the protein.

Future work will investigate the use of an automated BLAST routine in combination with the *de novo* sequencing algorithm.

BLAST is a registered trademark of the National Library of Medicine (USA)

Maxent is a trademark of Maxent Solutions Ltd.

**Author to whom all correspondence
should be addressed:**

Alan Millar
Waters Corporation
(Micromass UK Limited)
Floats Road, Wythenshawe
Manchester, M23 9LZ
Tel: + 44 (0) 161 946 2400
Fax: + 44 (0) 161 946 2480
e-mail: alan.millar@micromass.co.uk

WATERS CORPORATION
34 Maple St.
Milford, MA 01757 U.S.A.
T: 508 478 2000
F: 508 872 1990
www.waters.com

Made in the United Kingdom

Waters
RIGHT ON TIME.

