# Waters

## PRECURSOR ION DISCOVERY - GLYCOPEPTIDE DETECTION AND CHARACTERISATION FROM A SINGLE CAPLC Q-TOF EXPERIMENT

Mark A. Ritchie, Neil Johnson, John B. Hoyes, Alan Millar, Robin Carruthers, Chris Jones and James I. Langridge

Waters Corporation, Manchester, UK.

### Presented at ASMS ASILOMAR, USA, October, 2002

## OVERVIEW

- Parent ion discovery on a Micromass Q-Tof Ultima applied to glycopeptide analysis
- Detection and characterisation of glycopeptides from CapLCMS-MS/MS run of a tryptic digestion of fetuin
- Glycopeptide MS/MS data analysis using MaxEnt3 and CarboTools

### **INTRODUCTION**

The study of post-translational modifications (PTMs) is of ever increasing importance in the study of proteomics. Physiological changes in an organism can often leave their mark upon proteins in the form of altered PTMs most notably changes in phosphorylation and glycosylation.

Analytical approaches to glycosylation analysis often assume a "global" approach; glycans are released from the glycoprotein either enzymatically or chemically. This gives no information of glycosylation site occupancy, and will not distinguish the source of the glycan in a protein/glycoprotein mixture <sup>(1)</sup>.

In a proteomic experiment proteins are often identified by performing tryptic digests of either protein solutions or excised gel bands and analysing by electrospray (ESI) LC-MS/MS. Identification and characterisation of glycosylation in this same experiment would save the use of expensive exoglycosidases, glycan release protocols and would provide site specific information in conjunction with routine peptide analysis. N-linked glycopeptides are often large, the carbohydrate portion being typically 2-3KDa or more. Often these large and mobile moieties hinder protease hydrolysis, hence trypsin or other proteases may miss a cleavage site close to a glycosylated asparagine, resulting in an unusually large peptide portion <sup>(2)</sup>. Large species are only detected in the mass range as highly charged ions, hence 4,5 or 6+ ions are characteristic of glycopeptides <sup>(3)</sup>.

Parent ion discovery (PID) is a function of the Micromass Q-Tof (Waters Corp., Milford, MA) which allows specific targeting for MS/MS of species with a common structural motif during an LC-MS experiment. It has previous been successfully applied to detection of phosphopeptides <sup>(4)</sup>.

### **EXPERIMENTAL**

Bovine fetuin (Sigma) was reduced and akylated with DTT and IAA respectively, and digested with trypsin before being analysed by ESI LC MS and MS/MS. Samples were injected using an autosampler onto a C18 cartridge pre column and on-line separation was carried out on a 75µm column (C18 PepMap, LCPackings, Amsterdam, Netherlands). A gradient from 0-40% B (B=95% ACN 0.1% FA) over 40 minutes was provided using a Micromass CapLC.

The Micromass Q-Tof Ultima API mass spectrometer with a nano Z-spray source was run in Parent Ion Discovery (PID) mode, where the voltage on the gas collision cell is switched alternately between high (30v) and Iow (8v) every second <sup>(4)</sup>. This provides both a standard precursor MS spectrum and a spectrum of all the product ions seen in the standard scan.



# Waters

Upon detection of carbohydrate b ions (m/z 204, 366, 274, 292) the instrument switches to MS/MS mode and selects the most intense multiply charged ion for fragmentation. Only if the selected precursor ion yields the expected carbohydrate b type product ions will the MS/MS continue (scheme 1). The software was modified to specifically allow the selection of ions with 4 + to 7+ charges for MS/MS. A collision energy ramp of 20-40v was applied to obtain a diverse fragmentation spectra.



Scheme 1. Flow Chart of Glycopeptide PID experiment

For comparative purposes, the experiment was repeated but the instrument was prevented from switching into MS/MS. The low energy survey served as a view of the chromatography, while a reconstructed mass chromatogram for glycopeptide B ions gives an "pseudo-SIM" chromatogram for glycopeptide detection.

## RESULTS

### Data acquisition

The glycopeptides eluted during the last 10 minutes of the gradient, as seen on the pseudo-SIM scan (**figure 1, a + b**). From reconstructed mass chromatograms for glycopeptide B ions of the 4 MS/MS functions (**figure 1, c-f**) several glycopeptides appear to have been selected for MS/MS. The distinctive shape of the MS/MS chromatograms is attributed to the collision energy ramp applied.



Figure 1. Glycopeptide PID experiment chromatograms. Fetuin tryptic digest run by CapLCMS-MS/MS, in (a) and (b) the instrument was prevented from switching into MS/MS. (a) Low energy survey scan (b) Psuedo-SIM scan, reconstructed high energy mass chromatogram for carbohydrate b ions (m/z 204, 366, 274, 292). (c-f) MS/MS chromatograms reconstructed as mass chromatograms for carbohydrate b ions.

## Waters

#### Data processing and interpretation

The charge state and the m/z of the ion can be used to assign a mass to the glycopeptide, from which a putative carbohydrate composition may be assigned using Glycomod from the Expasy website (http://ca.expasy.org)

Glycopeptides were then characterised from the associated MS/MS data using MaxEnt3 and CarboTools in Masslynx 4.0 software (Waters Corp., Milford, MA). The MaxEnt3 de-iosotoped and deconvoluted MS/MS spectra were imported into CarboTools. The smallest y-type fragment ion (nomenclature as Domon and Costello <sup>(5)</sup>) detected is the usually the glycopeptide y1 ion, (the peptide with the reducing terminus N-acetyl hexosamine attached), Figure 2. Using Carbotools and working up from this intense ion a structure for the carbohydrate moiety on the glycopeptide can be established. In the example shown here, a tetrasialylated triantennary is characterised using CarboTools (Figure 3). In the case of insufficient MS/MS data for carbohydrate characterisation, glycan structure and peptides were assigned from Glycomod.

The glycopeptides characterised were largely trisialylated triantennaries, as expected <sup>(6)</sup>. Some tetra-sialylated triantennaries were also found.

A number of sodiated glycopeptide species were selected for MS/MS, but have been omitted from the summary in **table 1**.



Figure 2. MaxEnt3 deconvoluted deisotoped spectrum of a fetuin tryptic peptide bearing a trisialylated triantennary glycan. Fragmentation is preferential to the glycan, and no fragmentation of the peptide is observed. The y1 ion (peptide + HexNAc) at m/z 3758 is intense and can be used to easily determine the peptide mass.



Figure 3. MaxEnt3 data analysed using CarboTools. Carbohydrate differences between the fragment ions are measured and assigned manually with use of the software. Here a tetra-sialylated triantennary is sequenced entirely.

Glycopeptide we MW (Da)	Precursor Ion (m/z)	Peptide	Carbohydrate
6474.91	1296.66, 5+ 1619.73, 4+	72RPTGEVYDIEIDTLETTCHVLDPT PLANCSVR103 ICMC	
6419.00	1605.75, 4+	72RPTGEVYDIEIDTLETTCHVLDPT PLANCSVR103 No CMC	
5816.42	1459.23, 4+	72RPTGEVYDIEIDTLETTCHVLDPT PLANCSVR103 1CMC	┝⊒╼╺╱ ┝⊒╼┙
5760.91	1441.23, 4+	72RPTGEVYDIEIDTLETTCHVLDPT PLANCSVR103 No CMC	┝⊒╼╺╱ ┝⊒╼┙
6709.73	1342.95, 5+	72RPTGEVYDIEIDTLETTCHVLDPT PLANCSVR103 No CMC	
6533.67	1307.89, 5+ 1634.23, 4+	«VWPRRPTGEVYDIEIDTLET TCHVLDPTPLANCSVR103 2 CMC	
3904.61	1302.53, 3+	145LCPDCPLLAPLNDSR159 2CMC	н-
4033.032	1343.34, 3+	145LCPDCPLLAPLNDSR159 No CMC	
4486.93	1496.53, 4+	145LCPDCPLLAPL <u>NDS</u> R <sub>159</sub> No CMC	
7270.52	1455.10, 5+	143LCPDCPLLAPLNDSRVVHAV EVALATFNAESNGSYLQLVE ISR187 No CMC	
6167.90	1542.97, 4+	160VVHAVEVALATFNAESNGSY LQLVEISR187 No CMC	
5877.71	1470.43, 4+	160VVHAVEVALATFNAESNGSY LQLVEISR187 No CMC	

Vaters

Table 1. Tryptic glycopeptides of fetuin detected and characterised with the glycopeptide PID experiment. Where MS/MS data was insufficient for characterisation, most likely Glycomod assignments have been chosen.

 identification done on Glycomod alone, insufficient MS/MS data CMC - cysteines derivatised with IAA.

## CONCLUSION

From the fetuin digest a total of 12 glycopeptides were identified from all three N-linked glycosylation sites. Peptide heterogeneity from both derivatisation and cleavage location were observed, giving rise to a greater number of glycopeptides than expected. The use of an exclude list could possibly increase this number of glycopeptides detected after a subsequent run.

The method described is a sensitive, highly specific approach to glycopeptide detection and sequencing. In combination with other methods, this approach will allow full elucidation of N-linked glycosylaytion structure, including linkage information.

## REFERENCES

- (1) Harvey, D. J., Proteomics (2001) 1, 311-328
- Stimpson, E., Hope, J., Chong, A. L., Burlingame, A. L., *Biochemistry* (1999) 38, 4885-4895
- Ritchie, M. A., Deery, M. J., Lilley K., Gill, A.
  C., J Am Soc Mass Spectrom (2002) 13 (9): 1065-77
- (4) Langridge, J. I., Hoyes J. B., Bateman, R. H., Millar, A. L., Carruthers R., Jones C., Jensen O, N., poster at 49<sup>th</sup> ASMS, 2001.
- (5) Domon, B., Costello E. C., *Glycoconjugate J* (1988) 5:397-409
- Medzihradsky, K. F., Maltby, D. A., Hall, S. C., Settinerri, C. A., Burlingame, A. L., Am Soc Mass Spectrom. (1994) 5:350-358

Contact details for author: mark.richie@micromass.co.uk

#### **MS Sales Offices:**

BELGIUM 02-2534550 CANADA 514 694-1200 DENMARK 4657 4101 EUROPE +31 (0) 36-540 6000 FINLAND 02 284 56 11 FRANCE 0800-907016 GERMANY 0800-1817249 ITALY 02 2159 1415 NETHERLANDS 036-540 6160 NORDIC +46 (0) 8 555 115 10 SPAIN 93 440 71 30 SWEDEN 08 555 115 10 SWITZERLAND 62 88 92 030 UK 0161 435 4125 USA 978 524-8200

Author to whom all correspondence should be addressed: Mark A. Richie Waters Corporation (Micromass UK Limited) Floats Road, Wythenshawe Manchester, M23 9LZ Tel: + 44 (0) 161 946 2400 Fax: + 44 (0) 161 946 2480 e-mail: mark.richie@micromass.co.uk

WATERS CORPORATION 34 Maple St. Milford, MA 01757 U.S.A. T: 508 478 2000 F: 508 872 1990 www.waters.com

Waters RIGHT ON TIME.



Made in the United Kingdom

©2002 Waters Corporation October 2002 / WMP255 For research use only. Not for use in diagnostic procedures.



Certificate No: 951387