

Analysis of Steroid Isomers by Capillary Electrochromatography

Gordon Ross and
Rainer Schuster

Pharmaceutical

Abstract

Capillary electrochromatography (CEC) is a fusion of capillary electrophoresis and liquid chromatography which preserves and exploits the best aspects of both techniques. The separation is achieved for neutral molecules essentially through chromatographic selectivity mechanisms, while the mobile phase is propelled through the packed capillary by electromotive forces. Efficiencies can be increased by as much as ten-fold, compared with conventional LC systems. Such increases in efficiency can allow faster method development times and also separations of very closely related compounds. This application brief contrasts the first-pass separation of α and β isomers of 17-hydroxyoestradiol using HPLC and CEC. These isomers differ only by the position of the 17-OH group (figure 1).




Figure 2
Separation of alpha and beta 17-OH estradiol by HPLC




Figure 1
Structure of 17-OH estradiol

Conditions

Column 3 μ m ODS 125 mm \times 2 mm i.d.

Mobile Phase

NH₄Ac pH 7.0/acetonitrile (50/50)

Flow 0.26 ml/min

Temp 38 $^{\circ}$ C



Agilent Technologies

Innovating the HP Way




Figure 3
Separation of alpha and beta 17-OH estradiol isomers by CEC

Experimental

The CEC separation was performed using the Agilent Capillary Electrophoresis system equipped with diode array detector and computer controlled via Agilent ChemStation software. CEC capillary columns were Agilent CEC C18. The Agilent CE system is uniquely designed for operating CEC in that it can apply up to 12 bar pressure simultaneously to both vials in order to suppress bubble formation and maximize reliability.

HPLC of the isomers was carried out using a HP 1090 Series II liquid chromatograph equipped with diode-array detection and computer controlled via Agilent ChemStation software.


Figure 2 shows the separation obtained from HPLC analysis of these compounds. However in obtaining a rapid analysis the compounds are only partially resolved. When using CEC (figure 3) for separation, it should be noted that the retention time is markedly increased. This is due to the difference in stationary phase since the CEC phase has been specifically designed to meet the requirements of CEC. The two isomers are however well resolved and fully capable of quantitation.

Conditions

Column Agilent 3 μ m CEC C18, 250 mm (335 mm) \times 0.1 mm id
Mobile Phase 25 mM TRIS/acetonitrile (5/95)
Voltage 25 kV
Temperature 30 $^{\circ}$ C
Injection 10 s at 10 kV
Pressure 10 bar both sides

Equipment

- Agilent Capillary Electrophoresis System
- Agilent 1090 Series II LC
- Agilent ChemStation + software



Gordon Ross and Rainer Schuster are application chemists at Agilent Technologies, Waldbronn, Germany.

For more information on our products and services, visit our worldwide website at <http://www.agilent.com/chem>

© Copyright 1997 Agilent Technologies
 Released 06/97
 Publication Number 5965-9030E



Agilent Technologies
 Innovating the HP Way