

Purity Check of a Synthetic α -homo Polylysine (250) Preparation using Capillary Electrophoresis

Rudolf Grimm

Bioscience

Abstract

The smaller synthesis by-products can be resolved from the main component. The purity level of the 250-mer was estimated to be > 99.7 %.




Figure 1
Purity check of a synthetic α -homo polylysine (250)
preparation using capillary electrophoresis

Conditions

Buffer

50 mM phosphate, pH 2.0,

Sample

synthetic polylysine 250-mer

Capillary

effective length 72 cm

total length 80.5 cm

internal diameter 50 μ m

internal diameter at point of
detection is 150 μ m

Injection

200 mbars

Temperature

25 $^{\circ}$ C

Field strength

370 V/cm




Agilent Technologies

Innovating the HP Way

Equipment

Agilent Capillary Electrophoresis system



Rudolf Grimm is application chemist at Agilent Technologies, Waldbronn, Germany.

For more information on our products and services, visit our worldwide website at <http://www.agilent.com/chem>

© Copyright 1997 Agilent Technologies
Released 12/97
Publication Number 5966-2956E



Agilent Technologies
Innovating the HP Way