Analyzing Compounds of Environmental Interest Using an LC/Q-TOF Part 2: Fluorotelomer Unsaturated Acids Application

Environmental

Authors

Jim Lau, Chin-Kai Meng, and Jennifer Gushue Agilent Technologies, Inc. 2850 Centerville Road Wilmington, DE 19808 USA

Robert J. Letcher and Shaogang Chu Wildlife and Landscape Science Directorate Science and Technology Branch Environment Canada National Wildlife Research Centre Carleton University Ottawa, Ontario, K1A 0H3 Canada

Abstract

Perfluoroalkyl substances (PFASs) have been widely used in a variety of products due to their chemical inertness, resistance to heat, and ability to repel water and oils. PFASs exhibit a high propensity for persistence and bioaccumulation in wildlife, which is causing concern. Due to different manufacturing and degradation processes, different chain lengths and functional groups of PFAS exist in the environment. An LC/Q-TOF is best suited to screen, identify, and quantify many perfluoroalkyl compounds in different matrices.

This LC/Q-TOF application shows good mass accuracy (well below 3 ppm) and good resolution (> 13,500). Excellent quantification results for selected PFASs were obtained from extracts derived from blind liver samples.

Introduction

Environment Canada is tasked with risk assessment and the evaluation of impact of a variety of compounds in environmental matrices, including wildlife tissues, water, sediment, and air [1, 2]. A number of compounds classified as perfluoroalkyl substances (PFASs) have been widely used in a variety of products due to their chemical inertness, resistance to heat, and ability to repel water and oils. Some of the commercial PFAS products include lubricants, adhesives, stain and soil repellents, paper coatings, and fire-fighting foams. Due to their unique chemical and biological stability, some PFASs exhibit a high propensity for persistence and bioaccumulation in wildlife. In recent years, particular environmental concern has arisen as a number of PFASs have been reported in tissues of marine mammal, seabird, and fish species inhabiting various regions of the Arctic. More specifically, some PFASs, such as perfluorooctane sulfonate (PFOS) [CF₃(CF₂)₇SO₃H] and C8 to C15 chain length perfluorinated carboxylates (PFCAs), have been reported present at similar or higher concentrations than persistent organochlorines (OCs) in polar bears [3-5].

Many of the PFASs originate from the two manufacturing processes of electrochemical fluorination (perfluoroalkyl sulfonamido alcohols degrade to PFOS) and telomerization (fluorotelomer alcohols [FTOHs] are transformed to PFOA). The degradation pathway for telomers is:



From 8:2 fluorotelomer alcohol (8:2 FTOH) to:

Major products: 8:2 fluorotelomer aldehyde

(8:2 FTAL), 8:2 fluorotelomer carboxylate (8:2 FTCA), and perfluorooctanoic acid (PFOA) $[CF_3(CF_2)_6CO_2H]$

Minor products: 8:2 fluorotelomer unsatu-

> rated carboxylate (8:2 FTUCA, CF₃(CF₂)₆CF= COOH) and perfluorononanoic acid (PFNA)

Then from 8:2 FTCA and 8:2 FTUCA (degradation products from above) to:

• Major product: **PFOA** • Minor product: **PFNA**

This application uses fluorotelomer unsaturated carboxylates (FTUCAs - metabolites of telomer alcohols) to demonstrate the mass accuracy, resolution, and the quantification capability of the LC/Q-TOF.

Experimental

Samples

This study analyzed three standard solutions of the following three target compounds at 1, 100, and 250 ppb. Three internal standards (ISTDs) at 50 ppb each were also added to each of the three solutions. A fourth sample was a blind liver extract with an undisclosed concentration of the three target compounds.

Target Compounds:

FHUEA (C8, 6:2 FTUCA) 2H-Perfluoro-2-octenoic acid

CF₃(CF₂)₄CF=CH-COOH

FOUEA (C10, 8:2 FTUCA) 2H-Perfluoro-2-decenoic acid

CF3(CF2)6CF=CH-COOH

FDUEA (C12, 10:2 FTUCA) 2H-Perfluoro-2-dodecenoic acid

CF₃(CF₂)₈CF=CH-COOH

ISTDs:

 $^{2H\text{-Perfluoro-}[1,2-^{13}C_2]\text{-}2\text{-}octenoic}$ acid $\text{CF}_3(\text{CF}_2)_4\text{CF=}\text{C}^{13}\text{H-C}^{13}\text{OOH}$ **MFHUEA**

 $2 H\text{-Perfluoro-} [1,2^{-13}C_2]\text{-}2\text{-decenoic acid}$ **MFOUEA**

 $CF_3(CF_2)_6CF = C^{13}H - C^{13}OOH$

2H-Perfluoro-[1,2-¹³C₂]-2-dodecenoic acid **MFDUEA**

CF₃(CF₂)₈CF=C¹³H-C¹³OOH

Instrument Parameters

All sample analyses were performed on an Agilent 1200 SL Rapid Resolution LC coupled to an Agilent 6520 Q-TOF.

All sample analyses were performed under LC/ Q-TOF autotune conditions. Mass accuracy, sensitivity, and resolution for all samples were measured without any changes to 6520 Q-TOF instrument parameters, except ion source conditions appropriate for the spray chamber type, LC flow, and sample thermal stability.

3 mM NH₄OAc Mobile A

Mobile B Me0H

LC column ZORBAX XDB 2.1 mm × 50 mm,

C-18, 3.5-um particle size

Flow rate 0.5 mL/min Injection volume 10 µL

Scanned at 2.5 scans/sec, 50 to MS

 $1.100 \, m/z$

m/z 113, 1034 Negative ref. ions Q-TOF parameters Set by autotune Drying gas 12 L/min N2 at 300 °C

Nebulizer pressure 50 psi ESI (-) 3 KV Fragmentor 275 V

Results and Discussion

Table 1 shows the mass accuracy and resolution of the six compounds analyzed in this study. The accuracy for each compound was under 3 ppm and the resolution was > 13,500. Table 2 shows the quantitation ions and results. The quantitation ion was not the molecular ion but a fragment ion from the collision-induced dissociation (CID). The fragment ion is a lost fragment of "FCOOH" ([M-H] - 64) from the deprotonated molecular ion. The calibration range was 1 to 250 ppb (that is, 10 pg to 2.5 ng on column). Good signal-to-noise ratios (un-smoothed signal) were observed for the three FTUCA standards at 1 ppb (see Figure 1).

A blind mixture of the three standards in liver extract was quantified by Q-TOF. Excellent quantification results were obtained for all three targets:

| Compound | Measured | Actual |
|----------|----------|--------|
| FOUEA | 4.0 ppb | 5 ppb |
| FHUEA | 2.7 ppb | 5 ppb |
| FDUEA | 6.7 ppb | 5 ppb |

By comparing the ISTD responses, ion suppression (matrix effect) was observed in analyzing the liver extract. It has been shown using an LC/QQQ system with atmospheric pressure photoionization (APPI) is a more effective technique than electrospray ionization (ESI) to analyze FTOHs and perfluorinated sulfonamides [6].

Table 1. Mass Accuracy and Resolution of the Six FTUCA Analyzed in This Study

| | Formula | Mass accuracy (ppm) | Resolution |
|--------|--------------------------------|------------------------|------------|
| FHUEA | $C_8H_2F_{12}O_2$ | -0.98 | 13692 |
| MFHUEA | $[^{13}C]_2C_6H_2F_{12}O_2$ | - 0.39 | 13789 |
| FOUEA | $C_{10}H_2F_{16}O_2$ | - 0.19 | 14018 |
| MFOUEA | $[^{13}C]_2C_8H_2F_{16}O_2$ | - 2.08 | 13564 |
| FDUEA | $C_{12}H_2F_{20}O_2$ | - 0.62 | 15399 |
| MFDUEA | $[^{13}C]_2C_{10}H_2F_{20}O_2$ | – 1.33 | 14322 |

Table 2. Quantification lons and Linearity of the FTUCA
(The quantitation ion is a lost fragment of "FCOOH" from the molecular ion.)

| • | • | | | , | |
|--------|--|--------------------|--|-------------------------------|--|
| | Formula | [M–H] [–] | Quant lon ([M–H] [–] – 64) | R ² (1–250 ppb) | |
| FHUEA | C ₈ H ₂ F ₁₂ O ₂ | 357 | 293 | 0.9992 | |
| MFHUEA | $[^{13}C]_2C_6H_2F_{12}O_2$ | 359 | 294 | - | |
| FOUEA | $C_{10}H_2F_{16}O_2$ | 457 | 393 | 0.9998 | |
| MFOUEA | $[^{13}C]_2C_8H_2F_{16}O_2$ | 459 | 394 | - | |
| FDUEA | $C_{12}H_2F_{20}O_2$ | 557 | 493 | 0.9972 | |
| MFDUEA | $[^{13}C]_2C_{10}H_2F_{20}O_2$ | 559 | 494 | _ | |
| | | | | | |




Figure 1. Expansion of EICs for quant ions of 1 ppb FTUCA standards (S/N > 5, unsmoothed signal).

Conclusions

Some of the most impressive aspects about the LC/Q-TOF were the mass accuracy, the linear dynamic range, and the signal-to-noise ratios. For environmental applications, using Q-TOF can greatly reduce interference from complex matrices and improve the accuracy of the results. This application of Q-TOF shows good mass accuracy (well below 3 ppm) and good resolution (> 13,500). Excellent quantification results were obtained from a blind liver extract.

References

- Jim Lau, Chin-Kai Meng, Jennifer Gushue, Mehran Alaee, and Vimal Balakrishnan, "Analyzing Compounds of Environmental Interest Using an LC/Q-TOF – Part 1: Dyes and Pigments," Agilent Technologies publication 5989-9336EN, August 2008.
- Jim Lau, Chin-Kai Meng, Jennifer Gushue, Mark Hewitt, and Suzanne Batchelor, "Analyzing Compounds of Environmental Interest Using an LC/Q-TOF – Part 3: Imidacloprid and Manool," Agilent Technologies publication 5989-9129EN, August 2008.
- J. W. Martin, M. M. Smithwick, B. M. Braune, P. F. Hoekstra, D. C. G. Muir, and S. A. Mabury, "Identification of Long-Chain Perfluorinated Acids in Biota from the Canadian Arctic," *Environ. Sci. Technol.*, 2004, 38, 373–380.

www.agilent.com/chem

- M. M. Smithwick, D. C. G. Muir, S. A. Mabury, K. R. Solomon, J. W. Martin, C. Sonne, E. W. Born, R. J. Letcher, and R. Dietz, "Perfluoroalkyl Contaminants in Liver Tissue from East Greenland Polar Bears (*Ursus maritimus*)," *Environ. Toxicol. Chem.*, 2005, 24, 981–986.
- M. M. Smithwick, S. A. Mabury, K. R. Solomon, C. Sonne, J. W. Martin, E. W. Born, R. Dietz, A. E. Derocher, R. J. Letcher, T. J. Evans, G. W. Gabrielsen, J. Nagy, I. Stirling, M. K. Taylor, and D. C. G. Muir, "Circumpolar Study of Perfluoroalkyl Contaminants in Polar Bears (Ursus maritimus)," Environ. Sci. Technol., 2005, 39, 5517–5523.
- 6. S. G. Chu and R. J. Letcher, "Methodological Comparison of the Determination of Fluorotelomer Alcohols and Perfluorinated Sulfonamides in Biotic Samples by Liquid Chromatography-Atmospheric Pressure Photoionization and -Electrospray Mass Spectrometry," in preparation.

For More Information

For more information on our products and services, visit our Web site at www.agilent.com/chem.

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc. 2008

Printed in the USA August 11, 2008 5989-9132EN

