

# Software-assisted, high-throughput identification of main metabolites of pharmaceutical drugs

Rapid data acquisition by Agilent 1290 Infinity LC, TOF and Q-TOF instrumentation, and subsequent identification of metabolites by Agilent MassHunter Metabolite Identification software

# **Application Note**

Metabolite identification in drug discovery and drug development

# Author

Edgar Naegele Agilent Technologies Waldbronn, Germany



# Abstract

This Application Note describes:

- Rapid separation of metabolites generated from in-vitro experiments using the Agilent 1290 Infinity LC, system
- Fast acquisition of TOF mass spectra using Agilent 6530 Accurate-Mass Quadrupole Time-of-Flight LC/MS systems
- Fast, software-assisted identification of main metabolites from in-vitro experiments using Agilent MassHunter Metabolite Identification software
- Generation of reports for the identified metabolites using Agilent MassHunter software



# **Introduction**

In modern pharmaceutical drug development it is of crucial importance to analyze the adsorption, distribution, metabolism and excretion (ADME) properties of possible new drug candidates as quickly as possible in order to make decisions about further investments in the development of a special compound. To find compounds with the correct properties it is essential to screen a large number of compounds for their ADME properties, which requires to work in an high-throughput environment. This Application Note describes the application of the Agilent 1290 Infinity LC system, the Agilent 6530 Q-TOF MS system and the MassHunter Metabolite Identification software for fast, highthroughput identification of main metabolites of new pharmaceutical drug candidate compounds.

# **Experimental**

# Equipment

- Agilent 1290 Infinity LC system consisting of 1290 Infinity Binary Pump with integrated degasser, 1290 High Performance Autosampler with thermostat, and 1290 Infinity Thermostatted Column compartment
- Agilent 6530 Accurate-Mass Q-TOF LC/MS system
- Agilent MassHunter Metabolite Identification (MetID) software
- Column: ZORBAX SB-C18, 2.1 x 50 mm, 1.8 μm

# Sample preparation

The following stock solutions were used:

- 20 mg/mL microsomal S9 preparation
- 0.1 mg/mL buspirone in water
- 1.6 mg NADP in 1.6 mL 0.1 M phosphate buffer, pH 7.4

- 50 mM isocitrate/MgCl<sub>2</sub> (203 mg MgCl<sub>2</sub>.6H<sub>2</sub>O + 258.1 mg isocitrate in 20 mL H<sub>2</sub>O)
- Isocitrate dehydrogenase 0.33 unit/µL

NADPH regeneration system: 1.6 mL NADP solution + 1.6 mL Isocitrate solution + 100 µL IDH solution.

Incubation mixture: 3.85 µL substrate + 200 µL NADPH regeneration system + 746.15 µL phosphate buffer + 50 µL S9.

Incubation was carried out at 37 °C for 60 minutes. A 100  $\mu$ L aliquot was taken at the beginning (t=0) and at t=60 min. The reaction was stopped by adding 6  $\mu$ L perchloric acid and 100  $\mu$ L acetonitrile followed by centrifugation for 15 min at 14,000 rpm. The supernatant was evaporated to dryness using a SpeedVac concentrator and reconstituted with water containing 0.1 % formic acid for LC/MS analysis. The incubation sample stopped at 0 min was used as control.

### LC method

Solvent A:	Water + 0.1 % formic acid
Solvent B:	ACN + 0.1 % formic acid
Flow:	0.8 mL/min
Gradient	0 min, 5 %B; 0.10 min,
	5 %B; 1.10 min, 75 %B;
Stop time:	1.1.0 min
Post time:	1 min.
Injection:	Volume 5 µL, sample
	cooler at 4 °C, needle wash
	in 50 % methanol for 5 s,
	injection loop to bypass
	at 0.1 min with flush out
	factor 16
Column:	Temperature 60 °C

# **TOF MS** method

Source:	ESI positive
Capillary:	3500 V
Dry gas:	12 L/min
Nebulizer:	55 psi
Gas temp.:	350 °C

Skimmer:65 VFragmentor:200 VMass range:100-1000 m/zAcquisitionrate:rate:5 spectra/sReference121.0508 and 922.0080

# Data analysis method in the MetID software

The first step in the analysis comprised a comparison between the data file that contained the metabolite compounds (metabolite sample) and the data file that contained only the parent drug (control sample). All detectable mass signals were extracted from the MS level data using the Molecular Feature Extraction (MFE) algorithm. Related compound isotope masses and adduct masses were grouped together into discrete molecular features, and chemical noise was removed. The compounds lists of the metabolized sample and the control were then compared.

All new compounds or those that increased twofold in the metabolized sample were considered potential metabolites and were subjected to further analysis by different algorithms. The algorithms can identify and qualify new metabolites, or just qualify metabolites found by another algorithm. In this high-throughput experiment all algorithms' results were weighted equally and combined into a final identification relevance score. Metabolites were qualified when their final score was above the stringently defined relevance threshold. The results from all algorithms were collated in a results table, which could be inspected at-a-glance and reported<sup>1</sup>.

# **Results and discussion**

To achieve fast separation of the metabolites on a 50 mm, 1.8 µm particle size column, a 1 minute gradient was applied by the Agilent 1290 Infinity LC system. The metabolites were generated from the pharmaceutical test compound buspirone in an in-vitro assay. For adequate detection with the time-of-flight mass spectrometer the instrument was operated at a data rate of 5 Hz.

After generation the data was loaded into the MetID software and analyzed using a common method. The result was displayed by the MetID software in an at-a-glance table, in which the result for each metabolite could be examined in more detail (figure 1). From the results table a summary report was generated, which showed the available information for each metabolite (figure 2). The more extensive report contained the detailed results for each metabolite. As example the result for a mono-hydroxyl metabolite (figures 3 to 5) and a dihydroxy metabolite (figures 6 to 8) of buspirone are discussed here.



#### Figure 1

Result table showing an at-a-glance summary of buspirone metabolite analysis with overall identified metabolites, extracted ion chromatograms (EIC), extracted compound chromatograms (ECC), isotopic pattern analysis and calculated formulas.

Name	Mass	RT	Rel.	Qual.	User	SC	IPM	EIC	MDF	Form.	BioXF
2x Hydroxylation	417.2379	0.59	100.00	$\checkmark$	1	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	1	1
Hydroxylation	401.2423	0.63	100.00	$\checkmark$	1	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	1	1
Hydroxylation	401.2424	0.66	100.00	$\checkmark$	1						
2x Hydroxylation	417.2388	0.72	100.00	$\checkmark$							
Hydroxylation	401.2439	0.75	100.00	$\checkmark$							
Hydroxylation	401.2430	0.79	100.00	$\checkmark$							
Buspirone	385.2478	0.82	_		_	_	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	_
Hydroxylation	401.2429	0.84	75.00	×	1	×	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	1

#### Figure 2

Summary result report, including qualified metabolites sorted by their retention times (RT), with their metabolite names and relative score, molecular mass and the passed flag for individual algorithm results. SC=Sample-control comparison, IPM = Isotopic Pattern Matching, EIC = Extracted Ion Chromatogram, MDF = Mass Defect Filter, Form. = Calculated Formula, BioXF = Assigned Biotransformation, Qual. = Qualified by Score, User = Qualified by User.

The extensive report for the monohydroxyl metabolite, which eluted after 0.75 minutes at m/z 402.2511, showed the detailed information about the metabolite itself such as measured accurate mass, calculated formula, assigned biotransformation and ion species. Further, the report showed more detailed information about the result of each individual algorithm, for example, Molecular Feature Extraction (MFE), Extracted Ion Chromatogram (EIC) compound search and Mass Defect Filter Result (figure 3). For the hydroxyl metabolite the possible formula was calculated based not only on a defined mass error window but also on the measured isotopic pattern, which increased the quality of the calculated formula and limited the possible number of hits significantly. These results were also displayed in the detailed metabolite result report for the formula (figure 4).

Name Hy	droxylation	Bio	XF Name	Hydro	oxylation				
Formula C2	1H31N5O3	Mas	55	401.2	439				
<b>m/z</b> 40	2.2511	Spe	cies	(M+H	i)+				
RT 0.3	/54	San	nple Type	e Metab	boliteSample				
MFE Com	pound Search								
Mass	m/z	Spe	cies	RT	Start Time	End Time	Volume	Height	
401.2439	402.2511	(M+	H)+	0.754	0.739	0.774	192448	187344	
EIC Com	oound Search								
Mass	m/z	Spe	cies	RT	Start Time	End Time	Area	Area %	
401.2427	402.2500	(M+	H)+	0.755	0.739	0.774	149323	100.00	
Sample (	Comparison Res	sults							
<u>Sample (</u> Qualified ☑	Comparison Res Changed New	Resp. Ra	atio Cor	7. RT 1	Normalized Heig	ht			
Sample C Qualified ⊠ Isotopic	Comparison Res Changed New Pattern Matchi	Resp. Ra na Resu	atio Cor <i>Its</i>	r.RT I	Normalized Heig	ht			
<u>Sample (</u> Qualified ⊠ <u>Isotopic</u> Qualified	Comparison Res Changed New Pattern Matchi Score	Resp. Ra n <u>a Resu</u> Delta m	atio Cor <u>/ts</u> /z	r. RT 1	Normalized Heig	ht			
<u>Sample C</u> Qualified <u>∅</u> <u>Isotopic</u> Qualified ☑	Comparison Res Changed New Pattern Matchin Score 95.91	Resp. Ra n <u>a Resu</u> Delta m 0.00	atio Cor <u>Its</u> /z	r. RT I	Normalized Heig	ht			
<u>Sample C</u> Qualified ⊠ <u>Isotopic</u> Qualified ☑ Mass Def	Comparison Res Changed New Pattern Matchin Score 95,91 Sect Filter Result	Resp. Ra Resp. Ra <u>na Resu</u> Delta m, 0.00 <u>lts</u>	atio Cor <u>Its</u> /z	r. RT I	Normalized Heig	ht			
<u>Sample C</u> Qualified ⊠ <u>Isotopic</u> Qualified ⊠ <u>Mass Def</u> Qualified	Comparison Res Changed New Pattern Matchi Score 95.91 Fect Filter Resul Delta Mass [mDa	Resp. Ra Resp. Ra <u>na Resu</u> Delta m, 0.00 <u>lts</u> a]	atio Cor <u>Its</u> /z	r. RT I	Normalized Heig	ht			
Sample C Qualified ⊠ Isotopic Qualified ⊠ Mass Def Qualified ⊠	Changed New Pattern Matchi Score 95.91 Fitter Result Delta Mass [mDa -3.91	Resp. Ra n <u>a Resu</u> Delta m, 0.00 <u>Its</u> a]	atio Cor <u>Its</u> /z	r. RT I	Normalized Heig	ht			
Sample C Qualified ☑ Isotopic Qualified ☑ Mass Def Qualified ☑ Formula	Comparison Res Changed New Pattern Matchin Score 95,91 'ect Filter Result -3.91 Results	Resp. Ra n <u>a Resu</u> Delta m, 0.00 <u>lts</u> a]	atio Cor <u>Its</u> /z	r. RT	Normalized Heig	ht			
Sample C Qualified ☑ Isotopic Qualified ☑ Mass Def Qualified ☑ Formula	Changed New Pattern Matchin Score 95.91 <u>ect Filter Resul</u> Delta Mass [mDi -3.91 <u>Results</u> Neutral Formula	Resp. Ra Resp. Ra Delta m, 0.00 ( <u>ts</u> a]	atio Cor <u>Its</u> /z Calc. Ma	r. RT I	Normalized Heig Moss [mDa]	ht Delta M	fass [ppm]	Calcu	lation Base
Sample ( Qualified Zoualified Qualified Mass Def Qualified Formula Assigned	Changed New Pattern Matchin Score 95.91 Gett Filter Resul Delta Mass [mDi -3.91 Results Neutral Formula C21H31NSO3	Resp. Ra n <u>a Resu</u> Delta m, 0.00 I <u>ts</u> a]	atio Cor / <u>Its</u> /z Calc. Ma 401.2427	r. RT I nss Delta 7 -1.17	Normalized Heig Mass [mDa]	ht Delta M -2.92	fass [ppm]	<b>Calcu</b> MfeCo	<b>llation Base</b>
Sample ( Qualified Isotonic Qualified Mass Def Qualified Formula Assigned Biotransi	Comparison Res Changed New Pattern Matchin Score 95.91 Sect Filter Result Delta Mass [mDa -3.91 Results Neutral Formula C21H31N503 Sormation Result	Resp. R: Resp. R: Delta m, 0.00 <u>lts</u> 1 <u>ts</u>	atio Cor / <u>ts</u> /z Calc. Ma 401.2427	<b>r. RT I</b> <b>155 Delta</b> 7 -1.17	Normalized Heig a Mass [mDa]	ht Delta № -2.92	fass [ppm]	<b>Calcu</b> MfeCo	llation Base mpoundMsSpectrun
Sample C Qualified Zisotopic Qualified Mass Del Qualified Ziso Formula Assigned Biotranszi Assigned	Changed New Pattern Matchi, Score 95.91 fect Filter Result Delta Mass [mDi -3.91 Results Neutral Formula C21H31N503 formation Resul Name	Resp. Ra Resp. Ra Delta m, 0.00 ( <u>ts</u> Phase	atio Cor <u>Its</u> /z Calc. Ma 401.2427 Offset F	r. RT I ass Delta 7 -1.17 formula	Normalized Heig a Mass [mDa] Delta Mass	Delta M -2.92	łass [ppm] Delta Mass	Calcu MfeCo [ppm]	Nation Base mpoundMsSpectrum Calc. Mass

#### Figure 3

Detailed metabolite report for the buspirone hydroxy metabolite at retention time 0.75 min. This part of the report gives detailed information about the identified metabolite and the identifying algorithms. Other detailed information about formula (figure 4), chromatograms and isotopic pattern (figure 5) are also available.

Metab	olite I	nformation							
Name Formula m/z	Hydroxylat C21H31N5 402.2511	ion O3	BioXF Name Mass RT	Hydroxy 401.243 0.754	vlation 39				
Formul	a Summ	ary							
Select	ed So	ore Formula	Ion	Formul	a	Mass	Calc. M	ass	⊿ Mass [ppm]
TRUE	1	00.0 C21H31N5O3	C21	H32N5O3	3	401.2439	401.24	427	-2.92
Formula C21H31f Species (M+H)+ Formul	a (M) N5O3 la Results	Selecte TRUE m/z 402.251	<b>d</b> 1						
Ion For	mula	Score	e Mas	s ⊿N	lass [mDa]	⊿ Mass [p	pm] [	DBE	
C21H32	N5O3	100.	0 401.243	)	-1.17		2.92	9	
isotopi	C Peak In	Colo Normation			Color.				
	Abund %	Calc Abund%	o 0	m/z	Calc n	n/z ⊿m/	z [ppm]		
	100.00	100.0	0 402	.2511	402.2	500	-2.91		
	21.83	25.0	z 403	.2535	403.2	529	-1.35		
	2.59	3.0.	2 404	.2502	404.2	550	-1./0		

#### Figure 4

Detailed metabolite report about the formula including isotopic pattern, calculated for the buspirone hydroxy metabolite at retention time 0.75 min.



Figure 5

Detailed metabolite report for buspirone hydroxy metabolite at retention time 0.75 min:

A) Extracted Ion Chromatograms (EIC) of compounds with mass 402.25

B) Extracted Compound Chromatogram (ECC) of buspirone hydroxy metabolite at retention time 0.75 min

C) Measured isotopic pattern of buspirone hydroxy metabolite at retention time 0.75 min (blue lines) and caclulated isotopic pattern (CIP, green box).

Finally, the EIC, ECC and isotopic pattern were displayed (figure 5). The EIC of m/z 402.25 showed 5 peaks for possible hydroxyl metabolites of buspirone with the selected one at retention time 0.75 minutes (figure 5A). The ECC showed the extracted MFE compound for the molecular mass of 401.2439 at retention time 0.75 minutes identical to the EIC (figure 5B). The measured isotopic pattern of this compound showed an excellent fit to the calculated isotopic pattern as a basis for the formula calculation (figure 5C).

Within the same data analysis the dihydroxy metabolites at a level of two orders of magnitude below the monohydroxy metabolites were also identified. The extensive report showed detailed information about the dihydroxy metabolite, which elutes after 0.71 minutes at m/z 418.2461 and the detailed information about each algorithm (figure 6).

Metabolite	Information						
Name 2x	Hydroxylation	BioXF Na	ame 2x	Hydroxylation			
Formula C21	1H31N5O4	Mass	41	7.2388			
<b>m/z</b> 418	3.2461	Species	(M	1+H)+			
RT 0.7	16	Sample	Type Me	etaboliteSample			
MFE Comp	ound Search						
Mass	m/z	Species	RT	Start Time	End Time	Volume	Height
417.2388	418.2461	(M+H)+	0.716	0.700	0.726	3865	3889
EIC Compo	ound Search						
Mass	m/z	Species	RT	Start Time	End Time	Area	Area %
417.2376	418.2449	(M+H)+	0.713	0.703	0.739	3483	100.00
Sample Co	mparison Resu	lts					
Qualified	Changed	Resp. Ratio	Corr. RT	Normalized Hei	ght		
	New						
Isotopic Pa	attern Matching	Results					
Qualified	Score	Delta m/z					
	91.50	0.00					
Mass Defe	ct Filter Results						
Qualified	Delta Mass [mDa	]					
	-8.97						
Formula R	esults						
Assigned	Neutral Formula	Calc	. Mass De	elta Mass [mDa]	Delta M	Mass [ppm]	Calculation Base
	C21H31N5O4	417.2	2376 -1.	20	-2.87		MfeCompoundMsSpectrum
Biotransfor	mation Results	5					
Assigned	Name	Phase Offse	et Formula	Delta Mass	[mDA]	Delta Mass	[ppm] Calc. Mass
V	2x Hydroxylation	+02		1.20		2.87	417.2376

#### Figure 6

Detailed metabolite report for dihydroxy metabolite of buspirone at retention time 0.71 min. This part of the report gives detailed information about the identified metabolite and the identifying algorithms. Other detailed information about formula (see figure 7), chromatograms and isotopic pattern (see figure 8) are also available. The calculation of the formula was outlined in the detailed formula report (figure 7).

The EIC of m/z 418.24 showed about five significant peaks for possible dihydroxylated metabolites of buspirone with the selected peak at 0.71 minutes (figure 7A). The ECC showed the extracted MFE compound for the molecular mass of 417.2388 at retention time 0.71 identical to the EIC (figure 7B). The measured and calculated isotopic pattern of this compound is shown in figure 7C.

me 2x Hydro	oxylation	BioXF	Name	2x Hydroxy	lation					
rmula C21H31	N504	Mass		417.2388						
/z 418.246 ormula Sumi	1 mary	RT		0.716						
Selected	Score Formula	1	lon	Formula		Mass	Calc.	Mass	∆ Mass	[ppm]
TRUE	100.0 C21H31	N5O4	C21	132N5O4		417.2388	417	.2376		-2.87
Formula (M) C21H31N5O4	s <u>talls</u> S	elected								
Formula (M) C21H31N5O4 Species (M+H)+ Formula Rest	stans s T n 4 ults	elected RUE n/z 18.2461								
Formula (M) C21H31N5O4 Species (M+H)+ Formula Rest	s T n 4 Jlts	elected RUE n/z 18.2461 Score	Mas	s ∆Ma	iss [mDa]	∆ Mass	[ppm]	DBE		
Formula (M) C21H31N5O4 Species (M+H)+ Formula Resu Ion Formula C21H32N5O4 Isotopic Peak	s s n ults k Information	elected RUE 1/z 18.2461 Score 100.0	<b>Mas</b> 417.238	ss∆Ma ≋8	iss [mDa] -1.20	∆ Mass	<b>[ppm]</b> -2.87	DBE 9		
Formula (M) C21H31N5O4 Species (M+H)+ Formula Resu lon Formula C21H32N5O4 Isotopic Peak Abund	r <u>tans</u> S T 4 ults K Information % Calc A	elected RUE 1/z 18.2461 Score 100.0	<b>Mas</b> 417.238	s ∆Ma 88 mr/z	ss [mDa] -1.20 Calc m	∆Mass /z ∆	[ppm] -2.87 m/z [ppi	DBE 9 m]		
Formula (M) C21H31N504 Species (M+H)+ Formula Resu Ion Formula C21H32N504 Isotopic Peal Abund 100	rtans S T ults k Information % Calc A .00	elected RUE n/z 18.2461 \$Core 100.0 bund% 100.00	Mas 417.238 41	ss ∆Maa ≌8 mr/z 8.2461	ss [mDa] -1.20 Calc m 418.24	∆ Mass /z ∆ 19	[ppm] -2.87 m/z [ppn -2.2	<b>DBE</b> 9 <b>m]</b> 86		

#### Figure 7

Detailed metabolite report about the formula, including isotopic pattern, calculated for dihydroxy metabolite of buspirone at retention time 0.71 min.



#### Figure 8

Detailed metabolite report for dihydroxy metabolite buspirone at retention time 0.71 min:

A) Extracted Ion Chromatograms (EIC) of compounds with mass 418.24

B) Extracted Compound Chromatogram (ECC) of dihydroxy metabolite of buspirone at retention time 0.71 min

C) Measured and calculated isotopic pattern of dihydroxy buspirone metabolite at retention time 0.71 min.

# **Conclusion**

This Application Note demonstrated the use of the Agilent 1290 Infinity LC system with an Agilent Q-TOF LC/MS system for fast separation and accurate mass measurement of compounds in an in-vitro metabolite sample under high-throughput conditions. The metabolite compounds were separated in a run time below one minute and the width of the peaks extracted by the Metabolite ID software were below one second (FWHH). The major metabolites were identified quickly by means of the Agilent Metabolite Identification software. A summary report as well as detailed reports for each metabolite were generated.

# **References**

## 1.

E. Naegele, F. Wolf, U. Nassal, R. Jäger, H. Lehmann, F. Kuhlmann, K. Subramanian, "An interwoven, multi-algorithm approach for computerassisted identification of drug metabolites", *Agilent Technologies Application Note, publication number 5989-7375EN*, **2007**.

www.agilent.com/chem/metid

© Agilent Technologies, Inc., 2009 Published May 1, 2009 Publication Number 5989-9924EN

