

Polyisocyanate Analysis on Agilent PLgel MIXED-E with Gel Permeation Chromatography

Application Note

Materials Testing and Research, Polymers

Author

Graham Cleaver
Agilent Technologies, Inc.

Introduction

Polyisocyanates are made from multiple isocyanate functional groups. They are used as raw materials for the manufacture of polyurethane insulation, coatings, adhesives and sealants.

Analysis of polyisocyanate analysis is easy with gel permeation chromatography using Agilent PLgel 3 μ m MIXED-E columns. These columns are ideal for low molecular weight samples that contain oligomeric fractions, as well as polymers, up to 30,000 MW.





Analysis of a polyisocyanate

The chromatogram in Figure 1 shows the separation of a low molecular weight polyisocyanate using two Agilent PLgel 3 μ m MIXED-E columns in series.

Conditions

Column 2 × Agilent PLgel 3 µm MIXED-E, 300 × 7.5 mm

(p/n PL1110-6300)

Eluent Dichloromethane
Flow rate 1.0 mL/min
Detector RI

System Agilent PL-GPC 50




Figure 1. The separation of a low molecular weight polyisocyanate on an Agilent PLgel 3 µm MIXED-E two-column set.

For More Information

For more information on our products and services, visit our Web site at www.agilent.com/chem.

www.agilent.com/chem

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc., 2015 Printed in the USA April 30, 2015 5990-7984EN

