

Preparative Gel Permeation Chromatography of Phenol Distillates with Agilent PLgel

Application Note

Materials Testing and Research

Author

Graham Cleaver Agilent Technologies, Inc.

Introduction

Preparative gel permeation chromatography (GPC) is used to isolate components after suitable chromatographic conditions have been identified by running the same sample with analytical columns. A phenol sample is used here to demonstrate the method with analytical and preparative Agilent PLgel columns having the same packing.

Phenol Distillate Analysis

The sample of phenol (peak B) was found to contain contaminants (peaks A and C). For identification, analytical columns (Figure 1) were used to establish appropriate chromatographic conditions that were then used for the preparative separation (Figure 2).

Conditions - Analytical

Columns 2 × Agilent PLgel 10 µm 100Å, 7.5 × 300 mm

(p/n PL1110-6120)

Eluent Acetone

Flow rate 1.0 mL/min

Sample conc 1.0%

Inj vol 200 μL

Detector R

System Agilent PL-GPC 50








Figure 1. Preliminary analytical separation of a phenol distillate sample on Agilent PLgel 10 µm 100Å, 7.5 x 300 mm columns.

Conditions - Preparative

Columns $2 \times Agilent PLgel 10 \mu m 100 Å, 25 \times 300 mm$

(p/n PL1210-6120)

Eluent Acetone

Flow rate 1.0 mL/min

Sample conc 10%

Inj vol 2000 μL

Detector R

System PL-GPC 50




Figure 2. Preparative separation of the same phenol distillate sample on Agilent PLgel 10 μ m 100Å, 25 x 300 mm columns.

Conclusion

Agilent PLgel analytical and preparative columns can be used to investigate the composition of complex materials such as phenol distillates by gel permeation chromatography.

For More Information

These data represent typical results. For more information on our products and services, visit our Web site at www.agilent.com/chem.

www.agilent.com/chem

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc., 2015 Printed in the USA April 30, 2015 5990-8421EN

