

## **Analysis of Polystyrene**

### **Application Note**

Materials Testing and Research, Polymers

#### **Author**

Graham Cleaver Agilent Technologies, Inc.

#### Introduction

A set of three Agilent PLgel MiniMIX-B columns is ideal for routine analyses of polystyrene in tetrahydrofuran.





This sample has Mw=250,000 and Mn=100,000, and some low molecular weight components are detected with a UV at 254 nm.

PLgel 10 µm MiniMIX-B columns are designed for high MW polymer analysis and demanding eluent conditions. The PLgel 10 µm MiniMIX-B spans a wide range of molecular weights, up to 10 million, with a linear calibration curve. It is particularly useful for molecular weight distributions where slightly higher than average MWs are encountered. The 10 µm particle size provides good resolution with relatively low pressures for enhanced lifetimes in demanding conditions.

# emanding conditions. 10 min Figure 1. Analysis of polystyrene using PLgel 10 µm MiniMIX-B columns

#### **Conditions**

Column: 3 x PLgel 10 µm

Mini-MIX-B, 300 x 7.5 mm (part number PL1510-5100)

Eluent: THF
Flow Rate: 0.3 mL/min
Detection: UV, 254 nm

#### www.agilent.com/chem

This information is subject to change without notice.

© Agilent Technologies, Inc. 2015

Published in UK, April 30, 2015

5991-5810EN

30

