

WATERS OASIS® FAMILY OF SAMPLE EXTRACTION PRODUCTS

Expect high performance SPE using the right chemistry, format, and methodology. Now you can achieve robust SPE methods without having to endure irreproducible results and low recoveries caused by: undesirable silanol activity, sorbent drying, pH limitations, and breakthrough of polar compounds/metabolites.

Oasis® HLB sorbent enables a single, fast, and generic SPE protocol for rapid method development. Oasis® MCX and MAX sorbents provide highly selective and sensitive SPE methods for basic and acidic compounds, respectively. You can also clean up your samples on-line with the Oasis® on-line columns and Prospekt™ cartridges. The innovative design of the new Oasis® μ Elution plate enables you to elute analytes with as little as 25 μ L of solvent. This patented technology further increases your productivity and sensitivity.

No matter what your challenge, if you can do it by SPE, Oasis® can do it better.

