

# Polyester Polyol Analysis on Agilent PLgel 5 µm using Gel Permeation Chromatography

# **Application Note**

Materials Testing and Research, Polymers

### **Author**

**Graham Cleaver** 

Agilent Technologies, Inc.

# Introduction

Polyols are alcohols containing multiple hydroxyl groups. They are mainly used as reactants to make polymers such as polyester polyol, by condensation or step-growth polymerization of a diol polyol with dicarboxylic acid. Polyester polyols are themselves reacted with polyisocyanates to make polyurethanes for rigid-foam, flame-retardant building board. Natural oil polyester polyols from vegetable oils are beginning to replace some epoxide-based polyols.

Analysis of polyester polyols is easy with gel permeation chromatography and an Agilent PLgel column.





# **Analysis of Polyester Polyol**

Figure 1 shows the excellent oligomer resolution achieved for a relatively low molecular weight polyester polyol.




Figure 1. Excellent separation of oligomers in a polyester polyol using Agilent PLgel columns.

# **Conditions**

Column 2 × Agilent PLgel 5 µm 500Å, 300 × 7.5 mm

(p/n PL1110-6525)

Eluent THF

Flow rate 1.0 mL/min

Detector RI

System Agilent PL-GPC 50

### Conclusion

Low-pore-size PLgel columns can resolve oligomeric content from polymer samples, an important component that determines the behaviour of many polymers.

### For More Information

These data represent typical results. For more information on our products and services, visit our Web site at www.agilent.com/chem.

## wwww.agilent.com/chem

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc., 2015 Printed in the USA April 30, 2015 5990-8325EN

